

NB-326: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm2, 20 °C (reproduction NB-325)

Date: 2025-10-28

Tags: O2 Test Calibration NB Firing
O2 sensor H2 SrTiO3 Unisense
RhCrO3:Al:SrTiO3 H2 Sensor
temperature In situ Trace range robust
oxygen sensor photocatalysis Unisense
normal range

Category: SrTiO3


Status: Done

Created by: Nadzeya Brezhneva

Objectives

Reproduction of NB-325: simultaneous detection of H₂ and O₂ evolution in liquid phase for irradiated suspension of Rh_xCrO_y:Al:SrTiO₃ suspension (EA-358, 0.5 mg/mL), 365 nm LED, 20 mW/cm², 20 °C.

Reaction scheme



ChemDraw file linked: [NB-326-SrTiO3-photocatalytic H2O splitting.cdxml](#)

Literature/reference experiments

Literature	/
Reproduction	SrTiO3 - NB-325: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm2, 20 °C
Similar experiments	SrTiO3 - NB-316: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C

Reagents

Name	CAS Number / Experiment Number	Inventor y number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	Pressure [bar]	Concentration [mM]
milli-Q H ₂ O	/	/	/	/	/	/	/	0.998	25 + 25 (for calibration)	/	/
Al:SrTiO3 RhCrOx (EA-358)	SrTiO3 - EA-358: Modification of Al:SrTiO3 (EA-354) via deposition of Rh, Cr oxide co-catalyst, 350°C, 1h, Upscaling (3.33x)	/	/	/	12.50	12.72	/	/	/	/	/

Hydrogen	1333-74-0	/	/	/	/	/	/	/	2 balloons (ca. 2 L)	approx. 1	/
----------	-----------	---	---	---	---	---	---	---	-------------------------	--------------	---

Excel sheet for reagent calculation

/

Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-324: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup I I](#)

	Name
Used Set-up	Equipment - Advanced irradiation setup V1.0 I
Irradiation setup number	Equipment - Irradiation setup 4 (CEEC II, E002)

	Light Source Name	Power Source Name	Wavelength [nm]	Power Setting [mW]	Analog Setting [0.00 - 10.00]
First light source	Light Source - UHP LED 365 nm-4	Power Sources - BLS-18000-1 4	365	22	0.07

Used beam combiner [Name or None]	/
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	20
Stirring speed [rpm]	500
Irradiation start: 1. Firesting [relative to start log] 2. Unisense	1. 603 s 2. 17:50:27

Irradiation stop: 1. Firing [relative to start log] 2. Unisense	1. 1534 s 2. 18:05:57
--	--------------------------

O₂/H₂ sensor equipment

	Equipment	Used protocol
Used Firing	Equipment - Firing Fiber-Optic Oxygen Meter 2 Channel (Firing 2)	Protocol - Operation of Firing Fiber-Optic Oxygen Meter 2 Channel Software
Used O ₂ sensor	Equipment - Robust probe for liquid O ₂ measurement	Protocol - In-situ hydrogen and oxygen measurement in H ₂ /O ₂ reactor
Used H ₂ sensor	Equipment - H ₂ UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurement in H ₂ /O ₂ reactor

Procedure/observations

Date	Time	Step	Observations	Pictures/Files
28.10.2025		The experiment was done according to Protocol - In-situ hydrogen and oxygen measurement in H₂/O₂ reactor Important steps and deviations are listed below	/	/
	10:51-11:57	Conditioning of H ₂ sensor	NB-326-Logger1 1.8 mV at the end of polarization procedure	NB-326.ulong NB-326-Logger1.csv NB-326-Logger1.bmp
	ca. 14:30	Assembling the setup for calibration (25 mL of water was added using graduated cylinder), done according to Protocol - Liquid phase calibration of H₂ UniAmp sensor with H ₂ bubbling.	/	/
	14:53	Start of O ₂ logging.	NB-326-Ch2-1	2025-10-28_145357_NB-326-Ch2-1.txt 2025-10-28_145357_NB-326-Ch2-1.png

	14:54	Start of H2 logging.	NB-326-Logger2 offset -2 mV	NB-326.ulog NB-326-Logger2-calibration.csv NB-326-Logger2-calibration.bmp
	14:57	Degassing was started.	/	20251028_145926-degassing of water.jpg
	15:19	Introducing H2 sensor into the reactor under Ar flow.	/	/
	15:22	0 ppm was taken.	/	/
	15:30	H2 bubbling of the reactor was started	/	20251028_153053-purging with H2.jpg
	15:39	Introducing second H2 balloon.	/	20251028_154109-with 2nd H2 balloon.jpg
	15:50	1.000.000 ppm point was taken and calibration was saved	802 mV, slope: 1.026, 781.887 uM	20251028_155003-H2 table.jpg
	15:52	Stop of H2 logging.	/	/
	15:52	Stop of O2 logging.	/	/
	ca. 16:00	Deassembling the setup, drying the reactor with acetone and compressed air .	/	/
		Sample preparation		
	16:20	Weighing EA-358 photocatalyst in a 50 mL vial.	Creamy solid.	/
	16:25	Addition of 25 mL H2O to the vial via graduated cylinder.	/	/
	16:28-31	The suspension was vortexed for 3 min (Equipment - VWR® VV3, Vortex Mixer , stage 4/6), covered with Al foil before further use.	/	20251028_163114-suspension after vortex.jpg
		Continue in Protocol - In-situ hydrogen and oxygen measurment in H2/O2 reactor from step 6		
	16:40-50	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	/
	ca. 16:55	Assembling the setup.	Currently, stopper instead of H ₂ sensor, PT100, PT1000 and O ₂ robust probe are inside the reactor immersed in the liquid phase	20251028_165756-before degassing and irradiation.jpg

	17:00	Start of O2 logging.	NB-326-Ch2-2	2025-10-28_170045_NB-326-Ch2-2.txt 2025-10-28_170045_NB-326-Ch2-2.png
	17:03	The degassing was started	/	20251028_170354-degassing of the suspension.jpg
	17:32	Cannula was transferred to gas phase, above the suspension.	/	/
	17:35	H ₂ sensor was added in Ar counterflow.	/	/
	17:38	The degassing was stopped by removing the cannula and closing the valve.	/	20251028_174044-after introducing H2 sensor.jpg
	17:40	Stop of O2 logging.	/	/
	17:40	Start of O2 logging.	NB-326-Ch2-3	2025-10-28_174023_NB-326-Ch2-3.txt 2025-10-28_174023_NB-326-Ch2-3.png
	17:40	Start of H2 logging.	NB-326-Logger3	NB-326.ulong NB-326-Logger3-during irradiation.csv NB-326-Logger3-during irradiation.bmp
	17:40-50	Equilibration time.	/	/
	17:50	The irradiation was started	/	20251028_175536-after start of irradiation.jpg
	18:05	The irradiation was stopped.	/	/
	18:05-15	Equilibration time.	/	/
	18:16	Stop of O2 and H2 logging.	/	/

	ca. 18:20	Deassembling the setup, cleaning the reactor.	Tips of the sensors and reactor were covered with attached photocatalyst particles. Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode).	20251028_182055-after irradiation.jpg
--	-----------	---	--	---

Analysis

Used calibration for Firesting: [20250611-BOLA-fitting-liquid-phase-trace-oxygen-sensor-H2-O2 reactor.ini](#)

Used calibration for UniSense: NB-326-Logger2

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
28.10.2025	10:51	NB-326-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-326.ulong NB-326-Logger1.csv	/	NB-326-Logger1.bmp	/	
	14:54	NB-326-Logger2	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-326.ulong NB-326-Logger2-calibration.csv	/	NB-326-Logger2-calibration.bmp	/	
	17:40	NB-326-Logger3	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-326.ulong NB-326-Logger3-during irradiation.csv	NB-326-O2 and H2 curve.py	NB-326-Logger3-during irradiation.bmp NB-326-O2 and H2 curves.png	SrTiO3 - NB-325: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm2, 20 °C	
	14:53	NB-326-Ch2-1	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-10-28_145357_NB-326-Ch2-1.txt	/	2025-10-28_145357_NB-326-Ch2-1.png	/	
	17:00	NB-326-Ch2-2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-10-28_170045_NB-326-Ch2-2.txt	/	2025-10-28_170045_NB-326-Ch2-2.png	/	
	17:40	NB-326-Ch2-3	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-10-28_174023_NB-326-Ch2-3.txt	NB-326-O2 and H2 curve.py	2025-10-28_174023_NB-326-Ch2-3.png NB-326-O2 and H2 curves.png	SrTiO3 - NB-325: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm2, 20 °C	

Results

Simultaneous H₂ and O₂ measurements of irradiated suspension of EA-358 (0.5 mg/mL) in O₂/H₂ photoreactor under 365 nm irradiation (20 mW/cm², 20 °C, 15 min) were performed.

H₂ level at the end of irradiation - ca. 17.5 umol/L, O₂ level - ca. 17.1 umol/L.

Linked experiments

SrTiO₃ - NB-316: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

SrTiO₃ - NB-318: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm², 20 °C

SrTiO₃ - NB-322: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm², 20 °C (reproduction NB-318)

SrTiO₃ - NB-325: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm², 20 °C

Linked resources

Equipment - [VWR® VV3, Vortex Mixer](#)

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)


Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)

Protocol - [In-situ hydrogen and oxygen measurment in H2/O2 reactor](#)

Attached files

NB-326-O2 and H2 curves.png

sha256: 7c183cef06a97b7e0c57f5f3a2d87ce6db7197ff40d9ad37514dff7424383c52



NB-326-O2 and H2 curve.py


sha256: b582ddfb7aaf1ea15a8639ff018d59b9e2fd9074f127e638ec44d14b270d7a18

NB-326-Logger3-during irradiation.csv

sha256: d1a9d49223ceeff90feb69bffe06fb0df4260027d1c5032438fcadcb6da1f140

NB-326-Logger3-during irradiation.bmp

sha256: 42389795d1265dfd8964984841955cc196b666bc40068f4c749792987c9196c5



NB-326-SrTiO3-photocatalytic H2O splitting.cdxml

sha256: 4f2afbd4a5fd50272e5f8e543f1f9ae749d102033a1866e51d2f0744e92a0d58

NB-326-SrTiO3-photocatalytic H2O splitting.png

sha256: feff62799df5cea631f939e338cf4de88cbeeba0ca3457ed98c623d53659a829



20251028_145926-degassing of water.jpg

sha256: 763f528ebdd3ad681d721789fc4129fe21337c30b5b26fee7d3b85fa971ea231



20251028_153053-purging with H2.jpg

sha256: ebaa92eaab62a010ee6f752702be6a7384a17df967d3d1d2746ecd0a275c1b42



20251028_154109-with 2nd H2 balloon.jpg

sha256: 0b96f9e2b2d3f2aee65b90587961e20bc6281e3bb8453759a9f760aef83f045d



20251028_155003-H2 table.jpg

sha256: 69a334bccb957553e66de1e594db51e4716ef3911f0162e1fd4dbd18fb3cc4ea



20251028_170354-degassing of the suspension.jpg

sha256: b7816f5ede79edd36df1bacc783bd54783a3793af21464cfbc13ab215a512cd4



20251028_165756-before degassing and irradiation.jpg

sha256: 61f0ccf5f9de0a535807b48a3e25ed518efb0183243ece1b76409f49024edef8



20251028_163114-suspension after vortex.jpg

sha256: 04b78b51a5fb177e43eb04ab50d9595c893399b4e459901768e026b005cbb890



20251028_174044-after introducing H2 sensor.jpg

sha256: a9217ff279c9f6f1601e99007935a2ba84935ea717561ff5368929a5550a08f8



20251028_182055-after irradiation.jpg

sha256: 29c948ebb6e5fdaca240b8cf793e5cda7aa291910c7dce9fa3835e4766d64796



20251028_175536-after start of irradiation.jpg

sha256: de36d116cd425ed3904760ba25d62192b1b37c8c4e95a9fcd382a605cfc9f561




2025-10-28_174023_NB-326-Ch2-3.txt

sha256: b68deddce42f078fdca1378bb012300a5553e151860381336e045e1b214e1e76

2025-10-28_174023_NB-326-Ch2-3.png

sha256: 8ac03e0b7966380113f5f3967250179fa484b31b703ce22d6578e8faafe03cfa



2025-10-28_170045_NB-326-Ch2-2.txt

sha256: fd78a92b7dccc03498d04eb9cccd0b684fa6685a651c4f2d5241c91b341a8edf

2025-10-28_170045_NB-326-Ch2-2.png

sha256: ce87430b62ea9150517d257771a7a8b5242daaf81e6a075ff273e8df50717ab7



2025-10-28_145357_NB-326-Ch2-1.txt

sha256: 454231e42f208e24a0a19fb11c6652515c35d14cf49c177adee6fc2467d7e319

2025-10-28_145357_NB-326-Ch2-1.png

sha256: bcc988f288d4da2844dda0e4d1f103bc238f06aabc61500be3afdb59888ddcda




NB-326-Logger2-calibration.csv

sha256: e8f2a24b370ee174768ad23ceec48b507615878a6cd09fb31783afae52157622

NB-326-Logger2-calibration.bmp

sha256: 7957305628e0dd3d9700ed6fefe7e23e3964f369a659e55ac337f8767eb7c97




NB-326-Logger1.csv

sha256: ce041f86936e727d4352c019e41ef2a740f95fc6d44a39834a90b7045906f2a0

NB-326-Logger1.bmp

sha256: 595c965fb01a8bc64fdd6831aabfb7373caea653cad7001cb28b07bbbcabd428



NB-326.ulog

sha256: a95f6508c32e78edb4b5b9e01029a1452d1a9d54cfad4fb4ada20308e18615c2



Unique eLabID: 20251028-c26ef68cef85fae37f14d3b671fb3f6f232b0228

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3288>