

# NB-327: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrO<sub>x</sub>,Al:SrTiO<sub>3</sub> (NB-321, 0.05 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 20 °C

Date: 2025-10-28

Tags: O<sub>2</sub> Test Calibration Future NB  
Firesting O<sub>2</sub> sensor H<sub>2</sub> SrTiO<sub>3</sub> Unisense  
RhCrO<sub>3</sub>:Al:SrTiO<sub>3</sub> H<sub>2</sub> Sensor  
temperature In situ Trace range robust  
oxygen sensor photocatalysis Unisense  
normal range

Category: SrTiO<sub>3</sub>


Status: Need to be redone

Created by: Nadzeya Brezhneva

## Objectives

Simultaneous detection of H<sub>2</sub> and O<sub>2</sub> evolution in liquid phase for irradiated suspension of Rh<sub>x</sub>Cr<sub>y</sub>O<sub>3</sub>:Al:SrTiO<sub>3</sub> suspension (NB-321, **0.05 wt% Rh,Cr**, 0.5 mg/mL), 365 nm LED, 50 mW/cm<sup>2</sup>, 20 °C.

## Reaction scheme



ChemDraw file linked: [NB-327-SrTiO3-photocatalytic H2O splitting.cdxml](#)

## Literature/reference experiments

Literature	/
Reproduction	/
Similar experiments	<a href="#">SrTiO3 - NB-326: Liquid phase H2 and O2 of RhCrO<sub>x</sub>,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm<sup>2</sup>, 20 °C (reproduction NB-325)</a>

## Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass <sub>theo</sub> [mg]	Mass <sub>exp</sub> [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	Pressure [bar]	Concentration [mM]
milli-Q H <sub>2</sub> O	/	/	/	/	/	/	/	0.998	25 + 25 (for calibration)	/	/
Al:SrTiO <sub>3</sub> RhCrO <sub>x</sub> (NB-321)	<a href="#">SrTiO<sub>3</sub> - NB-321: Modification of EA-354 (SrTiO<sub>3</sub>:Al, upscaled batch) with Rh, Cr oxide cocatalyst (0.05 wt%)</a>	/	/	/	12.50	12.87	/	/	/	/	/
Hydrogen	1333-74-0	/	/	/	/	/	/	/	2 balloons (ca. 2 L)	ca. 1 bar	/

# Excel sheet for reagent calculation

/

## Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-314: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup](#)

	Name
Used Set-up	<a href="#">Equipment - Advanced irradiation setup V1.0 I</a>
Irradiation setup number	<a href="#">Equipment - Irradiation setup 4 (CEEC II, E002)</a>

	Light Source Name	Power Source Name	Wavelength [nm]	Power Setting [mW]	Analog Setting [0.00 - 10.00]
First light source	<a href="#">Light Source - UHP LED 365 nm-4</a>	<a href="#">Power Sources - BLS-18000-1 4</a>	365	56	0.19

Used beam combiner [Name or None]	/
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	20
Stirring speed [rpm]	500
Irradiation start: 1. Firesting [relative to start log] 2. Unisense	1. 603 s 2. 22:08:19
Irradiation stop: 1. Firesting [relative to start log] 2. Unisense	1. 1513 s 2. 22:23:30

## O<sub>2</sub>/H<sub>2</sub> sensor equipment

	Equipment	Used protocol
Used Firesting	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)	Protocol - Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software
Used O <sub>2</sub> sensor	Equipment - Robust probe for liquid O <sub>2</sub> measurment	Protocol - In-situ hydrogen and oxygen measurment in H <sub>2</sub> /O <sub>2</sub> reactor
Used H <sub>2</sub> sensor	Equipment - H <sub>2</sub> UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurment in H <sub>2</sub> /O <sub>2</sub> reactor

## Procedure/observations

Date	Time	Step	Observations	Pictures/Files
		Calibration from SrTiO <sub>3</sub> - NB-326: Liquid phase H <sub>2</sub> and O <sub>2</sub> of RhCrOx,Al:SrTiO <sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm <sup>2</sup> , 20 °C (reproduction NB-325) was used.	/	/
28.10.2025		<b>Sample preparation</b>		
	20:30	Weighing photocatalyst in a 50 mL vial (the vial was covered with Al foil before further use).	Creamy solid	/
	21:06	Addition of 25 mL H <sub>2</sub> O to the vial via graduated cylinder.	/	/
	21:06-09	The suspension was vortexed for 3 min (Equipment - VWR® VV3, Vortex Mixer, stage 4/6), covered with Al foil before further use.	Seems less optically dense than for the suspension with 0.1 wt% Rh,Cr catalyst EA-358	20251028_210937-suspension after vortex.jpg
		Continue in Protocol - In-situ hydrogen and oxygen measurment in H <sub>2</sub> /O <sub>2</sub> reactor from step 6		
	ca. 21:15-20	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	/
	21:20	Assembling the setup.	Currently, stopper instead of H <sub>2</sub> sensor, PT100, PT1000 and O <sub>2</sub> robust probe are inside the reactor immersed in the liquid phase	20251028_212335-before degassing.jpg

	21:25	Start of O2 logging.	<b>NB-327-Ch2-1</b>	2025-10-28_212536_NB-327-Ch2-1.txt 2025-10-28_212536_NB-327-Ch2-1.png
	21:30	The degassing was started	/	20251028_213111-suspension during degassing.jpg
	21:52	Cannula was transferred to gas phase, above the suspension.	/	/
	21:54	H <sub>2</sub> sensor was added in Ar counterflow.	/	/
	21:57	The degassing was stopped by removing the cannula and closing the valve.	/	/
	21:58	Stop of O2 logging.	/	/
	21:58	Start of O2 logging.	<b>NB-327-Ch2-2</b>	2025-10-28_215816_NB-327-Ch2-2.txt 2025-10-28_215816_NB-327-Ch2-2.png
	21:59	Start of H2 logging.	<b>NB-327-Logger1</b>	NB-327.ulong NB-327-Logger1.csv NB-327-Logger1.bmp
	21:59-22:09	Equilibration time.	/	/
	22:09	The irradiation was started	/	20251028_220831-after start of irradiation.jpg
	22:23	The irradiation was stopped.	/	/
	22:23-33	Equilibration time.	/	/
	22:33	Stop of O2 and H2 logging.	/	/

	ca. 22:40	Deassembling the setup, cleaning the reactor.	Tips of the sensors and reactor were covered with attached photocatalyst particles. Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode).	/
--	-----------	---	--	---

## Analysis

Used calibration for Firesting: [20250611-BOLA-fitting-liquid-phase-trace-oxygen-sensor-H2-O2 reactor.ini](#)

Used calibration for UniSense: NB-326-Logger2

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
28.10.2025	21:59	NB-327-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-327.ulong NB-327-Logger1.csv	NB-327-O2 and H2 curve.py	NB-327-Logger1.bmp NB-327-O2 and H2 curves.png	/	Clean response signal, slight decrease at 22:21. H2 value 34 uM H2
	21:25	NB-327-Ch2-1	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-10-28_212536_NB-327-Ch2-1.txt	/	2025-10-28_212536_NB-327-Ch2-1.png	/	Degassing of the suspension.
	21:58	NB-327-Ch2-2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-10-28_215816_NB-327-Ch2-2.txt	NB-327-O2 and H2 curve.py	2025-10-28_215816_NB-327-Ch2-2.png NB-327-O2 and H2 curves.png	/	A bit sharp increase in O2 value after start of irradiation, sharp decrease after stop of irradiation. O2 value 32 uM.

## Results

Simultaneous H<sub>2</sub> and O<sub>2</sub> measurements of irradiated suspension of NB-321 (0.05 wt% RhCr, 0.5 mg/mL) in O<sub>2</sub>/H<sub>2</sub> photoreactor under 365 nm irradiation (50 mW/cm<sup>2</sup>, 20 °C, 15 min) were performed.

H<sub>2</sub> level at the end of irradiation - 34 umol/L, O<sub>2</sub> level - 32 umol/L.

## Future recommendations

Old procedure	Problem	Suggested new procedure
/	Sharp increases in O2 and H2 signals - maybe related either with the position of the sensors in the reactor or the optical density.	<ul style="list-style-type: none"> <li>• Adjust the position of the sensors relatively to the light beam.</li> <li>• If necessary, destroy the lumps of the solid photocatalyst before adding water and vortexing to improve dispersion quality.</li> </ul>

## Linked experiments

SrTiO<sub>3</sub> - NB-316: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrO<sub>x</sub>,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 20 °C

SrTiO<sub>3</sub> - NB-318: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrO<sub>x</sub>,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm<sup>2</sup>, 20 °C

SrTiO<sub>3</sub> - NB-321: Modification of EA-354 (SrTiO<sub>3</sub>:Al, upscaled batch) with Rh, Cr oxide cocatalyst (0.05 wt%)

SrTiO<sub>3</sub> - NB-325: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrO<sub>x</sub>,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm<sup>2</sup>, 20 °C

SrTiO<sub>3</sub> - NB-326: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrO<sub>x</sub>,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm<sup>2</sup>, 20 °C (reproduction NB-325)

## Linked resources

Equipment - [VWR® VV3, Vortex Mixer](#)

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)

Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)


Protocol - [In-situ hydrogen and oxygen measurment in H<sub>2</sub>/O<sub>2</sub> reactor](#)

## Attached files

NB-327-O<sub>2</sub> and H<sub>2</sub> curve.py  
sha256: c76c6611dbbfe05cf7f921b2ece8982fb318b0a3e7b07134353113e71d9b00d0


NB-327-O<sub>2</sub> and H<sub>2</sub> curves.png

sha256: 698e5a642ce9a534ba48b15941e4cffd92107aeee3d418d5b4c8d2ab865810c9



NB-327-Logger1.bmp

sha256: f00f292145d1990a14a66fb2df5e8f8fe505c7c59e4e30233e35cddee01f4c36



NB-327-Logger1.csv

sha256: 78136f463d41603004857f1d1aa8c90064c60d2be0bd80cc3da8f4a334cd5751

NB-327-SrTiO3-photocatalytic H2O splitting.cdxml

sha256: a0a48b67ad0eff08e8d3187e4e962dc7a60f1575f047ff7a0ba610e2824d18d3

NB-327-SrTiO3-photocatalytic H2O splitting.png

sha256: 4ae3c4443899250e9b9f241542fe3f421c490189250c5c9fcc5fb5d79700e732



20251028\_212335-before degassing.jpg

sha256: 4bc2681585b5a6c8b428207dedbae18e186517dbc83d9a11f301250aa29df478



20251028\_213111-suspension during degassing.jpg

sha256: 5ef56ce8cdc32bc713e0bc9906bd984dbff5f4bf9d7182d3f99ba8b0ddc6fe62



20251028\_210937-suspension after vortex.jpg

sha256: a77974bfae808c8c6c05a519df75424bcb2290087a5f46123c502a75ee69adc3



20251028\_220831-after start of irradiation.jpg

sha256: 23394a56b25eb170503aafef5004625f9a8cdf42fa2263f2373b2a6382b68dff




NB-327.ulog

sha256: d2da7b55c77aa27fb739934720b1689b2c3c85220b3b4df0ccd5828557163ebc

2025-10-28\_215816\_NB-327-Ch2-2.png

sha256: 697cdc71c57b993da53fc43efc325ee3a462bafd338e4f05449bdd9e65463005




2025-10-28\_215816\_NB-327-Ch2-2.txt

sha256: e5fa4057b0115e976ca961e33046d33f550d17af925b21061230cfb8a382f44c

2025-10-28\_212536\_NB-327-Ch2-1.png

sha256: 88250e0a4d3802430c52b7f99bd04fe38a1521746039445c72294123a6d0dbf5



2025-10-28\_212536\_NB-327-Ch2-1.txt

sha256: 02373451af748dd2d83ca77bae2fa651c3318ab84d20bad20cf11eeebaa49053



Unique eLabID: 20251028-bce0ec20c997b46418aac31d6927a3690f055aa0  
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3289>