

# Gas phase calibration of H2 UniAmp sensor

**Date:** 2025-05-12  
**Tags:** AE H2 H2 Evolution H2 Sensor  
**Category:** Protocol  
**Created by:** Alexander Eith

## Goal

2 point gas phase calibration of [Equipment - H2 UniAmp Sensor](#) with [Equipment - H2 UniAmp Single Channel System](#)

## Prerequisites and preparation

Everything listed in [Protocol - Getting hydrogen from hydrogen bottle in CEEC II E014](#)

[Equipment - H2 UniAmp Sensor](#)

[Equipment - H2 UniAmp Single Channel System](#)

[Prep work - AE-482: Volume determination of GL14/NS14 5 mL Schlnek flask](#)

100 µL Hamilton syringe

BOLA fitting (2 mm hole)

Milli-Q water

## Safety

Everything stated in [Protocol - Getting hydrogen from hydrogen bottle in CEEC II E014](#)

Be careful with sensor, sensitive to touch

Polarization needs to be 100 mV (will be adjusted automatically when connected to [Equipment - H2 UniAmp Single Channel System](#) and [Equipment - H2 UniAmp Sensor](#)

## Steps

Step number	Step description	Pictures
0. Build setup I	Start the Laptop Plug it in (is always needed for grounding of the setup) Connect the <a href="#">Equipment - H2 UniAmp Single Channel System</a> with a USB-C USB-A to the laptop (Light on the device will turn on an green) Connect the <a href="#">Equipment - H2 UniAmp Sensor</a> to the device (Light on the device will turn off for 100 ms) Connect the PT100	

1. Start software	<p>If not installed: Install from <a href="https://unisense.com/product-category/software/">https://unisense.com/product-category/software/</a></p> <p>Open software (Unisense logger) on Laptop <a href="#">icon.png</a></p> <p>You get to a view of connected sensors. Should look somewhat like the picture <a href="#">start page of software.png</a></p> <ul style="list-style-type: none"> <li>• In column use: select if sensor reading should be measured during experiment</li> <li>• In column unit for the hydrogen sensor the unit can be changed (at the moment: <math>\mu\text{mol/L}</math> for liquid phase)</li> </ul> <p>Press start experiment</p> <p>Data is saved in Documents / Unisense Data --&gt;</p> <p>There save experiment in a folder with your abbreviation, give experiment a name <a href="#">save file.png</a></p>	<a href="#">icon.png</a> <a href="#">start page of software.png</a> <a href="#">save file.png</a>
2. Measurement configuration	<p>In software go to Calibration:</p> <ul style="list-style-type: none"> <li>• check in field UniAmp Channel Configuration that Polarisation is set to 100 mV - <b>do never change that value</b> if value is different immediately stop experiment, disconnect everything and inform AE</li> <li>• Preamp Gain is typically set to <math>1\text{pA}=1\text{mV}</math>, before changing it ask AE</li> <li>• Offset can be adjusted, so that signal is as close to 0 mV as possible</li> <li>• select solubility temperature compensation and sensor temperature compensation, when PT100 is connected</li> </ul>	<a href="#">Channel Configuration.png</a>
2a. Conditioning	<p>In software go to Data Logger:</p> <p>set sample interval to 10 s and start measurement (green triangle below sampling interval)</p> <p>In top of view sensors it can be selected which sensors are shown below, on the right some options are shown to play around with the visualization</p> <p>Wait until signal is stable (for start wait at least 20 min) - <b>do not plug out anything, since then conditioning has to be repeated</b></p>	<a href="#">Measurement 1.png</a>
3. Preparing Hydrogen saturated water	<p>Prepare a 5 mL flask with approx. 5 mL Milli-Q water and a septum. Use double balloons (speak balloon inside another balloon) and use two balloon volumes for <math>\text{H}_2</math> bubbling.</p> <p>Get hydrogen according to <a href="#">Protocol - Getting hydrogen from hydrogen bottle in CEEC II E014</a> - <b>follow stated safety measures</b></p> <p>Flush the flask with hydrogen by placing the balloon with cannula into the water and allow the gas to be purged with another cannula not connected to anything</p> <p>Let bubble for at least 20 min</p> <p>When balloon is nearly empty the cannula not connected to anything can be removed</p>	

4. Build setup II	<p>After step 2a and 3 are finished:  Stop experiment in software by pressing green square or by closing the software  Build the NS14/GI14 flask with a valve and a septum on the NS14  During that add milli-q water (1 mL) with a Eppendorf pipette to the flask. Also add a stir bar and start the stirring  Remove the protective plastic hull from the sensor  Add sensor to the flask using a 2.0 mm BOLA fitting - <b>be careful to not damage the BOLA fitting or the sensor</b>  Adjust sample interval to 1 s  Start new measurement by pressing on the green triangle  Wait until signal is stable (typically in less than 5 min)</p>	
5. Take calibration point at 0 ppm	<p>In software: go to calibration  In field "Sensor calibration &amp; experiment settings":  Adjust Temp. in drop down menu next to it to PT1000 when PT1000 is used, otherwise adjust temperature directly to desired temperature  Press on field "clear all points"  Press button "Add point"  Go to "Comment" and write there, when calibration point was taken or take note manually</p>	<a href="#">calibration.png</a>
6. Take calibration point with hydrogen inside	<p>Click on field "H2 table" and put in amount of hydrogen in ppm, as well as the temperature and pressure measured from the device  Add desired amount (<math>100\text{ }\mu\text{L} = 9381\text{ ppm} = 0.938\%</math>) of water saturated hydrogen (gas phase of flask) prepared in step 3 (should give concentration in the end similar to highest expected hydrogen concentration of experiment) with the gas tight hamilton syringe - close the syringe while transferring hydrogen  Wait till maximal hydrogen concentration is reached, then press add point. Note done when hydrogen was added and when calibration point was taken</p>	
7. Save calibration	<p>Press field "Save and use calibration"  Check at bottom of window, that correct calibration is now selected (most recent one)  Stop experiment in software by pressing green square in Data Logger</p>	
8. Dissassemble setup	<p>Remove BOLA fitting and sensor from flask  Do not disconnect from device and so on</p>	

9. Cleaning of sensor	<p>Sensor should not be dirty - so no cleaning needed, but it takes several minutes for the signal to go back to baseline value</p> <p>If no other experiments are planned for the day: Make sure that the sensor is actually dry and place it back in the protective hull and then in the storage box</p>	
10. Save data	<p>Data is automatically saved in format .ulog, which can always be opened with the software from unisense</p> <p>Data can be exported as excel file, .csv or picture from Data Logger</p> <p>In experiment note down relevant obtained values from calibration (slope and intercept) as well as H<sub>2</sub> concentration used for 2nd calibration point</p>	

## Linked resources

Equipment - [EPR, CEEC I lab K002](#)

Equipment - [H<sub>2</sub> UniAmp Sensor - Low range - 2.1 x 80 mm needle](#)

Equipment - [H<sub>2</sub> UniAmp Single Channel System](#)

Protocol - [EPR measurement on EPR, CEEC I lab K002](#)

Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)

Protocol - [Liquid phase calibration of H<sub>2</sub> UniAmp sensor](#)

## Attached files

calibration.png

sha256: 1aee43ae9c51f6ab16da305a1cb2514edac3d8835c3224cc131101789bb60baf



Measurment-1.png

sha256: 1ef332f79186f52c4f1787d4736abe44134077bdbaddef7022d43da623538887



Channel-Configuration.png

sha256: e6356a65c278ad48ea8c0064e75060804576a52121045317b4f18636ee752ddd



save-file.png

sha256: 17bf002499bf45a5990d17dec2be25270e397437359f030238ae4af8f2ebc75



start-page-of-software.png

sha256: d44ecd9a4cc534851b2236a78fb3deb14f6649813a13b94355a68f3a0bac84e7



icon.png

sha256: 20103d3125a59f726b00e54e5d70b13d6052dd22ed9390da9a4522e4d3da4e13



## Comment

On 2025-05-27 19:08:42 Nadzeya Brezhneva wrote:

Notes:

Future recommendation: placing water into NS14/GL14 Schlenk where H2 measurement is performed



Unique eLabID: 20250512-b5af754f312800aa983ea02c2d448b8753fc93e6

Link: <https://elab.water-splitting.org/database.php?mode=view&id=225>