

NB-356: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C IX (reproduction NB-316)

Date: 2025-11-20

Tags: O₂ Test Calibration Future NB
Firing O₂ sensor H₂ SrTiO₃
troubleshooting Unisense
RhCrO₃:Al:SrTiO₃ H₂ Sensor
temperature In situ Trace range robust
oxygen sensor photocatalysis Unisense
normal range

Category: SrTiO₃


Status: Done

Created by: Nadzeya Brezhneva

Objectives

Reproduction of NB-316: simultaneous detection of H₂ and O₂ evolution in liquid phase for irradiated suspension of Rh,CrO_x:Al:SrTiO₃ suspension (EA-358 sample, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C, 15 min (reference conditions).

Reaction scheme



ChemDraw file linked: [NB-316-SrTiO3-photocatalytic H2O splitting.cdxml](#)

Literature/reference experiments

Literature	/
Reproduction	<p>SrTiO₃ - NB-320: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C (reproduction NB-316) I</p> <p>SrTiO₃ - NB-329: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C (reproduction NB-316) I I</p> <p>SrTiO₃ - NB-331: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C (reproduction NB-316) I I I</p> <p>SrTiO₃ - NB-336: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C V (reproduction NB-316)</p> <p>SrTiO₃ - NB-339: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C VI (reproduction NB-316)</p> <p>SrTiO₃ - NB-348: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C VII (reproduction NB-316)</p>
Similar experiments	<p>SrTiO₃ - NB-315: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (NB-301, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C</p>

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	Pressure [bar]	Concentration [mM]
milli-Q H ₂ O	/	/	/	/	/	/	/	0.998	25 + 25	/	/
Al:SrTiO ₃ RhCrO _x (EA-358)	SrTiO ₃ - EA-358: Modification of Al:SrTiO ₃ (EA-354) via deposition of Rh, Cr oxide co-catalyst, 350°C, 1h, Upscaling (3.33x)	/	/	/	12.5	12.60	/	/	/	/	/
Hydrogen	1333-74-0	/	/	/	/	/	/	/	1 balloon (1 L)	ca. 1	/

Excel sheet for reagent calculation

/

Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-314: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup](#)

	Name
Used Set-up	Equipment - Advanced irradiation setup V1.0 I
Irradiation setup number	Equipment - Irradiation setup 4 (CEEC II, E002)

	Light Source Name	Power Source Name	Wavelength [nm]	Power Setting [mW]	Analog Setting [0.00 - 10.00]
First light source	Light Source - UHP LED 365 nm-4	Power Sources - BLS-18000-14	365	56	0.19

Used beam combiner [Name or None]	/
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	20
Stirring speed [rpm]	500
Irradiation start: 1. Firesting [relative to start log] 2. Unisense	1. 600 s 2. 23:10:16
Irradiation stop: 1. Firesting [relative to start log] 2. Unisense	1. 1516 s 2. 23:25:28

O₂/H₂ sensor equipment

	Equipment	Used protocol
Used Firesting	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)	Protocol - Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software
Used O ₂ sensor	Equipment - Robust probe for liquid O ₂ measurment	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor
Used H ₂ sensor	Equipment - H ₂ UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor

Procedure/observations

Date	Time	Step	Observations	Pictures/Files
20.11.2025		The experiment was done according to Protocol - In-situ hydrogen and oxygen measurment in H₂/O₂ reactor Important steps and deviations are listed below	/	/
	10:22-11:36	Conditioning of H ₂ sensor	NB-356-Logger1 1.69 mV at the end of polarization procedure	NB-356.ulog NB-356-Logger1-pre-polarization.csv NB-356-Logger1-pre-polarization.bmp

	16:00-05	Assembling the setup for calibration (25 mL of water was added using graduated cylinder), (LAUDA set to 20 °C) done according to Protocol - Liquid phase calibration of H2 UniAmp sensor with H2 bubbling.	/	/
	16:14	Start of O2 logging.	NB-356-Ch2-1	2025-11-20_161452_NB-356-Ch2-1.txt 2025-11-20_161452_NB-356-Ch2-1.png
	16:16	Start of H2 logging.	NB-356-Logger2 offset -2 mV	NB-356.ulog NB-356-Logger2-calibration.csv NB-356-Logger2-2point calibration.bmp NB-356-Logger2-calibration step.bmp
	16:18	Degassing was started.	/	20251120_161902-degassing of water.jpg
	17:12	Moving cannula to gas phase.	/	/
	17:15	Introducing H2 sensor into the reactor under Ar flow.	/	/
	17:24	0 ppm was taken.	/	/
	17:26	H2 bubbling of the reactor was started	/	20251120_172742-H2 bubbling.jpg
	17:45	1.000.000 ppm point was taken and calibration was saved	847 mV, slope: 1.076, 787 uM	20251120_174523-H2 table.jpg
	17:46	Stop of H2 logging.	/	/
	17:46	Stop of O2 logging.	/	/
	ca. 18:00	Deassembling the setup, drying the reactor with acetone and compressed air .	/	/
		Sample preparation		
	21:40	Weighing EA-358 photocatalyst in a 50 mL vial.	Creamy solid.	/
	21:44	Addition of 25 mL H2O to the vial via graduated cylinder.	/	/
	21:46-49	The suspension was vortexed for 3 min (Equipment - VWR® VV3, Vortex Mixer , stage 4/6), covered with Al foil before further use.	/	20251120_215000-suspension after vortex.jpg

		Continue in Protocol - In-situ hydrogen and oxygen measurment in H2/O2 reactor from step 6		
	22:00	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	/
	22:05	Assembling the setup.	Currently, stopper instead of H ₂ sensor, PT100, PT1000 and O ₂ robust probe are inside the reactor immersed in the liquid phase	/
	22:13	Start of O2 logging.	NB-356-Ch2-2	2025-11-20_221307_NB-356-Ch2-2.txt 2025-11-20_221307_NB-356-Ch2-2.png
	22:15	The degassing was started	/	20251120_221526-degassing of the suspension.jpg
	22:50	Cannula was transferred to gas phase, above the suspension.	/	/
	22:54	H ₂ sensor was added in Ar counterflow.	/	/
	22:58	The degassing was stopped by removing the cannula and closing the valve.	/	/
	22:59	Stop of O2 logging.	/	/
	23:00	Start of O2 logging.	NB-356-Ch2-3	2025-11-20_230014_NB-356-Ch2-3.txt 2025-11-20_230014_NB-356-Ch2-3.png
	23:00	Start of H2 logging.	NB-356-Logger3	NB-356.ulog NB-356-Logger3-during irradiation.csv NB-356-Logger3-during irradiation.bmp
	23:00-10	Equilibration time.	/	/
	23:10	The irradiation was started	/	20251120_231241-after start of irradiation.jpg
	23:25	The irradiation was stopped.	/	/
	23:25-35	Equilibration time.	/	/

	23:35	Stop of O2 and H2 logging.	/	/
	ca. 23:50	Deassembling the setup, cleaning the reactor.	Tips of the sensors and reactor were covered with attached photocatalyst particles. Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode).	20251121_000151-after irradiation.jpg

Analysis

Used calibration for Firesting: [20250611-BOLA-fitting-liquid-phase-trace-oxygen-sensor-H2-O2 reactor.ini](#)

Used calibration for UniSense: NB-356-Logger2

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
20.11.2025	10:22	NB-356-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-356.ulong NB-356-Logger1-pre-polarization.csv	/	NB-356-Logger1-pre-polarization.bmp	/	Pre-polarization of the sensor.
	16:16	NB-356-Logger2	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-356.ulong NB-356-Logger2-calibration.csv	/	NB-356-Logger2-2point calibration.bmp NB-356-Logger2-calibration step.bmp	/	Calibration of H2 sensor, signal at 10 ⁶ ppm - 847 mV, slope 1.076
	23:00	NB-356-Logger3	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-356.ulong NB-356-Logger3-during irradiation.csv	NB-356-O2 and H2 curve.py	NB-356-Logger3-during irradiation.bmp NB-356-O2 and H2 curves.png	SrTiO3 - NB-316: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C	H2 evolution during irradiation.
	16:14	NB-356-Ch2-1	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-20_161452_NB-356-Ch2-1.txt	/	2025-11-20_161452_NB-356-Ch2-1.png	/	Degassing of water followed by calibration of H2 sensor.
	22:13	NB-356-Ch2-2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-20_221307_NB-356-Ch2-2.txt	/	2025-11-20_221307_NB-356-Ch2-2.png	/	Degassing of the suspension.
	23:00	NB-356-Ch2-3	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-20_230014_NB-356-Ch2-3.txt	NB-356-O2 and H2 curve.py	2025-11-20_230014_NB-356-Ch2-3.png NB-356-O2 and H2 curves.png	SrTiO3 - NB-316: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C	O2 evolution during irradiation.

Results

Reproduction of NB-316, simultaneous H₂ and O₂ measurements of irradiated suspension of EA-358 (0.5 mg/mL) in O₂/H₂ photoreactor under 365 nm irradiation (50 mW/cm², 20 °C, 15 min, reference conditions) were performed.

Linked experiments

SrTiO₃ - NB-316: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

SrTiO₃ - NB-336: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C V (reproduction NB-316)

SrTiO₃ - NB-339: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C VI (reproduction NB-316)

SrTiO₃ - NB-348: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C VII (reproduction NB-316)

SrTiO₃ - NB-353: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C VIII (reproduction NB-316)

Linked resources

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)

Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)


Protocol - [Liquid phase calibration of H₂ UniAmp sensor](#)

Protocol - [In-situ hydrogen and oxygen measurment in H₂/O₂ reactor](#)

Attached files

NB-356-O₂ and H₂ curves.png

sha256: b88c3616b07aef9846109f298e02b33218f62bd346bc99a6c985fc8334c07792



NB-356-O2 and H2 curve.py

sha256: 16745d6e8c45ad9d3d920535e87fe621a75552d4038e7e133f7ae20c440fc034

20251120_161902-degassing of water.jpg

sha256: 5c30b62103aa95f6fa390ad8cf2b098c867a336866047354c1a08a184e634411



20251120_172742-H2 bubbling.jpg

sha256: 4e0be9c786f8d02a6b70daef59669e1fc390a64d80d2e4434c186e4758307c00




20251120_174523-H2 table.jpg

sha256: 6276ff984ac45d49fa6a2b5f82d52a914f39f0f7a83773407ec8d8e5a2cfd1bb



20251120_215000-suspension after vortex.jpg

sha256: b161cc3e07a9743e9cc95cc650feffd5a7f7022745bb96c010dfbde94aba20a3



20251120_221526-degassing of the suspension.jpg

sha256: cf49d1f044a4c4b34f17ec99981d2b70da9866c107ee7c8c5b0545b76589f17f



20251121_000151-after irradiation.jpg

sha256: 29cb13e6ae00cbf4ded79a64d2ff7c75764b818353012706502baac777d331d8



20251120_231241-after start of irradiation.jpg

sha256: f30466d9c59bb28cc2599347ac54599230443de66f7bb4b6e052111591190643




NB-356.ulog

sha256: eef826ae61bb04c2ded1b3b6af05031468a430e63d6a904e188df93de48eeec8

NB-356-Logger3-during irradiation.bmp

sha256: b0fa371ee67a3f2d227464ee289bc1397a61bce8f2b9f18b4dbf2497a2738057



NB-356-Logger3-during irradiation.csv

sha256: 598f3a411e9e4397ec10b9eed5f266af4239906ba31d20dd5197e11cf6671634

NB-356-Logger2-calibration.csv

sha256: c9d13f88a6a962cf50df3ea8b3b9925fcd02f5b3e40e79ba6ecd49a6c35ea6d

NB-356-Logger2-2point calibration.bmp

sha256: c2a6213cc9f96f589381156572ea0e820ac3bfe1d86db633c47b10cb38490cee



NB-356-Logger2-calibration step.bmp

sha256: ce93bb240f3ce0557e08fc495758ff577b8a82e6298fb6626750bb82f3177e40




NB-356-Logger1-pre-polarization.csv

sha256: fae44eb57cbbb6b3de31b69444c6d21738b91a771a37e0cc73e96cf461020d9a

NB-356-Logger1-pre-polarization.bmp

sha256: 8b2a9d1b2ded6ffec6a9f467fafebba2cb5a10eac8a37155f2db1de1ca93a7a




2025-11-20_230014_NB-356-Ch2-3.txt

sha256: 497dc363c0835a9fedb0538c9c2228bf20b68a00fc8945cd009d585496cc3ab3

2025-11-20_230014_NB-356-Ch2-3.png

sha256: b6407b60486378b06e2ecb9ca15ed032bf022c8e753e3b4aeccc06f9506bd415



2025-11-20_221307_NB-356-Ch2-2.txt

sha256: d79f20c4a4459377c86081e55f8779b2b5b5d3e6f2a1fcb9a25f4869828b8242


2025-11-20_221307_NB-356-Ch2-2.png

sha256: cdad2b6dc3b549ed279c0bacaf20d958fc2627aae68c57ae3e146063b29c3555



2025-11-20_161452_NB-356-Ch2-1.png

sha256: 12fa315862bce72fc11e15751ed66db43820bfcd72ba35346a50b766e800bcb9



2025-11-20_161452_NB-356-Ch2-1.txt

sha256: 1d68cfa16b7aba14d2bed4672b60245fc7c9bfb5dbff78a4466f64c2e990736e



Unique eLabID: 20251120-11d0950e68e9a0affdd96447c2302bb183c9769

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3516>