

NB-325: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm², 20 °C

Date: 2025-10-28

Tags: O₂ Test Calibration NB Firing
O₂ sensor H₂ SrTiO₃ Unisense
RhCrO₃:Al:SrTiO₃ H₂ Sensor
temperature In situ Trace range robust
oxygen sensor photocatalysis Unisense
normal range

Category: SrTiO₃


Status: Done

Created by: Nadzeya Brezhneva

Objectives

Simultaneous detection of H₂ and O₂ evolution in liquid phase for irradiated suspension of Rh,CrO_x:Al:SrTiO₃ suspension (EA-358, 0.5 mg/mL), 365 nm LED, **20 mW/cm²**, 20 °C.

Reaction scheme



ChemDraw file linked: [NB-325-SrTiO3-photocatalytic H2O splitting.cdxml](#)

Literature/reference experiments

Literature	/
Reproduction	/
Similar experiments	SrTiO3 - NB-316: Liquid phase H2 and O2 of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	Concentration [mM]
milli-Q H ₂ O	/	/	/	/	/	/	/	0.998	25	/
Al:SrTiO ₃ RhCrO _x (EA-358)	SrTiO3 - EA-358: Modification of Al:SrTiO₃ (EA-354) via deposition of Rh, Cr oxide co-catalyst, 350°C, 1h, Upscaling (3.33x)	/	/	/	12.50	12.54	/	/	/	/

Excel sheet for reagent calculation

/

Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-324: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup I I](#)

	Name
Used Set-up	Equipment - Advanced irradiation setup V1.0 I
Irradiation setup number	Equipment - Irradiation setup 4 (CEEC II, E002)

	Light Source Name	Power Source Name	Wavelength [nm]	Power Setting [mW]	Analog Setting [0.00 - 10.00]
First light source	Light Source - UHP LED 365 nm-4	Power Sources - BLS-18000-1 4	365	22	0.07

Used beam combiner [Name or None]	/
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	20
Stirring speed [rpm]	500
Irradiation start: 1. Firesting [relative to start log] 2. Unisense	1. 601 s 2. 2:13:54
Irradiation stop: 1. Firesting [relative to start log] 2. Unisense	1. 1531 s 2. 2:29:25

O₂/H₂ sensor equipment

	Equipment	Used protocol
Used Firesting	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)	Protocol - Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software
Used O ₂ sensor	Equipment - Robust probe for liquid O ₂ measurment	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor
Used H ₂ sensor	Equipment - H ₂ UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor

Procedure/observations

Date	Time	Step	Observations	Pictures/Files
		Calibration was used from experiment SrTiO ₃ - NB-320: Liquid phase H ₂ and O ₂ of RhCrOx,Al:SrTiO ₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm ² , 20 °C (reproduction NB-316) I		
28.10.2025		Sample preparation		
	ca. 0:50	Weighing EA-358 photocatalyst in a 50 mL vial.	Creamy solid.	/
	0:55	Addition of 25 mL H ₂ O to the vial via graduated cylinder.	/	/
	1:00	The suspension was vortexed for 3 min (Equipment - VWR® VV3, Vortex Mixer, stage 4/6), covered with Al foil before further use.	/	20251028_010451-suspension after vortex.jpg
		Continue in Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor from step 6		
	1:15	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	/
	1:20	Assembling the setup.	Currently, stopper instead of H ₂ sensor, PT100, PT1000 and O ₂ robust probe are inside the reactor immersed in the liquid phase	/

	1:22	Start of O2 logging.	NB-325-Ch2-1	2025-10-28_012227_NB-325-Ch2-1.txt 2025-10-28_012227_NB-325-Ch2-1.png
	1:24	The degassing was started	/	/
	1:56	Cannula was transferred to gas phase, above the suspension.	/	/
	1:59	H ₂ sensor was added in Ar counterflow.	/	/
	ca. 2:02	The degassing was stopped by removing the cannula and closing the valve.	/	/
	2:03	Stop of O2 logging.	/	20251028_020302-before irradiation.jpg
	2:03	Start of O2 logging.	NB-325-Ch2-2	2025-10-28_020353_NB-325-Ch2-2.txt 2025-10-28_020353_NB-325-Ch2-2.png
	2:05	Start of H2 logging.	NB-325-Logger1	NB-325.ulong NB-325-Logger1.csv NB-325-Logger1.bmp
	2:03-13	Equilibration time.	/	/
	2:13	The irradiation was started	/	20251028_021744-after start of irradiation.jpg
	2:29	The irradiation was stopped.	/	/
	2:29-39	Equilibration time.	/	/
	2:39	Stop of O2 and H2 logging.	Bubbles near the tips of the sensors	20251028_024050-after irradiation.jpg

	2:40-3:00	Deassembling the setup, cleaning the reactor.	Tips of the sensors and reactor were covered with attached photocatalyst particles. Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode).	/
--	-----------	---	--	---

Analysis

Used calibration for Firesting: [20250611-BOLA-fitting-liquid-phase-trace-oxygen-sensor-H2-O2 reactor.ini](#)

Used calibration for UniSense: NB-320-Logger6

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
28.10.2025	2:05	NB-325-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-325.ulong NB-325-Logger1.csv	NB-325-O2 and H2 curve.py	NB-325-Logger1.bmp NB-325-O2 and H2 curve.png	/	Clean response of H2 sensor, ca. 18.6 μM H2 at the end of irradiation.
	1:22	NB-325-Ch2-1	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-10-28_012227_NB-325-Ch2-1.txt	/	2025-10-28_012227_NB-325-Ch2-1.png	/	Degassing.
	2:03	NB-325-Ch2-2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-10-28_020353_NB-325-Ch2-2.txt	NB-325-O2 and H2 curve.py	2025-10-28_020353_NB-325-Ch2-2.png NB-325-O2 and H2 curve.png	/	Change in the slope of the O2 curve during irradiation, ca. 17.6 μM O2 at the end of irradiation.

Results

Simultaneous H₂ and O₂ measurements of irradiated suspension of EA-358 (0.5 mg/mL) in O₂/H₂ photoreactor under 365 nm irradiation (20 mW/cm², 20 °C, 15 min) were performed.

H₂ level at the end of irradiation - ca. 18.6 $\mu\text{mol/L}$, O₂ level - ca. 17.6 $\mu\text{mol/L}$.

Linked experiments

SrTiO₃ - NB-316: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

SrTiO₃ - NB-318: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm², 20 °C

SrTiO₃ - NB-320: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C (reproduction NB-316) I

SrTiO₃ - NB-322: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm², 20 °C (reproduction NB-318)

Linked resources

Equipment - VWR® VV3, Vortex Mixer

Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)

Equipment - Irradiation setup 4 (CEEC II, E002)

Protocol - Getting hydrogen from hydrogen bottle in CEEC II E014

Protocol - In-situ hydrogen and oxygen measurement in H₂/O₂ reactor

Attached files

NB-325-SrTiO₃-photocatalytic H₂O splitting.png

sha256: 84aba7b2eaa7e7c6c75ec4bc4edf260ee29035205566e821b254e3f00db612f8



NB-325-SrTiO₃-photocatalytic H₂O splitting.cdxml


sha256: 1e25f76a703054a4e3a1c95ca0b8c004e2e67c7a646758e715ac13662fa83a8c

NB-325-O₂ and H₂ curve.py

sha256: 4f3a701ca60f62fa6e4a9fd575d73ec5a989c2b1ece80795ed98ec303afe4df3

NB-325-O₂ and H₂ curve.png

sha256: ed31187f0af52b5b5a00247b655ec115b878f0a68ecf7a1d18d7759383c6cace



20251028_010451-suspension after vortex.jpg

sha256: e81f497fd4b555e5569291d23301ed1cc00547e7cfb02744625eefd7d63b5e78



20251028_021744-after start of irradiation.jpg

sha256: 0a9e297dac7cb833b4acff6770bb0bf97657aa36afa5a10483376daac217ac86



20251028_024050-after irradiation.jpg

sha256: ebe222b7176747583fd388813bbcf22f07178a282ae6dc2cbca126c20c960b7f




20251028_020302-before irradiation.jpg

sha256: b7eba86e19f39c1e7a1e46f1bfdecd2f2a7ff1ee979df270c41931e039e93dd1



NB-325-Logger1.bmp

sha256: 6f75b99e23ecf23e9348ae5c02ac6bf691fb9987e63b16039f0149eee0aac303




NB-325-Logger1.csv

sha256: 2a29a16e30cfafcdb99c199bf40deb976470b6bae6dad54a78ed749ee1f20b6d

2025-10-28_020353_NB-325-Ch2-2.png

sha256: 2e788ee24b678c7e274c2234861778735f27dc4cdcb6288e43d35307e84a099c



2025-10-28_020353_NB-325-Ch2-2.txt

sha256: b1e7319ad03b46672390ab28a3e3bdf278deae539a8e31b9ab9a9d42610d70ac

2025-10-28_012227_NB-325-Ch2-1.png

sha256: ef8c0f84f2ce9a61549246734de3ef0bf42af94d2e9321a303e6f5f86a6c6784



2025-10-28_012227_NB-325-Ch2-1.txt

sha256: 159305bfca94dd583954aaf399d97730d8cf9f55b574c036da2c72b6e22764b7

NB-325.ulog

sha256: 46defe82218f7652a217edac8e2dd2c5b862a639f3e543c8f71ccd0beee8f18a



Unique eLabID: 20251028-9fb5d102ee8a924df408cabb9a2b0363708dbc6a

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3267>