

NB-312: Gas phase H2 and O2 measurements with Unisense H2 sensor, Firesting O2 robust probe in irradiated Al:SrTiO3 RhCrOx (NB-289, 0.5 mg/mL), 365 nm, 50 mW, 1 h, degassing

Date: 2025-10-08

Tags: O2 Test Calibration NB Firesting
Irradiation O2 sensor H2 advanced irrada
setup Unisense H2 Sensor temperature
In situ Trace range robust oxygen sensor
photocatalysis

Category: SrTiO3


Status: Done

Created by: Nadzeya Brezhneva

Objectives

Simultaneous gas phase H2 and O2 measurements in irradiated Al:SrTiO3 RhCrOx (NB-289, 0.5 mg/mL), 365 nm, 50 mW, 1 h, degassing.

Reaction scheme



ChemDraw file linked: [NB-304-SrTiO3-photocatalytic H2O splitting.cdxml](#)

Literature/reference experiments

Literature	https://doi.org/10.1039/C9EE00310J
Reproduction	SrTiO3 - NB-304: Gas phase H2 and O2 measurements with Unisense H2 sensor, Firesting O2 robust probe in irradiated Al:SrTiO3 RhCrOx (NB-289, 0.5 mg/mL), 365 nm, 50 mW, 1 h, degassing
Similar experiments	SrTiO3 - NB-231: Gas phase H2 and O2 measurements with Unisense H2 sensor, Firesting O2 cap in irradiated Al:SrTiO3 RhCrOx (NB-162-4, 1 mg/mL), 365 nm, 50 mW, 1 h, degassing

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	pressure [bar]
milli-Q H ₂ O	/	/	/	/	/	/	/	/	25	/
Hydrogen	1333-74-0	/	/	/	/	/	/	/	1 ballon (approx. 2 L)	approx. 1

Al:SrTiO ₃ RhCrO _x (NB-289)	SrTiO ₃ - NB-289: Modification of NB-285 (SrTiO ₃ :Al (from self-made SrTiO ₃ , Osterloh, no Al ₂ O ₃ , PVDF filter) 1000 C, 10 h) with Rh, Cr oxide cocatalyst	/	/	/	12.50	12.62	/	/	/	/
---	--	---	---	---	-------	-------	---	---	---	---

Irradiation Parameters

Power measurement was performed in experiment [Prep work - NB-310: Measuring power output of UHP-365 nm #1 with 18A-2 in advanced irradiation setup](#)

	Name
Used set-up	Equipment - Advanced irradiation setup V1.0 I
Irradiation setup number	Equipment - Irradiation setup 4 (CEEC II, E002)

	Light Source Name	Used power source	Wavelength [nm]	Power Setting [mW]	Analogue settings [-]
First light source	Light Source - UHP LED 365 nm-1	Power Sources - BLS-18000-1 2	365	50	0.32

Used beam combiner [Name or None]	None
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	25
Stirring speed [1/min]	500
Start time [s], relative to start of log	930 s - Firesting 21:49:41 (not relative time) - Unisense

End time [s], relative to start of log	4650 s - Firing 22:51:41 - Unisense
---	---

O₂/H₂ sensor equipment

	Equipment	Used protocol
Used Firing	Equipment - Firing Fiber-Optic Oxygen Meter 2 Channel (Firing 2)	Protocol - Operation of Firing Fiber-Optic Oxygen Meter 2 Channel Software
Used O ₂ sensor	Equipment - Robust probe for liquid O ₂ measurement	Protocol - In-situ hydrogen and oxygen measurement in H ₂ /O ₂ reactor
Used H ₂ sensor	Equipment - H ₂ UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurement in H ₂ /O ₂ reactor

Procedure/observations

Date	Time	Step	Observations	Pictures/Files
09.10.2025		The experiment was done according to Protocol - In-situ hydrogen and oxygen measurement in H₂/O₂ reactor Important steps and deviations are listed below	/	/
		The PT100 was not placed inside the reactor, but outside the reactor inside the advanced irradiation setup, close to the double walled beaker.	/	/
	15:44-16:43	Conditioning of H₂ sensor	NB-312-Logger-1 1 s interval, 1.39 mV at the end of pre-polarization	NB-312-pre-measurements.bmp
	ca. 17:50	Assembling the set-up for calibration.	/	/
	ca. 18:00	Addition of 25 mL milli-Q H ₂ O to the photoreactor via funnel.	/	/
	ca. 18:10	Assembling sensors with appropriate BOLA fitting and GL14/NS14 adapters.	PT1000 - 4 mm BOLA fitting Robust oxygen probe - 3 mm BOLA fitting Unisense sensor - 2 mm fitting (used from the previous measurement, only removed from the protective holder)	/

		Calibration was done first according to Protocol - Liquid phase calibration of H2 UniAmp sensor with H2 bubbling	/	/
	18:19	The O2 log was started	NB-312-Ch2-1	2025-10-09_181947_NB-312-Ch2-1.txt 2025-10-09_181947_NB-312-Ch2-1.png
	18:19	The H2 log was started	NB-312-Logger-2 offset -1 mV	NB-312-calibration.bmp NB-312-calibration-1.bmp
	18:22	0 ppm point taken	0.51 mV	/
	18:23	Degassing was started by immersing cannula (from Ar supply setup) in water.	/	20251009_182518-calibration-during degassing.jpg
	18:41	Removing cannula from H ₂ O closing the valve.	/	/
	18:42	H2 bubbling of the reactor was started.	/	/
	18:50	1.000.000 ppm point was taken and calibration was saved	22.76 °C 1004.5 mbar 1.000.000 ppm 97674.034 Pa Signal of the sensor 1017 mV Slope: 0.010, intercept: 0.491	/
	18:51	The H2 logging was stopped.	/	/
	18:51	The O2 logging was stopped.	/	/
	18:55-20:00	Deassembling the setup. degreasing and drying the reactor at 120 °C.	/	/
		Sample preparation		
	20:10	Weighing 12.57 mg of NB-289 sample in a 50 mL vial.	Creamy solid.	20251009_201235-weighed photocatalyst.jpg
	20:17	Addition of 25 mL H2O to the vial via graduated cylinder.	/	/

	20:19-22	The suspension was vortexed for 3 min (Equipment - VWR® VV3, Vortex Mixer , stage 4/6), covered with Al foil before further use.	Milky white suspension.	20251009_202215-suspension after vortex.jpg
		Continue in Protocol - In-situ hydrogen and oxygen measurment in H2/O2 reactor from step 6	/	/
	ca. 20:40-50	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	20251009_205632-before irradiation.jpg 20251009_205642-positions of the sensors inside the reactor.jpg
	21:00	The O2 log was started	NB-312-Ch2-2	2025-10-09_210005_NB-312-Ch2-2.txt 2025-10-09_210005_NB-312-Ch2-2.png
	21:11	The degassing was started	/	/
	21:31	The degassing was stopped	No leakage	/
	21:33	The O2 log was stopped	Adjusting stirring, there were some problems with stirring of the suspension, but were eliminated.	/
	21:34	The O2 log was started	NB-312-Ch2-3	2025-10-09_213411_NB-312-Ch2-3.txt 2025-10-09_213411_NB-312-Ch2-3.png
	21:34	The H2 log was started	NB-312-Logger-3	NB-312-during irradiation-H2 evolution.bmp
	21:35-49	Equilibration time.	21:35-38 - making sure that the stirring was good	/
	21:49	The irradiation was started	0 h : 15 min : 30 s - Firesting 21:49:41 - Unisense (there is no lid on the setup, covered only by the side wall)	20251009_215035-after start of irradiation.jpg
	22:51	The irradiation was stopped	1 h : 17 min : 30 s - Firesting 22:51:41 - Unisense 0.75 vol.% O ₂ at the end of irradiation 1558 Pa H ₂ at the end of irradiation	/

	23:06	The H2 on O2 log were stopped	/	/
	23:10	Deassembling the setup, cleaning the reactor.	Adhesion of photocatalyst to the walls of the reactor.	20251009_230913-after irradiation.jpg

Analysis

Used calibration for Firesting: [Prep work - NB-280: Gas phase calibration of O2 robust oxygen sensor \(Firesting\) in a 4 neck H2/O2 photoreactor](#)

Used calibration for Unisense: NB-312-Logger2

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
09.10.2025	15:44	NB-312-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor	H2O	NB-313-Logger.ulong NB-312-Logger1-Data (1 H2 UniAmp (503157)).csv	/	NB-312-pre-measurements.bmp	/	Initial slight decrease in the signal, afterwards stabilizing (1.39 mV at the end of pre-polarization).
	18:19	NB-312-Logger2	electrochemical H2 detection	Equipment - H2 UniAmp Sensor	H2O	NB-313-Logger.ulong NB-312-Logger2-Data (1 H2 UniAmp (503157)).csv	/	NB-312-calibration.bmp NB-312-calibration-1.bmp	SrTiO3 - NB-304: Gas phase H2 and O2 measurements with Unisense H2 sensor, Firesting O2 robust probe in irradiated Al:SrTiO3 RhCrOx (NB-289, 0.5 mg/mL, 365 nm, 50 mW, 1 h, degassing	two point calibration (0 and 1.000.000 ppm), Slope = 0.010, intercept = 0.451, R ² = 1
	21:34	NB-312-Logger3	electrochemical H2 detection	Equipment - H2 UniAmp Sensor	H2O	NB-313-Logger.ulong NB-312-Logger3-Data (1 H2 UniAmp (503157)).csv	NB-312-with O2 and H2 axes.py NB-312-with O2, H2 and T axes.py	NB-312-during irradiation-H2 evolution.bmp NB-312-H2 and O2.jpeg NB-312-H2, O2 and T.jpeg	/	Gradual increase in H2 value (Up to ca. 1558 Pa at the end of irradiation).
	18:19	NB-312-Ch2-1	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)	H2O	2025-10-09_181947_NB-312-Ch2-1.txt	/	2025-10-09_181947_NB-312-Ch2-1.png	/	Degassing (during calibration) worked, till 0.00 vol.%
	21:00	NB-312-Ch2-2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)	H2O	2025-10-09_210005_NB-312-Ch2-2.txt	/	2025-10-09_210005_NB-312-Ch2-2.png	/	Degassing worked, no leakage after finishing degassing step in the configuration of the assembled setup for irradiation.
	21:34	NB-312-Ch2-3	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)	H2O	2025-10-09_213411_NB-312-Ch2-3.txt	NB-312-with O2 and H2 axes.py NB-312-with O2, H2 and T axes.py	2025-10-09_213411_NB-312-Ch2-3.png NB-312-H2 and O2.jpeg NB-312-H2, O2 and T.jpeg	/	Increase in O2 value during irradiation (up to 0.75 vol. % after 1 h irradiation).

Results

Gas phase measurements of O₂ and H₂ were performed during photocatalytic test of Al:SrTiO₃ suspension (reproduction of NB-304, sample **NB-289**, 0.5 mg/L, 365 nm, 50 mW, 1h) under degassed conditions were performed. Increase in O₂ and H₂ values during irradiation was observed (up to 0.75 vol.% O₂ and 1558 Pa

H₂ at the end of irradiation.

for comparison: [SrTiO₃ - NB-304: Gas phase H₂ and O₂ measurements with Unisense H₂ sensor, Firesting O₂ robust probe in irradiated Al:SrTiO₃ RhCrOx \(NB-289, 0.5 mg/mL\), 365 nm, 50 mW, 1 h, degassing](#)

O₂: 0.71 vol.% at the end of irradiation,

H₂: 1370 Pa at the end of irradiation.

Linked experiments

[SrTiO₃ - NB-289: Modification of NB-285 \(SrTiO₃:Al \(from self-made SrTiO₃, Osterloh, no Al₂O₃, PVDF filter\) 1000 C, 10 h\) with Rh, Cr oxide cocatalyst](#)

[SrTiO₃ - NB-304: Gas phase H₂ and O₂ measurements with Unisense H₂ sensor, Firesting O₂ robust probe in irradiated Al:SrTiO₃ RhCrOx \(NB-289, 0.5 mg/mL\), 365 nm, 50 mW, 1 h, degassing](#)

Linked resources

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Robust probe for liquid O₂ measurment](#)

Equipment - [Advanced irradiation chamber V1.0 I](#)

Equipment - [H₂ UniAmp Sensor - Normal range - 2.1 x 80 mm needle](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)

Protocol - [Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software](#)

Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)


Protocol - [Gas phase calibration of H₂ UniAmp sensor](#)

Protocol - [In-situ hydrogen and oxygen measurment in H₂/O₂ reactor](#)

Attached files


NB-312-H₂, O₂ and T.jpeg

sha256: 878b306abfef70f0e830b93c9db5e3c99282faa2affa0c5b4bd8e6ed96b9eb41



NB-312-H2 and O2.jpeg

sha256: 0b2cfcead1f03f08d19fee2ebe972f9b5784d109634b271391b83203da7814de



NB-312-with O2 and H2 axes.py

sha256: c384ddf9dcaf3fb5f8a28db102432bb7790d6ca97b19f6f8a2a1cf2285c4a854

NB-312-with O2, H2 and T axes.py

sha256: 75883973adb3c7ee0dcc03dcf936fe9c7051078aa7e7ddb0020e639e9280a8e7

20251009_182518-calibration-during degassing.jpg

sha256: 227b7d514987c4e224465f2bf6f8215d190ebc51c5b984fbc5977cbce72257a1



20251009_184335-H2 bubbling.jpg

sha256: b5c26179c001acb2b7abe8b27c122e83c2bcb0e210c09d6cc43c7d93cc90cafd



20251009_201235-weighed photocatalyst.jpg

sha256: e0b0495f4aae2f80f981a9acf2aac40991231ac1df2891c792339283b78cb212



20251009_202215-suspension after vortex.jpg

sha256: 89764e8845d0fce6d93c787211cf67f8bfa748a327e41b529761c730ae1b21f9



20251009_205632-before irradiation.jpg

sha256: a46630bde01b34339c342f748dbb6e3419e091ae9a455fad1f05cb9d8aeddec9



20251009_211211-degassing of the suspension.jpg

sha256: 7c87650bc8481e33b3502a4ff5a2f011d1fea6f8c4247b3d3ef0f024ca0bba28



20251009_205642-positions of the sensors inside the reactor.jpg

sha256: 4b597ca6f5d81938bffee3374a2478d6b3157b9d5bb02a5b773e84e849fe6952



20251009_230913-after irradiation.jpg

sha256: 0e63e341818f93ccaafc3b539a749f4bf4f2060e59b0379eeee36b26734a8c8e3



20251009_215035-after start of irradiation.jpg

sha256: 7e685c09831cf77c8f39bdaacef4f2c374ee04261142f5193c60d2f3703d9d57




2025-10-09_213411_NB-312-Ch2-3.txt

sha256: b51c8c772d5be3c02c77b9419b7685101620db170f0698eb9e003b26cc691d4d

2025-10-09_213411_NB-312-Ch2-3.png

sha256: 531b00bd8f61e7f417aaf37afbd7d7088601e69677c3c0479cb670c00a3a594c




2025-10-09_210005_NB-312-Ch2-2.txt

sha256: 3962b6bd09347347c36993edb8d274917d0719337ed3d0494c38851da16406db


2025-10-09_210005_NB-312-Ch2-2.png

sha256: f5a6e956362e6cbddc8bff159a1462005c8dc9f2f8bdb7313acda9fa3f1e8706



2025-10-09_181947_NB-312-Ch2-1.png

sha256: a1714955b0618df1c3cec67121f57a77cbefb0c6602c2e91c087e5611722ee61



2025-10-09_181947_NB-312-Ch2-1.txt

sha256: ed2b805f227d970a3576e1c8adec8b63d4448f000def3ec37e113a9f21887157

NB-312-Logger3-Data (1 H2 UniAmp (503157)).csv

sha256: 6e3fd08276d07720b526b605b3e0bb8033812ebb8e3ac6391831607e7cbd442f

NB-312-Logger2-Data (1 H2 UniAmp (503157)).csv


sha256: 0065b9235ee41f6d5be787eab2d76a34a020602685320d717c3e8dc4233d8c0f

NB-312-Logger1-Data (1 H2 UniAmp (503157)).csv

sha256: c00be40c86e04ad5a03e5ba876c103139c29162e5a3096f9c8d35f6a3362e793

NB-312-calibration.bmp

sha256: 46e554b6bd93066ed697b976bef3a0bd5496a059812c9e30f44e099ddfde6800




NB-313-Logger.ulog

sha256: d695475fe30cb64653bdc2d7f51dd6a815243fef4f5debfda53901598b745e0b


NB-312-pre-measurements.bmp

sha256: 713489ccaae91f3438c56020a89e05b9289353e4032de21e858bc3da9efe2eaf




NB-312-during irradiation-H2 evolution.bmp

sha256: 6fd914118b4c4c2767a63645924347e0879289028b2ff809b11aa8e5913a6805



NB-312-calibration-1.bmp

sha256: b40cc1f766326c17328e73b19d006bcaa80b436a1201aabaaf820219d49dcf84



Unique eLabID: 20251008-23b4af4a73c64cb2b5765ff58dacf260e68aa366

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3129>