

NB-363: Gas phase H2 and O2 measurements of Al:SrTiO3 RhCrOx (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C, 70 min, degassing (reproduction NB-362)

Date: 2025-11-26

Tags: O2 Test Calibration NB Firesting
Irradiation O2 sensor H2 advanced irrada
setup Unisense H2 Sensor temperature
In situ Trace range robust oxygen sensor
photocatalysis

Category: SrTiO3


Status: Done

Created by: Nadzeya Brezhneva

Objectives

Reproduction of NB-362: simultaneous detection of O₂ and H₂ evolution in gas phase for irradiated suspension of Rh_xCrO_y:Al:SrTiO₃ suspension (EA-358, 0.5 mg/mL), 365 nm LED, 50 mW/cm², 20 °C during 70 min.

Reaction scheme



ChemDraw file linked: [NB-362-SrTiO3-photocatalytic H2O splitting.cdxml](#)

Literature/reference experiments

Literature	/
Reproduction	SrTiO3 - NB-362: Gas phase H2 and O2 measurements of Al:SrTiO3 RhCrOx (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C, 70 min, degassing
Similar experiments	SrTiO3 - NB-361: Gas phase H2 and O2 measurements of Al:SrTiO3 RhCrOx (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C, 15 min, degassing

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	pressure [bar]
milli-Q H ₂ O	/	/	/	/	/	/	/	/	25	/

Al:SrTiO3 RhCrOx (EA-358)	SrTiO3 - EA-358: Modification of Al:SrTiO3 (EA-354) via deposition of Rh, Cr oxide co- catalyst, 350°C, 1h, Upscaling (3.33x)	/	/	/	12.50	12.47	/	/	/	/
---------------------------------	---	---	---	---	-------	-------	---	---	---	---

Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-314: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup](#)

	Name
Used Set-up	Equipment - Advanced irradiation setup V1.0 I
Irradiation setup number	Equipment - Irradiation setup 4 (CEEC II, E002)

	Light Source Name	Power Source Name	Wavelength [nm]	Power Setting [mW]	Analog Setting [0.00 - 10.00]
First light source	Light Source - UHP LED 365 nm-4	Power Sources - BLS-18000-14	365	56	0.19

Used beam combiner [Name or None]	/
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	20
Stirring speed [rpm]	500
Irradiation start: 1. Firesting [relative to start log] 2. Unisense	1. 610 s 2. 23:30:27

Irradiation stop: 1. Firesting [relative to start log] 2. Unisense	1. 4855 s 2. 0:41:12
---	-------------------------

O₂/H₂ sensor equipment

	Equipment	Used protocol
Used Firesting	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)	Protocol - Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software
Used O ₂ sensor	Equipment - Robust probe for liquid O ₂ measurment	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor
Used H ₂ sensor	Equipment - H ₂ UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor

Procedure/observations

Date	Time	Step	Observations	Pictures/Files
25.11.2025		Calibration from experiment SrTiO₃ - NB-361 : Gas phase H ₂ and O ₂ measurements of Al:SrTiO₃ RhCrOx (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm ² , 20 °C, 15 min, degassing was used.		
		Sample preparation		
	22:20	Weighing EA-358 photocatalyst in a 50 mL vial.	Creamy solid.	/
	22:24	Addition of 25 mL H ₂ O to the vial via graduated cylinder.	/	/
	22:25-28	The suspension was vortexed for 3 min (Equipment - VWR® VV3, Vortex Mixer , stage 4/6), covered with Al foil before further use.	/	/
		Continue in Protocol - In-situ hydrogen and oxygen measurment in H₂/O₂ reactor from step 6		
	22:30	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	20251125_222810-suspension after vortex.jpg
	22:35	Assembling the setup.	/	/

	22:45	Start of O2 logging.	NB-363-Ch2-1	2025-11-25_224558_NB-363-Ch2-1.txt 2025-11-25_224558_NB-363-Ch2-1.png
	22:48	The degassing was started	/	20251125_225034-degassing of the suspension.jpg
	23:20	The degassing was stopped by removing the cannula and closing the valve.	/	/
	23:19	Stop of O2 logging.	/	/
	23:20	Start of O2 logging.	NB-363-Ch2-2	2025-11-25_232017_NB-363-Ch2-2.txt 2025-11-25_232017_NB-363-Ch2-2.png
	23:20	Start of H2 logging.	NB-363-Logger1	NB-363.ulog NB-363-Logger1.csv NB-363-Logger1.bmp
	23:20-30	Equilibration time.	/	/
	23:30	The irradiation was started	After 20:10 the O ₂ and H ₂ curves reached plateau	20251125_233106-after start of irradiation.jpg
26.11.2025	0:40	The irradiation was stopped.	/	/
	0:40-51	Equilibration time.	/	/
	0:51	Stop of O2 and H2 logging.	/	/
	ca. 1:00	Deassembling the setup, cleaning the reactor.	Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode).	/

Analysis

Used calibration for Firesting: [20250910-BOLA fitting-gas phase-4-neck photoreactor-trace oxygen robust probe-Ch2.ini](#)

Used calibration for UniSense: NB-361-Logger2

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
25.11.2025	23:20	NB-363-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-363.ulong NB-363-Logger1.csv	NB-363-O2 and H2 curve.py	NB-363-Logger1.bmp NB-363-O2 and H2 curves.png	SrTiO3 - NB-362: Gas phase H2 and O2 measurements of Al:SrTiO3 RhCrOx (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C, 70 min, degassing	H2 evolution during irradiation
	22:45	NB-363-Ch2-1	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-25_224558_NB-363-Ch2-1.txt	/	2025-11-25_224558_NB-363-Ch2-1.png	/	Degassing of the suspension.
	23:20	NB-363-Ch2-2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-25_232017_NB-363-Ch2-2.txt	NB-363-O2 and H2 curve.py	2025-11-25_232017_NB-363-Ch2-2.png NB-363-O2 and H2 curves.png	SrTiO3 - NB-362: Gas phase H2 and O2 measurements of Al:SrTiO3 RhCrOx (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C, 70 min, degassing	O2 evolution during irradiation

Results

Reproduction of NB-362: simultaneous H₂ and O₂ measurements (gas phase) of irradiated suspension of EA-358 (0.5 mg/mL) in O₂/H₂ photoreactor under 365 nm irradiation (50 mW/cm², 20 °C, 70 min) were performed.

Linked experiments

SrTiO3 - NB-312: Gas phase H2 and O2 measurements with Unisense H2 sensor, Firesting O2 robust probe in irradiated Al:SrTiO3 RhCrOx (NB-289, 0.5 mg/mL), 365 nm, 50 mW, 1 h, degassing (reproduction NB-304)

SrTiO3 - NB-361: Gas phase H2 and O2 measurements of Al:SrTiO3 RhCrOx (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C, 15 min, degassing

SrTiO3 - NB-362: Gas phase H2 and O2 measurements of Al:SrTiO3 RhCrOx (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C, 70 min, degassing

Linked resources

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Robust probe for liquid O2 measurment](#)

Equipment - [Advanced irradiation chamber V1.0 I](#)

Equipment - [H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)

Protocol - [Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software](#)

Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)

Protocol - [Gas phase calibration of H2 UniAmp sensor](#)

Protocol - [In-situ hydrogen and oxygen measurment in H2/O2 reactor](#)


Attached files

NB-363-O2 and H2 curve.py

sha256: c1727f05410721a2361f9b3f784f209f9c1ac8ae1508eac38f31c82d1e96e6ce


NB-363-O2 and H2 curves.png

sha256: 8de361b71694fc219ce174b8604df6e7c266a4c53b720024428c88ac6f95dc8d



2025-11-25_232017_NB-363-Ch2-2.png

sha256: b1340b9897b2ceee77f20eded659041c3441f63a6d499da77ffb9e16a45c2d98



2025-11-25_232017_NB-363-Ch2-2.txt


sha256: 6f84838931b34117d614bffae16edf24f361f045af4fda970a26aca956101759

2025-11-25_224558_NB-363-Ch2-1.txt

sha256: eb0096c43747ced63c18ff87da997ead5c9e922f7ef3d96a6c4daead3a1b52c9


2025-11-25_224558_NB-363-Ch2-1.png

sha256: 087e7aed4ed0cf9d8bb181f7b7741fff317d01922e97c23fad73354e43dc1131



NB-363-Logger1.bmp

sha256: 9519996e77e7ec5925a2e6ed0fd3e6ed43e992e4a46adb97bdb2af3e2af6a74b



NB-363.ulog

sha256: 11f622e8f2f453acb7a11eeef4e7375a2fb149dea99b89aa0d06457036d2b28

NB-363-Logger1.csv

sha256: e19f768f14bad54e78197816084e8e8fbb799f341afd151dd0ea6dece2f7e930

20251125_225034-degassing of the suspension.jpg

sha256: d836ff4720f37b537bbf8915266a30b9dcf98ecfaa56f979ce67283a566634b9



20251125_222810-suspension after vortex.jpg

sha256: 9733cb802de8c2f4e455d2d1e912ddebada79784c05b56156164a193fb61d9540



20251125_233106-after start of irradiation.jpg

sha256: 94307d9356d73f68b81c740efded67bf9f692cf09c49ae0be84737f4dbbb326b



20251126_005240-after irradiation.jpg

sha256: 440ac83918b1156b71254c566be45e0507e8f9664be10b2c276f4f6de0a9dcde



Unique eLabID: 20251126-4def0d3764a01a0bd94090ae58d63a360acae764
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3558>