

# NB-351: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrO<sub>x</sub>,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 10 °C

Date: 2025-11-18

Tags: O<sub>2</sub> Test Calibration NB Firing  
O<sub>2</sub> sensor H<sub>2</sub> SrTiO<sub>3</sub> Unisense  
RhCrO<sub>3</sub>:Al:SrTiO<sub>3</sub> H<sub>2</sub> Sensor  
temperature In situ H<sub>2</sub>/O<sub>2</sub> reactor Trace  
range robust oxygen sensor  
photocatalysis Unisense normal range

Category: SrTiO<sub>3</sub>


Status: Done

Created by: Nadzeya Brezhneva

## Objectives

Simultaneous detection of H<sub>2</sub> and O<sub>2</sub> evolution in liquid phase for irradiated suspension of Rh<sub>x</sub>CrO<sub>y</sub>:Al:SrTiO<sub>3</sub> suspension (EA-358, 0.5 mg/mL), 365 nm LED, 50 mW/cm<sup>2</sup>, 10 °C.

## Reaction scheme



ChemDraw file linked: [NB-351-SrTiO3-photocatalytic H2O splitting.cdxml](#)

## Literature/reference experiments

Literature	/
Reproduction	/
Similar experiments	<a href="#">SrTiO3 - NB-316: Liquid phase H2 and O2 of RhCrO<sub>x</sub>,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 20 °C</a>

## Reagents

Name	CAS Number / Experiment Number	Inventor y number	Amount [mmol]	Equival ents	Mass <sub>theo</sub> [mg]	Mass <sub>exp</sub> [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	Pressure [bar]	Concentratio n [mM]
milli-Q H <sub>2</sub> O	/	/	/	/	/	/	/	0.998	25.00 + 25 (for calibration)	/	/
Al:SrTiO <sub>3</sub> RhCrO <sub>x</sub> (EA-358)	<a href="#">SrTiO<sub>3</sub> - EA-358: Modification of Al:SrTiO<sub>3</sub> (EA-354) via deposition of Rh, Cr oxide co-catalyst, 350°C, 1h, Upscaling (3.33x)</a>	/	/	/	12.50	12.67	/	/	/	/	/
Hydrogen	1333-74-0	/	/	/	/	/	/	/	2 balloons (ca. 2 L)	approx. 1	/

# Excel sheet for reagent calculation

/

## Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-314: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup](#)

	Name
Used Set-up	<a href="#">Equipment - Advanced irradiation setup V1.0 I</a>
Irradiation setup number	<a href="#">Equipment - Irradiation setup 4 (CEEC II, E002)</a>

	Light Source Name	Power Source Name	Wavelength [nm]	Power Setting [mW]	Analog Setting [0.00 - 10.00]
First light source	<a href="#">Light Source - UHP LED 365 nm-4</a>	<a href="#">Power Sources - BLS-18000-1 4</a>	365	56	0.19

Used beam combiner [Name or None]	/
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	10
Stirring speed [rpm]	500
Irradiation start: 1. Firesting [relative to start log] 2. Unisense	1. 603 s 2. 21:26:03
Irradiation stop: 1. Firesting [relative to start log] 2. Unisense	1. 1518 s 2. 21:41:17

## O<sub>2</sub>/H<sub>2</sub> sensor equipment

	Equipment	Used protocol
Used Firesting	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)	Protocol - Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software
Used O <sub>2</sub> sensor	Equipment - Robust probe for liquid O <sub>2</sub> measurment	Protocol - In-situ hydrogen and oxygen measurment in H <sub>2</sub> /O <sub>2</sub> reactor
Used H <sub>2</sub> sensor	Equipment - H <sub>2</sub> UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurment in H <sub>2</sub> /O <sub>2</sub> reactor

## Procedure/observations

Date	Time	Step	Observations	Pictures/Files
18.11.2025		The experiment was done according to <a href="#">Protocol - In-situ hydrogen and oxygen measurment in H<sub>2</sub>/O<sub>2</sub> reactor</a> Important steps and deviations are listed below	/	/
	9:13-10:28	Conditioning of H <sub>2</sub> sensor	<b>NB-351-Logger1</b> 1.6 mV at the end of polarization procedure	<a href="#">NB-351.ulong</a> <a href="#">NB-351-Logger1-pre-polarization.csv</a> <a href="#">NB-351-Logger1-pre-polarization.bmp</a>
	16:49	Assembling the setup for calibration (25 mL of water was added using graduated cylinder), (LAUDA set to 10 °C) done according to <a href="#">Protocol - Liquid phase calibration of H<sub>2</sub> UniAmp sensor</a> with H <sub>2</sub> bubbling.	/	/
	16:50-17:03	Keeping the system under stirring to reach 10 °C (T control by PT1000).	/	/
	17:04	Start of O <sub>2</sub> logging.	<b>NB-351-Ch2-2</b>	<a href="#">2025-11-18_170449_NB-351-Ch2-2.txt</a> <a href="#">2025-11-18_170449_NB-351-Ch2-2.png</a>
	17:06	Degassing was started.	/	<a href="#">20251118_170710-degassing of water.jpg</a>

	17:07	Start of H2 logging.	<b>NB-351-Logger2</b> offset -2 mV	NB-351.ulog NB-351-Logger2-calibration.csv NB-351-2point calibration.bmp NB-351-Logger2-calibration step.bmp
	17:26	Moving cannula to gas phase above the liquid level.	/	/
	17:29	Introducing H2 sensor into the reactor under Ar flow.	/	/
	17:30	0 ppm was taken.	/	/
	17:37	Degassing was stopped.	/	/
	17:38	H2 bubbling of the reactor was started	/	/
	17:46	Introducing second H2 balloon.	Drop in H2 value after addition of the 2nd balloon, but afterwards increased and could reach steady state.	20251118_174602-after introducing 2nd H2 balloon.jpg
	17:55	1.000.000 ppm point was taken and calibration was saved	640 mV, slope: 0.750, 856 uM	/
	18:02	Stop of H2 logging.	/	/
	18:02	Stop of O2 logging.	/	/
	ca. 18:10	Deassembling the setup, drying the reactor with acetone and compressed air .	/	/
		<b>Sample preparation</b>		
	19:15	Weighing EA-358 photocatalyst in a 50 mL vial.	Creamy solid.	/
	19:20	Addition of 25 mL H2O to the vial via graduated cylinder.	/	/

	19:22-25	The suspension was vortexed for 3 min ( <a href="#">Equipment - VWR® VV3, Vortex Mixer</a> , stage 4/6), covered with Al foil before further use.	/	<a href="#">20251118_192609-suspension after vortex.jpg</a>
		Continue in <a href="#">Protocol - In-situ hydrogen and oxygen measurment in H2/O2 reactor</a> from step 6		
	ca. 19:35	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	/
	19:40	Assembling the setup.	Currently, stopper instead of H <sub>2</sub> sensor, PT100, PT1000 and O <sub>2</sub> robust probe are inside the reactor immersed in the liquid phase	/
	19:40-50	Keeping suspension at 10 °C, T control by PT1000.	/	/
	19:50	<a href="#">Start of O2 logging.</a>	<b>NB-351-Ch2-3</b>	<a href="#">2025-11-18_195017_NB-351-Ch2-3.txt</a> <a href="#">2025-11-18_195017_NB-351-Ch2-3.png</a>
	19:52	The degassing was started	/	<a href="#">20251118_195316-degassing of the suspension.jpg</a>
	20:30	Cannula was transferred to gas phase, above the suspension.	/	/
	20:34	H <sub>2</sub> sensor was added in Ar counterflow.	/	/
	20:36	The degassing was stopped by removing the cannula and closing the valve.	/	/
	20:38	<a href="#">Stop of O2 logging.</a>	/	/
	20:38	<a href="#">Start of O2 logging.</a>	<b>NB-351-Ch2-4</b>	<a href="#">2025-11-18_203830_NB-351-Ch2-4.txt</a> <a href="#">2025-11-18_203830_NB-351-Ch2-4.png</a>
	20:38	<a href="#">Start of H2 logging.</a>	<b>NB-351-Logger3</b> High H2 signal after introducing sensor to the reactor.	<a href="#">NB-351.ulong</a> <a href="#">NB-351-Logger3.csv</a> <a href="#">NB-351-Logger3.bmp</a>

	20:43	Stop of H2 and O2 loggings.	Removing the sensor from the reactor, introducing once again, checking cable management.	/
	20:48	Start degassing once again.	<b>NB-351-Ch2-5</b>	<a href="#">2025-11-18_204849_NB-351-Ch2-5.txt</a> <a href="#">2025-11-18_204849_NB-351-Ch2-5.png</a>
	21:07	Moving cannula to gas phase above the liquid level.	/	/
	21:10	Introducing H2 sensor under Ar flow.	/	/
	21:13	Removing cannula, closing the valve.	/	/
	21:15	Stop of O2 logging.	/	/
	21:16	Start of O2 logging.	<b>NB-351-Ch2-6</b>	<a href="#">2025-11-18_211601_NB-351-Ch2-6.txt</a> <a href="#">2025-11-18_211601_NB-351-Ch2-6.png</a>
	21:16	Start of H2 logging.	<b>NB-351-Logger4</b>	<a href="#">NB-351.ulong</a> <a href="#">NB-351-Logger4-during irradiation.csv</a> <a href="#">NB-351-Logger4-during irradiation.bmp</a>
	21:16-26	Equilibration time.	/	/
	21:26	The irradiation was started	/	<a href="#">20251118_212721-after start of irradiation.jpg</a>
	21:41	The irradiation was stopped.	/	/
	21:41-51	Equilibration time.	/	/
	21:51	Stop of O2 and H2 logging.	/	/

	ca. 22:00	Deassembling the setup, cleaning the reactor.	Tips of the sensors and reactor were covered with attached photocatalyst particles. Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode).	<a href="#">20251118_215222-after irradiation.jpg</a>
--	-----------	---	--	---

## Analysis

Used calibration for Firesting: [20250611-BOLA-fitting-liquid-phase-trace-oxygen-sensor-H2-O2 reactor.ini](#)

Used calibration for UniSense: NB-351-Logger2

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
18.11.2025	9:13	NB-351-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	<a href="#">NB-351.ulong</a> <a href="#">NB-351-Logger1-pre-polarization.csv</a>	/	<a href="#">NB-351-Logger1-pre-polarization.bmp</a>	/	Pre-polarization of H2 sensor.
	17:07	NB-351-Logger2	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	<a href="#">NB-351.ulong</a> <a href="#">NB-351-Logger2-calibration.csv</a>	/	<a href="#">NB-351-2point calibration.bmp</a> <a href="#">NB-351-Logger2-calibration step.bmp</a>	/	Calibration of H2 sensor, at 10 <sup>6</sup> ppm signal 640 mV, slope 0.750
	20:38	NB-351-Logger3	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	<a href="#">NB-351.ulong</a> <a href="#">NB-351-Logger3.csv</a>	/	<a href="#">NB-351-Logger3.bmp</a>	/	High signal after introducing H2 snsor to the suspension --> needs to be eliminated.
	21:16	NB-351-Logger4	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	<a href="#">NB-351.ulong</a> <a href="#">NB-351-Logger4-during irradiation.csv</a>	<a href="#">NB-351-O2 and H2 curve.py</a>	<a href="#">NB-351-Logger4-during irradiation.bmp</a> <a href="#">NB-351-O2 and H2 curves.png</a>		H2 evolution during irradiation.
	17:04	NB-351-Ch2-2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	<a href="#">2025-11-18_170449_NB-351-Ch2-2.txt</a>	/	<a href="#">2025-11-18_170449_NB-351-Ch2-2.png</a>	/	Degassing of water followed by calibration of H2 sensor.
	19:50	NB-351-Ch2-3	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	<a href="#">2025-11-18_195017_NB-351-Ch2-3.txt</a>	/	<a href="#">2025-11-18_195017_NB-351-Ch2-3.png</a>	/	First degassing.

	20:38	NB-351-Ch2-4	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-18_203830_NB-351-Ch2-4.txt	/	2025-11-18_203830_NB-351-Ch2-4.png	/	Check of O2 value during 1st logging - first attempt to start pre-reaction baseline.
	20:48	NB-351-Ch2-5	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-18_204849_NB-351-Ch2-5.txt	/	2025-11-18_204849_NB-351-Ch2-5.png	/	Degassing after second introduction of H2 sensor in the suspension.
	21:16	NB-351-Ch2-6	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-18_211601_NB-351-Ch2-6.txt	NB-351-O2 and H2 curve.py	2025-11-18_211601_NB-351-Ch2-6.png NB-351-O2 and H2 curves.png	/	O2 evolution during irradiation.

## Results

Simultaneous H<sub>2</sub> and O<sub>2</sub> measurements of irradiated suspension of EA-358 (0.5 mg/mL) in O<sub>2</sub>/H<sub>2</sub> photoreactor under 365 nm irradiation (50 mW/cm<sup>2</sup>, 10 °C, 15 min) were performed.

## Linked experiments

SrTiO<sub>3</sub> - NB-316: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrOx,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 20 °C

SrTiO<sub>3</sub> - NB-318: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrOx,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm<sup>2</sup>, 20 °C

SrTiO<sub>3</sub> - NB-322: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrOx,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm<sup>2</sup>, 20 °C (reproduction NB-318)

SrTiO<sub>3</sub> - NB-325: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrOx,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm<sup>2</sup>, 20 °C

SrTiO<sub>3</sub> - NB-330: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrOx,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 30 °C

SrTiO<sub>3</sub> - NB-334: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrOx,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 30 °C (reproduction NB-330)



## Linked resources

Equipment - [VWR® VV3, Vortex Mixer](#)

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Advanced power measurment chamber V1.0 I](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)

Light Source - [UHP LED 365 nm-4](#)

Power Sources - [BLS-18000-1 4](#)


Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)

Protocol - [In-situ hydrogen and oxygen measurment in H2/O2 reactor](#)

## Attached files

NB-351-O2 and H2 curves.png

sha256: b32963d36d3dfc263c5fc52f3048a48428ab4b200e0326e20958e18bb39052aa



NB-351-O2 and H2 curve.py

sha256: 111901ace35b9834fc3831f4ca512acce87aa02a313c012cc267ef16c6b47287

NB-351-SrTiO3-photocatalytic H2O splitting.cdxml

sha256: 53cc5ac17c1daa623cac433566de7e5c8d37e57d743d4e1fc906582213d107c9

NB-351-SrTiO3-photocatalytic H2O splitting.png

sha256: c5fed95bcf7614912c6697125be7fbbf8c073e07ccbab0c5516f39dff01bf71a




20251118\_170710-degassing of water.jpg

sha256: 32617286e202023bc95de2719dd865934e5b555897ae216bee8d7629d3114ba3



20251118\_174602-after introducing 2nd H2 balloon.jpg

sha256: c70a496e1a1a4532563a3041c43ef19690058a8d9d84e6152a43899a9dbe9703



20251118\_195316-degassing of the suspension.jpg

sha256: e814b635c8fd867c0e6f96ae9b8505a4bc8cedec369328f2b6acea8cb9712af2



20251118\_192609-suspension after vortex.jpg

sha256: 49b93a4d0a94d477eeb2e5855239d948b5e171480b7c586efe3afa36197d5099



20251118\_215222-after irradiation.jpg

sha256: b8a8c06ce0cd465e2f6e65df49fc3eda8b435364527d631efb5761ebcd0b88a1



20251118\_212721-after start of irradiation.jpg

sha256: 8726f844e90506532841db2a770e4e3ceb09ab2b71d9454be63dc0756c97a474



NB-351.ulog


sha256: 3f0cc580cad02bcd979fff276382cdae42c049aaec1716eefdaedd7fb46ed6b

NB-351-Logger4-during irradiation.csv

sha256: 9af569a2aeabb4aceb3f232b65ad704accd48fb9933dcc44f2cf719c480e7d7e

NB-351-Logger4-during irradiation.bmp

sha256: 656c1838ac699321ff27928e4ecb570c9eedf95482ee186b9ae87cfa02900c66



NB-351-Logger3.bmp

sha256: 05a096c8ac1a2de91a099d898c706a2550ba0cce7e72c4c2f61a03779f9f9c91



NB-351-Logger3.csv

sha256: f5d02c2169fa60aaba4a651449fde413e44b4db2783a5b816dff2795739aa7a5

NB-351-Logger2-calibration.csv

sha256: a1513e779e6fcd09bd07d0b200f40674a4652d1e99d8bbb5997e0de27125fdab

NB-351-2point calibration.bmp

sha256: 8485642c00c99066d352e9d32fb15e269125cf5a326845994cb80ceb3af5461e




NB-351-Logger2-calibration step.bmp

sha256: 5718189704bbb94d88986b44f2199b10fd24d80c7074ae95be8ce97433f70495



NB-351-Logger1-pre-polarization.bmp

sha256: 88431ae22e5e4ac3150fb86f2489144c9583bc576b4507004a6bbbbb71a966ca6



NB-351-Logger1-pre-polarization.csv


sha256: 031f39f400b0c50db4243745a19f8953fbbda8c7d0c585e9a692acd97ed66b26

2025-11-18\_211601\_NB-351-Ch2-6.txt

sha256: e792dabde3fdadbac8a3e3595bb26261415b7f801032344cffe0b94ad9072d0

2025-11-18\_211601\_NB-351-Ch2-6.png

sha256: f28085422f632509b32204e4ae86b476d5029949c1e563ae36842199202c75ed




2025-11-18\_204849\_NB-351-Ch2-5.txt

sha256: 753a6cfbf3429f478d69c4e2202682fb27683c18674c29dad14ac6ec96e0c569

2025-11-18\_204849\_NB-351-Ch2-5.png

sha256: 21ae35b1615c6b633af39ba5e8a29d668a9f58f32297b9f0000979dc3e1c6d8f



2025-11-18\_203830\_NB-351-Ch2-4.txt

sha256: 99c5ff37246ff5c51d30becff7b04e23b91d4a8a4ae6e250d0bb0fef08231237

2025-11-18\_203830\_NB-351-Ch2-4.png

sha256: a20ca39baac008debabeb0dcbf6a6098d89622ae0b7aea740dfac4ebd3a2ebdb



2025-11-18\_195017\_NB-351-Ch2-3.txt

sha256: da0913bb3c57031e1b7cfe8f9c8a09bd9d6517294bd7c858a780c47205747f2a

2025-11-18\_195017\_NB-351-Ch2-3.png

sha256: 6a18162fee9a30b4741f1b7972e13814a0a5a4a38d586c323e388b7d6f7587a6



2025-11-18\_170449\_NB-351-Ch2-2.png

sha256: 441cc6e931cee0956f7810f9fc24b92e7024660c1ad907b8a4ffb6778284982b



2025-11-18\_170449\_NB-351-Ch2-2.txt

sha256: 1a06ed1c901f192b5f411a607abcc0d30aa503e2a6fd3bb662329f5bed1a6b64



Unique eLabID: 20251118-c8339e87c70c80570c8da352fa09bc971d1de7e2  
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3476>