

NB-332: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (NB-321, 0.05 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C (reproduction NB-327)

Date: 2025-11-03

Tags: [O2](#) [Test](#) [Calibration](#) [NB](#) [Firing](#) [O2 sensor](#) [H2](#) [SrTiO3](#) [Unisense](#) [RhCrO3:Al:SrTiO3](#) [H2 Sensor](#) [temperature](#) [In situ](#) [Trace range robust](#) [oxygen sensor](#) [photocatalysis](#) [Unisense](#) [normal range](#)

Category: SrTiO₃


Status: Done

Created by: Nadzeya Brezhneva

Objectives

Reproduction of NB-327: Simultaneous detection of H₂ and O₂ evolution in liquid phase for irradiated suspension of Rh_xCr_yO₃:Al:SrTiO₃ suspension (NB-321, **0.05 wt% Rh,Cr**, 0.5 mg/mL), 365 nm LED, 50 mW/cm², 20 °C.

Reaction scheme



ChemDraw file linked: [NB-327-SrTiO3-photocatalytic H2O splitting.cdxml](#)

Literature/reference experiments

Literature	/
Reproduction	SrTiO3 - NB-327: Liquid phase H2 and O2 of RhCrO_x,Al:SrTiO3 (NB-321, 0.05 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C
Similar experiments	SrTiO3 - NB-326: Liquid phase H2 and O2 of RhCrO_x,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm², 20 °C (reproduction NB-325)

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	Concentration [mM]
milli-Q H ₂ O	/	/	/	/	/	/	/	0.998	25	/
Al:SrTiO ₃ RhCrO _x (NB-321)	SrTiO3 - NB-321: Modification of EA-354 (SrTiO3:Al, upscaled batch) with Rh, Cr oxide cocatalyst (0.05 wt%)	/	/	/	12.50	12.51	/	/	/	/

Excel sheet for reagent calculation

/

Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-314: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup](#)

	Name
Used Set-up	Equipment - Advanced irradiation setup V1.0 I
Irradiation setup number	Equipment - Irradiation setup 4 (CEEC II, E002)

	Light Source Name	Power Source Name	Wavelength [nm]	Power Setting [mW]	Analog Setting [0.00 - 10.00]
First light source	Light Source - UHP LED 365 nm-4	Power Sources - BLS-18000-1 4	365	56	0.19

Used beam combiner [Name or None]	/
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	20
Stirring speed [rpm]	500
Irradiation start: 1. Firesting [relative to start log] 2. Unisense	1. 615 s 2. 22:29:35
Irradiation stop: 1. Firesting [relative to start log] 2. Unisense	1. 1525 s 2. 22:44:45

O₂/H₂ sensor equipment

	Equipment	Used protocol
Used Firesting	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)	Protocol - Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software
Used O ₂ sensor	Equipment - Robust probe for liquid O ₂ measurment	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor
Used H ₂ sensor	Equipment - H ₂ UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor

Procedure/observations

Date	Time	Step	Observations	Pictures/Files
		Calibration from SrTiO ₃ - NB-331: Liquid phase H ₂ and O ₂ of RhCrOx,Al:SrTiO ₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm ² , 20 °C (reproduction NB-316) I I I was used.	/	/
03.11.2025		Sample preparation		
	ca.21:05	Weighing photocatalyst in a 50 mL vial.	Slightly creamy solid	/
		Important note: before weighing, the solid was dispersed inside the vial manually with Smartspatula - to break lumps manually, but to avoid mortaring.		
	21:13	Addition of 25 mL H ₂ O to the vial via graduated cylinder.	/	/
	21:18-21	The suspension was vortexed for 3 min (Equipment - VWR® VV3, Vortex Mixer, stage 4/6), covered with Al foil before further use.	/	20251103_212138-suspension after vortex.jpg
		Continue in Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor from step 6		
	ca. 21:25	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	/

	21:30	Assembling the setup.	Currently, stopper instead of H ₂ sensor, PT100, PT1000 and O ₂ robust probe are inside the reactor immersed in the liquid phase	20251103_212855-before degassing and irradiation.jpg
	21:30	Start of O2 logging.	NB-332-Ch2-1	2025-11-03_213054_NB-332-Ch2-1.txt 2025-11-03_213054_NB-332-Ch2-1.png
	21:32	The degassing was started	/	20251103_213303-degassing of the suspension.jpg
	22:09	Cannula was transferred to gas phase, above the suspension.	/	/
	22:11	H ₂ sensor was added in Ar counterflow.	/	/
	ca. 22:17	The degassing was stopped by removing the cannula and closing the valve.	/	/
	22:18	Stop of O2 logging.	/	/
	22:19	Start of O2 logging.	NB-332-Ch2-2	2025-11-03_221920_NB-332-Ch2-2.txt 2025-11-03_221920_NB-332-Ch2-2.png
	22:19	Start of H2 logging.	NB-332-Logger1	NB-332.ulong NB-332-Logger1.csv NB-332-Logger1.bmp
	22:19-29	Equilibration time.	/	/
	22:29	The irradiation was started	/	20251103_223003-after start of irradiation.jpg
	22:44	The irradiation was stopped.	/	/
	22:44-54	Equilibration time.	/	/
	22:54	Stop of O2 and H2 logging.	/	/

	ca. 23:00	Deassembling the setup, cleaning the reactor.	Tips of the sensors and reactor were covered with attached photocatalyst particles, bubbles attached to the surface of the sensors tips. Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode).	/
--	--------------	---	---	---

Analysis

Used calibration for Firesting: [20250611-BOLA-fitting-liquid-phase-trace-oxygen-sensor-H2-O2 reactor.ini](#)

Used calibration for UniSense: NB-331-Logger3

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
03.11.2025	22:19	NB-332-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-332.ulong NB-332-Logger1.csv	NB-332-O2 and H2 curve.py	NB-332-Logger1.bmp NB-332-O2 and H2 curve.png	SrTiO3 - NB-327: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (NB-321, 0.05 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C	Several artefacts during measurements.
	21:30	NB-332-Ch2-1	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-03_213054_NB-332-Ch2-1.txt	/	2025-11-03_213054_NB-332-Ch2-1.png	/	Degassing of the suspension.
	22:19	NB-332-Ch2-2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-03_221920_NB-332-Ch2-2.txt	NB-332-O2 and H2 curve.py	2025-11-03_221920_NB-332-Ch2-2.png NB-332-O2 and H2 curve.png	SrTiO3 - NB-327: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (NB-321, 0.05 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C	Slight step at ca. 22:41.

Results

Simultaneous H₂ and O₂ measurements of irradiated suspension of NB-321 (0.05 wt% RhCr, 0.5 mg/mL) in O₂/H₂ photoreactor under 365 nm irradiation (50 mW/cm², 20 °C, 15 min) were performed, reproduction NB-327.

H₂ level at the end of irradiation - ca. 34 umol/L, O₂ level - ca. 33 umol/L.

Linked experiments

SrTiO₃ - NB-316: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

SrTiO₃ - NB-318: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm², 20 °C

SrTiO₃ - NB-321: Modification of EA-354 (SrTiO₃:Al, upscaled batch) with Rh, Cr oxide cocatalyst (0.05 wt%)

SrTiO₃ - NB-322: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm², 20 °C (reproduction NB-318)

SrTiO₃ - NB-325: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm², 20 °C

SrTiO₃ - NB-326: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm², 20 °C (reproduction NB-325)

SrTiO₃ - NB-327: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (NB-321, 0.05 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

SrTiO₃ - NB-331: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C (reproduction NB-316) I I I

Linked resources

Equipment - [VWR® VV3, Vortex Mixer](#)

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)

Protocol - [Hydrogen Measurement Using GC](#)

Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)

Protocol - [Liquid phase calibration of H₂ UniAmp sensor](#)

Protocol - [In-situ hydrogen and oxygen measurment in H₂/O₂ reactor](#)


Attached files

NB-332-O2 and H₂ curve.py

sha256: c398cb8b06b7201511f21ad1b609732315d460defa453e6410efae9c018770cf

NB-332-O2 and H2 curve.png

sha256: b064232db282d423a3799860f72804c4242aedaf062cb04f06cfb9e85be38998



20251103_213303-degassing of the suspension.jpg

sha256: 3b43522c4d074f133ca4324986296fa94a1a37a18b929eccc152044647e07d00



20251103_212855-before degassing and irradiation.jpg

sha256: 11b680fadb61907cede394bce1678a928b9dbd498618b9ee822a745196de56a9



20251103_212138-suspension after vortex.jpg

sha256: 27b7b0d36e2c94c24c05c3e37fff36459b9b348840443ffb8e6a73bde998adb2



20251103_223003-after start of irradiation.jpg

sha256: b8ca5afd65e7359ee6fa38bdd566db2a198812407dd00056c32f711acbf12dda




20251103_225558-after irradiation.jpg

sha256: 383fc4fff88e3c9a75b8fe5bcc9f6de488c2d5a16b2f71526e762874d1e6ceaa



2025-11-03_221920_NB-332-Ch2-2.png

sha256: 1deeb6769a4f7d3f433c7f770ca2ef891a1df3ccfcc24fccda9c1237a719dfb7




2025-11-03_221920_NB-332-Ch2-2.txt

sha256: 86c87cea99a8b1baf6177d12d1b98c5bba82400285a2ec61a101e015564b73f0

2025-11-03_213054_NB-332-Ch2-1.png

sha256: 03cd3d092b4261cca78357dd2435b4369451c7e13e2c0d39d72c4e03bc35b532




2025-11-03_213054_NB-332-Ch2-1.txt

sha256: 087a14286a9c01fbdaaffcb621f3864f12335a1864970a4aca8ee3b4723baf06f

NB-332-Logger1.bmp

sha256: 670967f7ac98296202c822e619745008c0492e5ce2b7dafa19f655d52be4b844



NB-332.ulog

sha256: bf190960bece0b6daf8e2fda4fb26096dfdd91e377cd1986380928ae85fa16c6

NB-332-Logger1.csv

sha256: bc5f29ff72252166fedd61ce2fba4e652a3b4b15825f71a2074c4be3ca976809



Unique eLabID: 20251103-999b8c746d43d0b953efbacd9be026aa38c79807
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3324>