

NB-334: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 30 °C (reproduction NB-330)

Date: 2025-11-05

Tags: O₂ Test Calibration NB Firing
O₂ sensor H₂ SrTiO₃ Unisense
RhCrO₃:Al:SrTiO₃ H₂ Sensor
temperature In situ H₂/O₂ reactor Trace
range robust oxygen sensor
photocatalysis Unisense normal range

Category: SrTiO₃


Status: Done

Created by: Nadzeya Brezhneva

Objectives

Reproduction of NB-330: simultaneous detection of H₂ and O₂ evolution in liquid phase for irradiated suspension of Rh_xCrO_y:Al:SrTiO₃ suspension (EA-358, 0.5 mg/mL), 365 nm LED, 50 mW/cm², 30 °C.

Reaction scheme



ChemDraw file linked: [NB-330-SrTiO3-photocatalytic H2O splitting.cdxml](#)

Literature/reference experiments

Literature	/
Reproduction	SrTiO₃ - NB-330: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 30 °C
Similar experiments	SrTiO₃ - NB-316: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	Pressue [bar]	Concentration [mM]
milli-Q H ₂ O	/	/	/	/	/	/	/	0.998	25 + 25 (for calibration)	/	/
Al:SrTiO ₃ RhCrO _x (EA-358)	SrTiO₃ - EA-358: Modification of Al:SrTiO₃ (EA-354) via deposition of Rh, Cr oxide co-catalyst, 350°C, 1h, Upscaling (3.33x)	/	/	/	12.50	12.51	/	/	/	/	/

Hydrogen	1333-74-0	/	/	/	/	/	/	/	2 ballons (ca. 2 L)	ca. 1	/
----------	-----------	---	---	---	---	---	---	---	------------------------	-------	---

Excel sheet for reagent calculation

/

Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-314: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup](#)

	Name
Used Set-up	Equipment - Advanced irradiation setup V1.0 I
Irradiation setup number	Equipment - Irradiation setup 4 (CEEC II, E002)

	Light Source Name	Power Source Name	Wavelength [nm]	Power Setting [mW]	Analog Setting [0.00 - 10.00]
First light source	Light Source - UHP LED 365 nm-4	Power Sources - BLS-18000-1 4	365	56	0.19

Used beam combiner [Name or None]	/
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	30
Stirring speed [rpm]	500
Irradiation start: 1. Firesting [relative to start log] 2. Unisense	1. 601 s 2. 17:16:13

Irradiation stop: 1. Firing [relative to start log] 2. Unisense	1. 1513 s 2. 17:41:21
--	--------------------------

O₂/H₂ sensor equipment

	Equipment	Used protocol
Used Firing	Equipment - Firing Fiber-Optic Oxygen Meter 2 Channel (Firing 2)	Protocol - Operation of Firing Fiber-Optic Oxygen Meter 2 Channel Software
Used O ₂ sensor	Equipment - Robust probe for liquid O ₂ measurement	Protocol - In-situ hydrogen and oxygen measurement in H ₂ /O ₂ reactor
Used H ₂ sensor	Equipment - H ₂ UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurement in H ₂ /O ₂ reactor

Procedure/observations

Date	Time	Step	Observations	Pictures/Files
05.11.2025		The experiment was done according to Protocol - In-situ hydrogen and oxygen measurement in H₂/O₂ reactor Important steps and deviations are listed below	/	/
	10:08-11:15	Conditioning of H ₂ sensor	NB-334-Logger1 1.7 mV at the end of polarization procedure	NB-334.ulong NB-334-Logger1-pre-polarization.csv NB-334-Logger1-pre-polarization.bmp
	11:20-45	Assembling the setup for calibration (25 mL of water was added using graduated cylinder), (LAUDA set to 30 °C) done according to Protocol - Liquid phase calibration of H₂ UniAmp sensor with H ₂ bubbling.	/	/
	11:45-12:14	Keeping the system at 30 °C (control by PT1000 sensor immersed in water inside the reactor).	/	/
	12:14	Start of O ₂ logging.	NB-334-Ch2-1	2025-11-05_121425_NB-334-Ch2-1.txt 2025-11-05_121425_NB-334-Ch2-1.png

	12:17	Degassing was started.	ca. 13:00 - increasing Ar flow a bit	20251105_121826-during degassing of water.jpg
	12:18	Start of H2 logging.	NB-330-Logger2 offset -2 mV	NB-334.uog NB-334-Logger2-calibration.csv NB-334-Logger-2point-calibration.bmp NB-334-calibration-step.bmp
	13:04	Moving cannula to the gas phase above the liquid	/	/
	13:05	Introducing H2 sensor into the reactor under Ar flow.	Offset -3 mV	/
	13:08	0 ppm was taken.	/	/
	13:10	Stop degassing by removing the cannula, closing the valve.	/	/
	13:17	H2 bubbling of the reactor was started by introducing the first H2 balloon.	ca. 13:20 - very slight moving of H2 sensor, slight change in the slope of the curve	20251105_131800-introducing 1 H2 balloon.jpg
	13:25	Introducing second H2 balloon.	/	20251105_132700-introducing 2nd H2 balloon.jpg
	13:31	1.000.000 ppm point was taken and calibration was saved	999 mV, slope: 1.315, 759.873 μ M	20251105_133102-H2 table.jpg
	13:32	Stop of H2 logging.	/	/
	13:32	Stop of O2 logging.	/	/
	13:35	Deassembling the setup, drying the reactor with acetone and compressed air .	/	/
		Sample preparation		
	15:45	Weighing EA-358 photocatalyst in a 50 mL vial.	Creamy solid, manually breaking lumps inside the vial using Smartspatula	/
	15:50	Addition of 25 mL H2O to the vial via graduated cylinder.	/	/

	15:51-54	The suspension was vortexed for 3 min (Equipment - VWR® VV3, Vortex Mixer , stage 4/6), covered with Al foil before further use.	/	20251105_155807-suspension after vortex.jpg
		Continue in Protocol - In-situ hydrogen and oxygen measurment in H2/O2 reactor from step 6		
	ca. 16:00	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	20251105_160440-suspension inside the reactor.jpg
	ca. 16:05-10	Assembling the setup.	Currently, stopper instead of H ₂ sensor, PT100, PT1000 and O ₂ robust probe are inside the reactor immersed in the liquid phase /	
	16:10-25	Stirring of the suspension at 30 °C.	/	20251105_162423-assembling setup.jpg
	16:25	Start of O2 logging.	NB-334-Ch2-2	2025-11-05_162517_NB-334-Ch2-2.txt 2025-11-05_162517_NB-334-Ch2-2.png
	16:28	The degassing was started	16:43 - The glass connector popped up from the tubing in Ar line - -> problem eliminated using screw driver and heat gun (to provide better fixing of the connector)	20251105_162942-degassing of the suspension.jpg
	16:54	Continuation of degassing.	/	/
	17:08	Cannula was transferred to gas phase, above the suspension.	/	/
	17:11	H ₂ sensor was added in Ar counterflow.	/	20251105_175313-after introducing H2 sensor.jpg
	17:14	The degassing was stopped by removing the cannula and closing the valve.	/	/
	17:15	Stop of O2 logging.	/	/
	17:16	Start of O2 logging.	NB-334-Ch2-3	2025-11-05_171608_NB-334-Ch2-3.txt 2025-11-05_171608_NB-334-Ch2-3.png

	17:16	Start of H2 logging.	NB-334-Logger3	NB-334.ulong NB-334-Logger3-during irradiation.csv NB-334-Logger3-during irradiation.bmp
	17:16-26	Equilibration time.	/	/
	17:26	The irradiation was started	/	20251105_172634-after start of irradiation.jpg
	17:41	The irradiation was stopped.	/	/
	17:41-51	Equilibration time.	/	/
	17:51	Stop of O2 and H2 logging.	/	/
	ca. 17:55	Deassembling the setup, cleaning the reactor.	Tips of the sensors and reactor were covered with attached photocatalyst particles. Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode).	/

Analysis

Used calibration for Firesting: [20250611-BOLA-fitting-liquid-phase-trace-oxygen-sensor-H2-O2 reactor.ini](#)

Used calibration for UniSense: NB-334-Logger2

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
05.11.2025	10:08	NB-334-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-334.ulong NB-334-Logger1-pre-polarization.csv	/	NB-334-Logger1-pre-polarization.bmp	/	1.7 mV at the end of polarization
	12:18	NB-334-Logger2	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-334.ulong NB-334-Logger2-calibration.csv	/	NB-334-calibration-step.bmp NB-334-Logger-2point-calibration.bmp	/	Calibration with 1.000.000 ppm pint: slope of 1.315

17:16	NB-334-Logger3	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-334.ulog NB-334-Logger3-during irradiation.csv	NB-334-O2 and H2 curve.py	NB-334-Logger3-during irradiation.bmp NB-334-O2 and H2 curve.png	SrTiO3 - NB-330: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 30 °C	worked, clean signal response, ca. 58 uM at the end of irradiation
12:14	NB-334-Ch2-1	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-05_121425_NB-334-Ch2-1.txt	/	2025-11-05_121425_NB-334-Ch2-1.png	/	O2 measurements during H2 calibration (degassing with Ar and purging with H2 afterwards).
16:25	NB-334-Ch2-2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-05_162517_NB-334-Ch2-2.txt	/	2025-11-05_162517_NB-334-Ch2-2.png	/	Degassing of the photocatalyst suspension.
17:16	NB-334-Ch2-3	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-05_171608_NB-334-Ch2-3.txt	NB-334-O2 and H2 curve.py	2025-11-05_171608_NB-334-Ch2-3.png NB-334-O2 and H2 curve.png	SrTiO3 - NB-330: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 30 °C	worked, ca. 47 umol/L at the end of irradiation

Results

Successful reproduction of NB-330: simultaneous H₂ and O₂ measurements of irradiated suspension of EA-358 (0.5 mg/mL) in O₂/H₂ photoreactor under 365 nm irradiation (50 mW/cm², **30 °C**, 15 min) were performed.

H₂ level at the end of irradiation - ca. 58 umol/L, O₂ level - ca. 47 umol/L.

Linked experiments

SrTiO3 - NB-316: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C

SrTiO3 - NB-318: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm2, 20 °C

SrTiO3 - NB-322: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm2, 20 °C (reproduction NB-318)

SrTiO3 - NB-325: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm2, 20 °C

SrTiO3 - NB-330: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 30 °C

Linked resources

Equipment - [VWR® VV3, Vortex Mixer](#)

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Advanced power measurment chamber V1.0 I](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)

Light Source - [UHP LED 365 nm-4](#)

Power Sources - [BLS-18000-1 4](#)

Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)

Protocol - [In-situ hydrogen and oxygen measurment in H2/O2 reactor](#)


Attached files

NB-334-O2 and H2 curve.py

sha256: fced42a69b86ea7e617afaa60dc65710d344b036e1cfeef9d21e2bae89cc8a26

NB-334-O2 and H2 curve.png

sha256: b0c13d4bd13becf4b4fc643017209d4fd96d9920a0bd32af01a7099dc64749f2



20251105_121826-during degassing of water.jpg

sha256: 5768da2fcbfe2bfa3c15b4b99115c038583992b530a7ce9b891205dd6f56769f



20251105_131800-introducing 1 H2 balloon.jpg

sha256: 8f07235e77d08b88e82ea982baec7f9f01584720d885bf0bc697ec39ef46fc0f




20251105_133102-H2 table.jpg

sha256: c7840f645b119ce71831c50930695bed9fae452f56587ccff0561a079a793f5e



20251105_132700-introducing 2nd H2 balloon.jpg

sha256: 376ca4b649cfc067096f6e372a433c098e7924709b385ed793fce69c3c97adcc



20251105_155807-suspension after vortex.jpg

sha256: 0e9ef14e2614852ccd31127ad270ede074d65cd7f5719917bdd6a3582eacec32



20251105_160440-suspension inside the reactor.jpg

sha256: a5c0908388648da72615d32dfe9323c146e14c41a24e5c9f5227240ab215dc55



20251105_162423-assembling setup.jpg

sha256: bd54f92af902a42149ff89cae5a1789396656bb6050778b15dc606f351b0e937



20251105_162942-degassing of the suspension.jpg

sha256: ae66fc3c7d7c529d112cfb8b511c2a080513c9fa21611eeb96c985b2e5830904



20251105_172634-after start of irradiation.jpg

sha256: c67a46113f30aae96599af78af6d5d94642fcf4608ab5e3edea5cb09b21d87fd



20251105_175313-after introducing H2 sensor.jpg

sha256: 6391a3857adfcf4af26b4d1143ed78257cfefba37ea6444c3dbd018a43371742




2025-11-05_171608_NB-334-Ch2-3.txt

sha256: d413583b99a300a15a3ca4ae35577d4f8330a4d25045972e368e7583a67e8800

2025-11-05_171608_NB-334-Ch2-3.png

sha256: 72199cab35d1cb179c12d38ea6bf3bf5708764f10cc60255337557b3276ee2c1




2025-11-05_162517_NB-334-Ch2-2.txt

sha256: 34b8eb7d7938bf4d8d81fd8d0fb94e3413a2f177ca0996d9b5392048027695c0

2025-11-05_162517_NB-334-Ch2-2.png

sha256: f047a7aa684146959f586c820df13232ce80f241eda4b2e5dd93208676f024b3



2025-11-05_121425_NB-334-Ch2-1.txt

sha256: 4fc3f3cda53aafd2d6c8ce15ef2e830bbaaac6165f64d737736091f2f75d8ec8

2025-11-05_121425_NB-334-Ch2-1.png

sha256: 0c41df2c46392f3adc7722d001ac1a00ba4c8498ab9018b95a9024cc98720b9d




NB-334-Logger3-during irradiation.csv

sha256: 6b5a37ce5b71d7165800d29deae521d14b0cef238f0fbb46d7e0fcecdd524686

NB-334-Logger3-during irradiation.bmp

sha256: 9b4dd803926a03ddf78d9cdc4ae04e3f997d4fd4e630fa3702e473200d6b604e



NB-334-Logger-2point-calibration.bmp

sha256: b416bdd76aab69a890f1b6f2f06229f32e42116e967283e6a3af194d509f80e6



NB-334-Logger2-calibration.csv

sha256: 1bdb6f1136946cb6329a3426cfbb14c00da53ed2db079f17f92c7682fc79485f

NB-334-calibration-step.bmp

sha256: 520cf28a3dfa90afeee25120349593bc14ff208f6b71eff6111a0c098d747c1a




NB-334-Logger1-pre-polarization.csv

sha256: c7fc838f5d2a1af454ab2504d616d9a4151fdb9440989ed885d944f41926db97

NB-334-Logger1-pre-polarization.bmp

sha256: b538222f3a5074525c2e83129b4f47fc52e00dde697626452817974a39116872



NB-334.ulog

sha256: 2bb0e7d38032ce20ea015811ae9965a9dfbeb0e4f961bc0d7f4fbaa14eb2eb5a



Unique eLabID: 20251105-15ba5f80e68133baf651b116b6fcb886b5790ae
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3345>