

NB-333: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (NB-323, 0.2 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C (reproduction NB-328)

Date: 2025-11-03

Tags: O₂ Test Calibration Future NB
Firesting O₂ sensor H₂ SrTiO₃ Unisense
RhCrO₃:Al:SrTiO₃ H₂ Sensor
temperature In situ Trace range robust
oxygen sensor photocatalysis Unisense
normal range

Category: SrTiO₃


Status: Done

Created by: Nadzeya Brezhneva

Objectives

Reproduction of NB-328: simultaneous detection of H₂ and O₂ evolution in liquid phase for irradiated suspension of Rh_xCr_yO₃:Al:SrTiO₃ suspension (NB-323, **0.2 wt% Rh,Cr**, 0.5 mg/mL), 365 nm LED, 50 mW/cm², 20 °C.

Reaction scheme



ChemDraw file linked: [NB-328-SrTiO3-photocatalytic H2O splitting.cdxml](#)

Literature/reference experiments

Literature	/
Reproduction	SrTiO ₃ - NB-328: Liquid phase H ₂ and O ₂ of RhCrO _x ,Al:SrTiO ₃ (NB-323, 0.2 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm ² , 20 °C
Similar experiments	SrTiO ₃ - NB-316: Liquid phase H ₂ and O ₂ of RhCrO _x ,Al:SrTiO ₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm ² , 20 °C SrTiO ₃ - NB-327: Liquid phase H ₂ and O ₂ of RhCrO _x ,Al:SrTiO ₃ (NB-321, 0.05 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm ² , 20 °C

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	Concentration [mM]
milli-Q H ₂ O	/	/	/	/	/	/	/	0.998	25	/
Al:SrTiO ₃ RhCrO _x (NB-323)	SrTiO ₃ - NB-323: Modification of EA-354 (SrTiO ₃ :Al, upscaled batch) with Rh, Cr oxide cocatalyst (0.2 wt%)	/	/	/	12.50	12.68	/	/	/	/

Excel sheet for reagent calculation

/

Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-314: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup](#)

	Name
Used Set-up	Equipment - Advanced irradiation setup V1.0 I
Irradiation setup number	Equipment - Irradiation setup 4 (CEEC II, E002)

	Light Source Name	Power Source Name	Wavelength [nm]	Power Setting [mW]	Analog Setting [0.00 - 10.00]
First light source	Light Source - UHP LED 365 nm-4	Power Sources - BLS-18000-1 4	365	56	0.19

Used beam combiner [Name or None]	/
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	20
Stirring speed [rpm]	500
Irradiation start: 1. Firesting [relative to start log] 2. Unisense	1. 605 s 2. 1:37:41
Irradiation stop: 1. Firesting [relative to start log] 2. Unisense	1. 1525 s 2. 1:53:01

O₂/H₂ sensor equipment

	Equipment	Used protocol
Used Firesting	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)	Protocol - Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software
Used O ₂ sensor	Equipment - Robust probe for liquid O ₂ measurment	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor
Used H ₂ sensor	Equipment - H ₂ UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor

Procedure/observations

Date	Time	Step	Observations	Pictures/Files
		Calibration from SrTiO ₃ - NB-331: Liquid phase H ₂ and O ₂ of RhCrOx,Al:SrTiO ₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm ² , 20 °C (reproduction NB-316) I I I was used.	/	/
04.11.2025		Sample preparation		
	ca. 0:10	Weighing photocatalyst in a 50 mL vial (the vial was covered with Al foil before further use) (preliminary, the lumps of the solid were broken with Smartspatula inside the vial).	Slightly grey solid	/
	0:20	Addition of 25 mL H ₂ O to the vial via graduated cylinder.	/	/
	0:22-25	The suspension was vortexed for 3 min (Equipment - VWR® VV3, Vortex Mixer, stage 4/6), covered with Al foil before further use.	Dispersion seems to look better than in the previous experiment SrTiO ₃ - NB-328: Liquid phase H ₂ and O ₂ of RhCrOx,Al:SrTiO ₃ (NB-323, 0.2 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm ² , 20 °C	20251104_002643-suspension after vortex.jpg
		Continue in Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor from step 6		
	0:30-35	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	/

	ca. 0:37	Assembling the setup.	Currently, stopper instead of H ₂ sensor, PT100, PT1000 and O ₂ robust probe are inside the reactor immersed in the liquid phase	20251104_003621-before degassing and irradiation.jpg
	0:39	Start of O2 logging.	NB-333-Ch2-1	2025-11-04_003927_NB-333-Ch2-1.txt 2025-11-04_003927_NB-333-Ch2-1.png
	0:41	The degassing was started	/	
	1:19	Cannula was transferred to gas phase, above the suspension.	/	
	1:21	H ₂ sensor was added in Ar counterflow.	/	
	1:25	The degassing was stopped by removing the cannula and closing the valve.	/	
	1:27	Stop of O2 logging.	/	
	1:27	Start of O2 logging.	NB-333-Ch2-2	2025-11-04_012735_NB-333-Ch2-2.txt 2025-11-04_012735_NB-333-Ch2-2.png
	1:27	Start of H2 logging.	NB-333-Logger1	NB-333.ulog NB-333-Logger1.csv NB-333-Logger1.bmp
	1:27-37	Equilibration time.	/	/
	1:37	The irradiation was started	/	20251104_013911-after start of irradiation.jpg
	1:53	The irradiation was stopped.	/	/
	1:53-2:03	Equilibration time.	/	/
	2:03	Stop of O2 and H2 logging.	/	/

	2:10	Deassembling the setup, cleaning the reactor.	Tips of the sensors and reactor were covered with attached photocatalyst particles. Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode).	20251104_021611-after irradiation.jpg
--	------	---	--	---

Analysis

Used calibration for Firesting: [20250611-BOLA-fitting-liquid-phase-trace-oxygen-sensor-H2-O2 reactor.ini](#)

Used calibration for UniSense: NB-331-Logger3

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
04.11.2025	1:27	NB-333-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-333.ulog NB-333-Logger1.csv	NB-333-O2 and H2 curve.py	NB-333-Logger1.bmp NB-333-O2 and H2 curve.png	SrTiO3 - NB-328: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (NB-323, 0.2 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C	Clean H2 response signal, H2 value ca. 29 uM
	0:39	NB-333-Ch2-1	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-04_003927_NB-333-Ch2-1.txt	/	2025-11-04_003927_NB-333-Ch2-1.png	/	Degassing of the suspension.
	1:27	NB-333-Ch2-2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-04_012735_NB-333-Ch2-2.txt	NB-333-O2 and H2 curve.py	2025-11-04_012735_NB-333-Ch2-2.png NB-333-O2 and H2 curve.png	SrTiO3 - NB-328: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (NB-323, 0.2 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C	Sharp increase/decrease in O2 value after start and stop of irradiation. O2 value ca. 29 uM.

Results

Reproduction of NB-328, simultaneous H₂ and O₂ measurements of irradiated suspension of NB-323 (0.5 mg/mL) in O₂/H₂ photoreactor under 365 nm irradiation (50 mW/cm², 20 °C, 15 min) were performed.

H₂ level at the end of irradiation - 29 umol/L, O₂ level - 29 umol/L (but sharp change in signal after start and finish of irradiation in O2 curve - needs to be eliminated).

Future recommendations

Old procedure	Problem	Suggested new procedure
/	Sharp increases in O ₂ signal - maybe related either with the position of the sensors in the reactor.	Adjust the position of the sensors relatively to the light beam.

Linked experiments

SrTiO₃ - NB-316: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

SrTiO₃ - NB-318: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm², 20 °C

SrTiO₃ - NB-322: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm², 20 °C (reproduction NB-318)

SrTiO₃ - NB-323: Modification of EA-354 (SrTiO₃:Al, upscaled batch) with Rh, Cr oxide cocatalyst (0.2 wt%)

SrTiO₃ - NB-325: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm², 20 °C

SrTiO₃ - NB-326: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm², 20 °C (reproduction NB-325)

SrTiO₃ - NB-327: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (NB-321, 0.05 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

SrTiO₃ - NB-328: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (NB-323, 0.2 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

SrTiO₃ - NB-331: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C (reproduction NB-316) I I I

Linked resources

Equipment - [VWR® VV3, Vortex Mixer](#)

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)

Protocol - [Hydrogen Measurement Using GC](#)

Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)

Protocol - [In-situ hydrogen and oxygen measurment in H2/O2 reactor](#)


Attached files

NB-333-O2 and H2 curve.py

sha256: 2a9a7e0682be7d43f1eba3e624d06391babd220065b66d5fc6e8dd2c0978713e

NB-333-O2 and H2 curve.png

sha256: 2386802e5e0c0e8e857ddfe04c6aaec0f5676f2cba5c7e4c7f5778693ba71808



20251104_002643-suspension after vortex.jpg

sha256: 84e1ed5e7b7a3a645ce1eda84cb3e8fa962b676b0420675136dffdcf368786e0



20251104_021611-after irradiation.jpg

sha256: 2598b411351b79f29a88c4aa89680faa17199d88a68bc1952d1db80b4f78b1b2



20251104_004141-during degassing.jpg

sha256: 64a27b3a5202f98639bbfdbc0aa2ab00d280b0e7146df6b60b35b2f4ded68b7f



20251104_003621-before degassing and irradiation.jpg

sha256: 10c507128fc4e7ba08f4f1855a72664eed17cac9ebb2d852a5a5e8b3e599a795



20251104_013911-after start of irradiation.jpg

sha256: 40711856d718a432919645cfcc012ca003867d80df72a03f6ad0b600400bdad1



2025-11-04_012735_NB-333-Ch2-2.png

sha256: aa701e3ebf502d5e784ce0da5c94c1658c648613092560666e9127af8155d4f0



2025-11-04_012735_NB-333-Ch2-2.txt

sha256: 07a165852f6508ffdb8d4deeb6952202294d3d0735e152b69b5e6094d585511

2025-11-04_003927_NB-333-Ch2-1.png

sha256: 5eb4b9dd918578eef63b35dda273a55f25b0fc758949412553a41a4f0bb29e25



2025-11-04_003927_NB-333-Ch2-1.txt


sha256: c05d29c2daf21d5c3751276c1d3710a4e9268bd2500605d64bdd246ba311e612

NB-333-Logger1.csv

sha256: 1cd6ca2f344a2e3890733f7b0b753cbc105319577117e827d7b9427e735a7881

NB-333-Logger1.bmp

sha256: be81f2224b3e91208d250f6e960df2818c7e1e346f6fcbe931accf166fc5f511



NB-333.ulog

sha256: 26193abef77dadfc6e0622aeb955569c0db1826cf69531389d898c990742fca0



Unique eLabID: 20251103-f1add9728c1eea04ca4a757345d2d612a88750f1
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3325>