


# AE-594: Design and printing of 2nd generation of holder for power measurement

Date: 2025-10-16  
Tags: AE CAD Autodesk Inventor Additive manufacturing 3D print advanced irrad setup Prusa Slicer  
Category: Photoreactor  
Status: Done  
Created by: Alexander Eith

## CAD and Image of 3D Print

CAD	Photo
	

## Literature/reference experiments

Literature	/
Reproduction	/
Similar experiments	Photoreactor - AE-431: Print and Setup of advanced power measurement setup V1.0

## Filament

Name	CAS Number / Experiment Number	Mass [g]
PLA (PrimaValue PLA+ black, 1.75 mm diameter)	26680-10-4	72

# 3D printing parameters

Used printer	Equipment - 3D printer - PRUSA MK4S - CEEC I Lab 208
Printing temperature [°C]	215
Bed temperature [°C]	60
Layer height [mm]	0.2
Infill (%)	15
Infill type	"Gitternetz"
Support structure type	organic

## Printing procedure

Part identifier	Date	Time	CAD file	STL file	Description	Reference part	Modifications	Print Result	Future modifications
power_sesnor_aperture_V1.0	14.10	13:00	<a href="#">power_sensor_aperture_V1.0.ipt</a>	<a href="#">power_sensor_aperture_V1.0.stl</a>	Aperture for holding power sensor	/	/	worked	Sensor does not fit into holder, make it a bit wider
power_sensor_aperture_V1.1	22.10	14:00	<a href="#">power_sensor_aperture_V1.1.ipt</a>	<a href="#">power_sensor_aperture_V1.1.stl</a>	Aperture for holding power sensor	power_sensor_aperture_v1.0	Holder 1 mm shorter small sides and 3 mm longer longer sides	worked	fits well

## Results

power\_sensor\_aperture\_V1.1 can and should be used with [Equipment - Advanced power measurement setup V1.0](#) |

## Linked experiment

Photoreactor - [AE-431: Print and Setup of advanced power measurement setup V1.0](#)

## Linked resources

Equipment - [Advanced irradiation setup V1.0](#) |

Equipment - [Advanced power measurement setup V1.0](#) |

Equipment - [3D printer - PRUSA MK4S - CEEC I Lab 208](#)

## Attached files

power\_sensor\_aperture\_V1.1.stl

sha256: 22cae2f452eccba690c6ba73ef4114235968e1e063f5eac81d67d10c5977a292

power\_sensor\_aperture\_V1.1.ckpt

sha256: df99e8654d9818574c8b9a58be563ecb9ae4a0486022d473fddd987bb397375f

power\_sensor\_aperture\_V1.0.stl


sha256: 47509540277b39f076f56935aeb1a1d7c607342d61569ba5e09c6d49c6cadb73

power\_sensor\_aperture\_V1.0.ckpt

sha256: 407299f36a5a335d09c2e51a8a935df09c08e714832fd2cf36e8cc4a08cc1795

CAD file.png

sha256: 9008a75b4ba04b2682fbf7d7602aa02a29decd15fefe57821cc5b4065041e593



20251016\_100018.jpg

sha256: c47e4526a62d26e58f1796f37612066fb4e5215bf4783205f891621f6947ad39



CAD file.png

sha256: 9008a75b4ba04b2682fbf7d7602aa02a29decd15fefe57821cc5b4065041e593



picture.jpg

sha256: c47e4526a62d26e58f1796f37612066fb4e5215bf4783205f891621f6947ad39



Unique eLabID: 20251016-22b68d3e8311fc48104b42d304f4761da352a5ad

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3202>