

NB-340: Preparation of $\text{RhCl}_3 \cdot 3\text{H}_2\text{O}$ and $\text{Cr}(\text{NO}_3)_3 \cdot 9\text{H}_2\text{O}$ stock solutions

Date: 2025-11-11

Tags: Stocksolution KRA NB
 $\text{RhCl}_3 \cdot 3\text{H}_2\text{O}$ $\text{Cr}(\text{NO}_3)_3 \cdot 9\text{H}_2\text{O}$ Solution

Category: Prep work

Status: Done

Created by: Nadzeya Brezhneva

Objective

Preparation of stock solutions of RhCl_3 and $\text{Cr}(\text{NO}_3)_3$ salts.

Literature/reference experiments

Literature	https://doi.org/10.1039/D4SC03978E
Reproduction	/
Similar experiments	Prep work - EA-373: Preparation of RhCl_3 and $\text{Cr}(\text{NO}_3)_3$ stock solutions Prep work - NB-275: Preparation of RhCl_3 and $\text{Cr}(\text{NO}_3)_3$ solutions I Prep work - EA-MEJ-357: Preparation of RhCl_3 and $\text{Cr}(\text{NO}_3)_3$ stock solutions

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume _{theo} [ml]	Volume _{exp} [ml]	Concentration [mM]
$\text{RhCl}_3 \cdot 3\text{H}_2\text{O}$, BLD Pharmatech GmbH	13569-65-8	C138335	0.0189	/	5	5	263.31	/	/	/	9.87
milli-Q water for $\text{RhCl}_3 \cdot 3\text{H}_2\text{O}$ solution	/	/	/	/	/	/	/	/	1.923	1.923	/
$\text{Cr}(\text{NO}_3)_3 \cdot 9\text{H}_2\text{O}$	7789-02-8	/	0.0509	/	20	20.38	400.15	/	/	/	28.81
milli-Q water for $\text{Cr}(\text{NO}_3)_3 \cdot 9\text{H}_2\text{O}$ solution	/	/	/	/	/	/	/	/	1.767	1.767	/

Excel sheet for reagent calculation

[NB-340- stock solutions calculations.xlsx](#)

Glovebox parameters

Glovebox used	Equipment - Glovebox CEEC I K004
Date	11.11.2025
Time	9:50-10:17
O2 value before use	93.0 ppm

H2O value before use	<0.1 ppm
Picture of log book page	20251111_101740-glovebox-logbook.jpg

Procedure/observations

For transfer of precise liquid amount, Eppendorf pipettes were used (above 100 ul: 1000 ul Eppendorf).

Date	Time	Step	Observations	Pictures
11.11.2025	ca. 10:00	Weighing $\text{RhCl}_3 \cdot 3\text{H}_2\text{O}$ inside the glovebox Equipment - Glovebox CEEC I K004 , transfer it to a 10 ml vial (performed by KRA).	Dark red crystals	/
	10:40	Weighing $\text{Cr}(\text{NO}_3)_3 \cdot 9\text{H}_2\text{O}$ in a 10 ml vial.	Dark blue crystals	/
	10:55	The calculated amount of water to reach the specific concentration was added to the vial with $\text{RhCl}_3 \cdot 3\text{H}_2\text{O}$ crystals.	Orange-red solution NB-340-RhCl₃ solution	/
	11:00	The calculated amount of water to reach the specific concentration was added to the vial with $\text{Cr}(\text{NO}_3)_3 \cdot 9\text{H}_2\text{O}$.	Blue solution NB-340-Cr(NO₃)₃ solution	/
	ca. 11:05	Solution of NB-340-RhCl₃ was transferred into a 4ml screw cap vial.	Red solution NB-340-RhCl₃ 9.87 mM 11.11.2025	20251111_110630-final solutions.jpg
	ca. 11:10	Solution of NB-340-Cr(NO₃)₃ vial was transferred into a 4ml screw cap vial using an Eppendorf pipette.	Blue solution NB-340-Cr(NO₃)₃ 28.81 mM 11.11.2025	20251111_110630-final solutions.jpg
	ca. 11:15	The vials were covered with Al foil and stored under fume hood.	/	/

Product characterization

Sample	Concentration mM	Mass [mg]	Purity	Volume, mL	Amount [μmol]	Yield [%]	Description	Image	Storage location
NB-340-RhCl ₃ 11.11.2025	9.87	/	/	1.923	/	/	Red solution	20251111_110630-final solutions.jpg	Laboratory - Lab E004 - CEEC II, NB fume hood
NB-340-Cr(NO ₃) ₃ 11.11.2025	28.81	/	/	1.767	/	/	Blue solution	20251111_110630-final solutions.jpg	Laboratory - Lab E004 - CEEC II, NB fume hood

Results

Fresh stock solution of RhCl₃ (9.87 mM in H₂O) and Cr(NO₃)₃ (28.81 mM in H₂O) salts were prepared.

Linked experiments

Prep work - [NB-275: Preparation of RhCl₃ and Cr\(NO₃\)₃ solutions I](#)

Prep work - [EA-373: Preparation of RhCl₃ and Cr\(NO₃\)₃ stock solutions](#)

Linked resources

Equipment - [Glovebox CEEC I K004](#)

Laboratory - [Lab E004 - CEEC II](#)

Attached files

NB-340- stock solutions calculations.xlsx

sha256: e06e911254ebc01e51438317fe28176cc61c42317ee4eef86d134585f7928e09

20251111_101740-glovebox-logbook.jpg

sha256: 77201e60bd7c95cd87fc4a987af338a985835351ae27d16841e3fe766ad0a2ae



20251111_110630-final solutions.jpg

sha256: 753ff99278bec0877377a9aae8fcfd80c3d9dbf24089d4f14b483be9196b9a7a





Unique eLabID: 20251111-db383ff0fd11deee2542aa280f5ece9849fd7620

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3412>