

# NB-320: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrO<sub>x</sub>,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 20 °C (reproduction NB-316) I

Date: 2025-10-27

Tags: O<sub>2</sub> Test Calibration Future NB  
Firing O<sub>2</sub> sensor H<sub>2</sub> SrTiO<sub>3</sub>  
troubleshooting Unisense  
RhCrO<sub>3</sub>:Al:SrTiO<sub>3</sub> H<sub>2</sub> Sensor  
temperature In situ Trace range robust  
oxygen sensor photocatalysis Unisense  
normal range

Category: SrTiO<sub>3</sub>


Status: Done

Created by: Nadzeya Brezhneva

## Objectives

Second reproduction of NB-316 (reference experiment): simultaneous detection of H<sub>2</sub> and O<sub>2</sub> evolution in liquid phase for irradiated suspension of Rh<sub>x</sub>CrO<sub>y</sub>:Al:SrTiO<sub>3</sub> suspension (EA-358, 0.5 mg/mL), 365 nm LED, 50 mW/cm<sup>2</sup>, 20 °C.

## Reaction scheme



ChemDraw file linked: [NB-320-SrTiO3-photocatalytic H2O splitting.cdxml](#)

## Literature/reference experiments

Literature	/
Reproduction	SrTiO <sub>3</sub> - NB-316: Liquid phase H <sub>2</sub> and O <sub>2</sub> of RhCrO <sub>x</sub> ,Al:SrTiO <sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm <sup>2</sup> , 20 °C SrTiO <sub>3</sub> - NB-319: Liquid phase H <sub>2</sub> and O <sub>2</sub> of RhCrO <sub>x</sub> ,Al:SrTiO <sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm <sup>2</sup> , 20 °C (reproduction NB-316)
Similar experiments	SrTiO <sub>3</sub> - NB-315: Liquid phase H <sub>2</sub> and O <sub>2</sub> of RhCrO <sub>x</sub> ,Al:SrTiO <sub>3</sub> (NB-301, 0.5 mg/mL), 365 nm, 50 mW/cm <sup>2</sup> , 20 °C

## Reagents

Name	CAS Number / Experiment Number	Inventor y number	Amount [mmol]	Equivalen ts	Mass <sub>theo</sub> [mg]	Mass <sub>exp</sub> [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	Pressure [bar]	Concentrati on [mM]
milli-Q H <sub>2</sub> O	/	/	/	/	/	/	/	0.998	25 + 25 (for calibration)	/	/

Al:SrTiO3 RhCrOx (EA-358)	SrTiO3 - EA-358: Modification of Al:SrTiO3 (EA-354) via deposition of Rh, Cr oxide co-catalyst, 350°C, 1h, Upscaling (3.33x)	/	/	/	12.50	12.80	/	/	/	/	/
Hydrogen	1333-74-0	/	/	/	/	/	/	/	2 balloons (ca. 2 L)	approx. 1	/

## Excel sheet for reagent calculation

/

## Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-314: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup](#)

	Name
Used Set-up	<a href="#">Equipment - Advanced irradiation setup V1.0 I</a>
Irradiation setup number	<a href="#">Equipment - Irradiation setup 4 (CEEC II, E002)</a>

	Light Source Name	Power Source Name	Wavelength [nm]	Power Setting [mW]	Analog Setting [0.00 - 10.00]
First light source	<a href="#">Light Source - UHP LED 365 nm-4</a>	<a href="#">Power Sources - BLS-18000-1 4</a>	365	56	0.19

Used beam combiner [Name or None]	/
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	20

<b>Stirring speed [rpm]</b>	500
<b>Irradiation start:</b> <b>1. Firesting [relative to start log]</b> <b>2. Unisense</b>	1. 600 s 2. 17:29:20
<b>Irradiation stop:</b> <b>1. Firesting [relative to start log]</b> <b>2. Unisense</b>	1. 1513 s 2. 17:44:33

## O<sub>2</sub>/H<sub>2</sub> sensor equipment

	Equipment	Used protocol
Used Firesting	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)	Protocol - Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software
Used O <sub>2</sub> sensor	Equipment - Robust probe for liquid O <sub>2</sub> measurment	Protocol - In-situ hydrogen and oxygen measurment in H <sub>2</sub> /O <sub>2</sub> reactor
Used H <sub>2</sub> sensor	Equipment - H <sub>2</sub> UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurment in H <sub>2</sub> /O <sub>2</sub> reactor

## Procedure/observations

Date	Time	Step	Observations	Pictures/Files
27.10.2025		The experiment was done according to <a href="#">Protocol - In-situ hydrogen and oxygen measurment in H<sub>2</sub>/O<sub>2</sub> reactor</a> Important steps and deviations are listed below	/	/
	9:27-10:27	Conditioning of H <sub>2</sub> sensor	<b>NB-320-Logger1</b> ca. 1.60 mV at the end of polarization procedure	<a href="#">NB-320.ulog</a> <a href="#">NB-320-Logger1.csv</a> <a href="#">NB-320-Logger1.bmp</a>
		Assembling the setup for calibration (25 mL of water was added using graduated cylinder), done according to <a href="#">Protocol - Liquid phase calibration of H<sub>2</sub> UniAmp sensor with H<sub>2</sub> bubbling</a> .		

	11:11	Start of O <sub>2</sub> logging.	<b>NB-320-Ch2-1</b>	<a href="#">2025-10-27_111158_NB-320-Ch2-1.txt</a> <a href="#">2025-10-27_111158_NB-320-Ch2-1.png</a>
	11:12	Start of H <sub>2</sub> logging.	<b>NB-320-Logger2</b> -2 mV offset	<a href="#">NB-320.ulog</a> <a href="#">NB-320-Logger2.bmp</a>
		Important note: BOLA fitting for PT1000 sensor did not fit well, so the measurement was stopped and the setup was deassembled to eliminate the problem.		
	11:21	Start of O <sub>2</sub> logging.	<b>NB-320-Ch2-2</b>	<a href="#">2025-10-27_112151_NB-320-Ch2-2.txt</a> <a href="#">2025-10-27_112151_NB-320-Ch2-2.png</a>
	11:22	Start of H <sub>2</sub> logging.	<b>NB-320-Logger3</b>	<a href="#">NB-320.ulog</a> <a href="#">NB-320-Logger3-not saved calibration.csv</a> <a href="#">NB-320-Logger3-calibration curve step.bmp</a>
	11:22	Degassing was started.	/	/
	ca. 12:03	Removal cannula from water, keeping it in the gas phase above the liquid.	/	/
	12:05	Introducing H <sub>2</sub> sensor into the reactor under Ar flow.	/	/
	12:15	H <sub>2</sub> bubbling of the reactor was started	/	/
	12:21	Removing cannula from the reactor, closing the valve.	/	<a href="#">20251027_122208-first attempt of calibration.jpg</a>
	12:31	Introducing second H <sub>2</sub> balloon.	/	/
	ca. 12:40	1.000.000 ppm point was taken	853 mV, 775 $\mu$ M	/
		Important note: But the 0 point was not taken, calibration needs to be redone.		
	ca. 12:50-13:00	Deassembling the setup, cleaning with acetone and compressed air.	/	/
	ca. 13:10-15	Assembling the setup for calibration once again.	/	/

	13:19	Start of O <sub>2</sub> logging.	NB-320-Ch2-3	2025-10-27_131923_NB-320-Ch2-3.txt 2025-10-27_131923_NB-320-Ch2-3.png
	13:20	Start of H <sub>2</sub> logging.	NB-320-Logger4 (used calibration)	NB-320.ulong NB-320-Logger4-used calibration.csv NB-320-Logger4-calibration info.png NB-320-Logger4-calibration curve step.bmp NB-320-Logger4-calibration 2 points.bmp
	13:24	Start degassing.	/	20251027_132610-degassing of water.jpg
	13:51	Moving cannula to the gas phase above the liquid.	/	/
	13:53	Introducing H <sub>2</sub> sensor under Ar flow.	/	/
	13:56	0 point was taken.	0.4 mV	/
	14:02	Introducing H <sub>2</sub> balloon.	/	20251027_140423-introducing H2 balloon.jpg
	14:10	Introducing 2nd H <sub>2</sub> balloon.	/	20251027_141044-2nd balloon.jpg
	14:21	1.000.000 ppm point was taken and calibration was saved.	814 mV, 775 uM, slope 1.049	20251027_142159-H2 table.jpg
	14:23	Stop of O <sub>2</sub> and H <sub>2</sub> loggings.	/	/
	ca. 14:30	Deassembling the setup, drying the reactor with acetone and compressed air .	/	/
		<b>Sample preparation</b>	/	/
	14:55	Weighing <b>EA-358</b> photocatalyst in a 50 mL vial.	Creamy solid.	/
	15:00	Addition of 25 mL H <sub>2</sub> O to the vial via graduated cylinder.	/	/
	15:07-10	The suspension was vortexed for 3 min ( Equipment - VWR® VV3, Vortex Mixer, stage 4/6), covered with Al foil before further use.	/	20251027_151031-suspension after vortex.jpg
		Continue in Protocol - In-situ hydrogen and oxygen measurment in H2/O2 reactor from step 6		

	ca. 15:20	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	/
	15:30-50	Eliminating leakage in the water cooling system.	Tubing was loosened and not completely sealed to the olive, leakage of water appeared when LAUDA was on --> problem was eliminated.	/
	ca. 15:55	Assembling the setup.	Currently, stopper instead of H <sub>2</sub> sensor, PT100, PT1000 and O <sub>2</sub> robust probe are inside the reactor immersed in the liquid phase	/
	15:57	The degassing was started	/	<a href="#">20251027_160356-degassing of the suspension.jpg</a>
	16:05	<a href="#">Start of O<sub>2</sub> logging.</a>	<b>NB-320-Ch2-4</b>	<a href="#">2025-10-27_160525_NB-320-Ch2-4.txt</a> <a href="#">2025-10-27_160525_NB-320-Ch2-4.png</a>
	16:26	Cannula was transferred to gas phase, above the suspension.	/	/
	16:28	H <sub>2</sub> sensor was added in Ar counterflow.	/	<a href="#">20251027_163001-introducing H2 sensor under Ar flow.jpg</a>
	16:31	The degassing was stopped by removing the cannula and closing the valve.	/	/
	16:35	<a href="#">Stop of O<sub>2</sub> logging.</a>	/	/
	16:35	<a href="#">Start of O<sub>2</sub> logging.</a>	<b>NB-320-Ch2-5</b>	<a href="#">2025-10-27_163556_NB-320-Ch2-5.txt</a> <a href="#">2025-10-27_163556_NB-320-Ch2-5.png</a>
	16.35	<a href="#">Start of H<sub>2</sub> logging.</a>	<b>NB-320-Logger5</b>	<a href="#">NB-320.ulog</a>
	16:35-45	Equilibration time.	16:39 - putting side cover in the box --> change in O <sub>2</sub> value, the stirring was stopped	/
	16:47	Switching stirring of the suspension ON again.	/	/

	16:51	Stop of O <sub>2</sub> and H <sub>2</sub> loggings.	Too high O <sub>2</sub> leakage, BOLA fittings of all sensors need to be checked.	/
	17:00	Start O <sub>2</sub> logging.	<b>NB-320-Ch2-6</b>	2025-10-27_170003_NB-320-Ch2-6.txt 2025-10-27_170003_NB-320-Ch2-6.png
	17:01	Start degassing once again.	/	/
	17:12	Removing cannula from the liquid phase, placing above the suspension in the gas phase.	/	/
	17:13	Introducing H <sub>2</sub> sensor under Ar.	/	/
	17:16	Stop degassing, closing the valve.	/	20251027_171842-after introducing H2 sensor to the suspension.jpg
	17:18	Stop O <sub>2</sub> logging.	/	/
	17:19	Start H <sub>2</sub> logging.	<b>NB-320-Logger6 (H<sub>2</sub> logging for irradiation)</b>	NB-320.ulog NB-320-Logger6-during irradiation.csv NB-320-Logger6-during irradiation.bmp
	17:19	Start O <sub>2</sub> logging.	<b>NB-320-Ch2-7 (O<sub>2</sub> logging for irradiation)</b>	2025-10-27_171920_NB-320-Ch2-7.txt 2025-10-27_171920_NB-320-Ch2-7.png
	17:19-29	Equilibration time.	Leakage rate ca. 0.5 um in 10 min	/
	17:29	The irradiation was started	Step in H <sub>2</sub> curve at ca. 17:42	20251027_172931-after start of irradiation.jpg
	17:44	The irradiation was stopped.	/	/
	17:44-54	Equilibration time.	/	/
	17:54	Stop of O <sub>2</sub> and H <sub>2</sub> logging.	/	/

	ca. 18:00	Deassembling the setup, cleaning the reactor.	Tips of the sensors and reactor were covered with attached photocatalyst particles. Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode).	/
--	-----------	---	--	---

## Analysis

Used calibration for Firesting: [20250611-BOLA-fitting-liquid-phase-trace-oxygen-sensor-H2-O2 reactor.ini](#)

Used calibration for UniSense: NB-320-Logger4

(for feasibility, the succesful O<sub>2</sub> and H<sub>2</sub> loggings used for calibration and photocatalytic test are highlighted with green)

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
27.10.2025	9:27-10:27	NB-320-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-320.ulong NB-320-Logger1.csv	/	NB-320-Logger1.bmp	/	Pre-polarization
	11:12	NB-320-Logger2	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-320.ulong	/	NB-320-Logger2.bmp	/	Unsuccessful.
	11:22	NB-320-Logger3	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-320.ulong NB-320-Logger3-not saved calibration.csv	/	NB-320-Logger3-calibration curve step.bmp	/	0 point was not taken
	13:20	NB-320-Logger4	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-320.ulong NB-320-Logger4-used calibration.csv	/	NB-320-Logger4-calibration info.png NB-320-Logger4-calibration curve step.bmp NB-320-Logger4-calibration 2 points.bmp	SrTiO3 - NB-316: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C	<b>Used calibration for H2 sensor</b> 10 <sup>-6</sup> ppm - 814 mV, slope 1.049
	16:35	NB-320-Logger5	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-320.ulong	/	/	/	Unsuccessful.



	17:19	NB-320-Logger6	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-320.ulong NB-320-Logger6-during irradiation.csv	NB-320-O2 and H2 curve.py	NB-320-Logger6-during irradiation.bmp NB-320-O2 and H2 curves.png	/	<b>H2 logging for the photocatalytic test</b> Clean response of H2 sensor, step in signal at ca. 17.42, ca. 47 uM H2 at the end of irradiation.
	11:11	NB-320-Ch2-1	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-10-27_111158_NB-320-Ch2-1.txt	/	2025-10-27_111158_NB-320-Ch2-1.png	/	First attempt, BOLA for PT1000 was not tight, needed to be redone.
	11:21	NB-320-Ch2-2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-10-27_112151_NB-320-Ch2-2.txt	/	2025-10-27_112151_NB-320-Ch2-2.png	/	First attempt of degassing in the first attempt of calibration.
	13:19	NB-320-Ch2-3	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-10-27_131923_NB-320-Ch2-3.txt	/	2025-10-27_131923_NB-320-Ch2-3.png	/	Degassing for the calibration of H2 sensor used further in the experiment.
	16:05	NB-320-Ch2-4	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-10-27_160525_NB-320-Ch2-4.txt	/	2025-10-27_160525_NB-320-Ch2-4.png	/	First degassing of the photocatalyst suspension.
	16:35	NB-320-Ch2-5	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-10-27_163556_NB-320-Ch2-5.txt	/	2025-10-27_163556_NB-320-Ch2-5.png	/	Change in O2 value during stirring and after stop of stirring. High O2 leakage rate afterwards, elimination of the problem was required.
	17:00	NB-320-Ch2-6	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-10-27_170003_NB-320-Ch2-6.txt	/	2025-10-27_170003_NB-320-Ch2-6.png	/	2nd degassing of the photocatalyst suspension after eliminating leakage problem.
	17:19	NB-320-Ch2-7	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-10-27_171920_NB-320-Ch2-7.txt	NB-320-O2 and H2 curve.py	2025-10-27_171920_NB-320-Ch2-7.png NB-320-O2 and H2 curves.png	SrTiO3 - NB-316: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C	<b>O2 logging for the photocatalytic test.</b> Increase in O2 value during irradiation, clean response, O2 value ca. 37 uM at the end of irradiation

## Results

Simultaneous H<sub>2</sub> and O<sub>2</sub> measurements of irradiated suspension of EA-358 (0.5 mg/mL) in O<sub>2</sub>/H<sub>2</sub> photoreactor under 365 nm irradiation (50 mW/cm<sup>2</sup>, 20 °C, 15 min, reproduction of [SrTiO3 - NB-316: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 \(EA-358, 0.5 mg/mL\), 365 nm, 50 mW/cm2, 20 °C](#)) were performed.

H<sub>2</sub> level at the end of irradiation - ca. 47 umol/L, O<sub>2</sub> level - ca. 37 umol/L.

## Future recommendations

Old procedure	Problem	Suggested new procedure
/	High leakage rate in O <sub>2</sub> measurements.	Check BOLA fitting of all sensors (should be tight and fit well to all sensors), check how well the GL14/NS14 and GL18/NS14 adapters are greased.

/	0 point during H <sub>2</sub> calibration was not saved.	Don't forget to take 0 point after introducing H <sub>2</sub> sensor in the reactor and before introducing H <sub>2</sub> balloon.
/	Stop stirring of the stirring plate when closing the irradiation box with the front cover.	Check the position of the stirring plate in the box. When closing the box before irradiation, no touching of the wall with the press button of the stirrer should appear.

## Linked experiments

SrTiO<sub>3</sub> - NB-316: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrO<sub>x</sub>,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 20 °C

SrTiO<sub>3</sub> - NB-319: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrO<sub>x</sub>,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 20 °C (reproduction NB-316)

## Linked resources

Equipment - [VWR® VV3, Vortex Mixer](#)

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)

Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)


Protocol - [Liquid phase calibration of H<sub>2</sub> UniAmp sensor](#)

Protocol - [In-situ hydrogen and oxygen measurment in H<sub>2</sub>/O<sub>2</sub> reactor](#)

## Attached files

NB-320-O2 and H<sub>2</sub> curves.png

sha256: d9d7b2c215feeb0f98b8a672f92447f6783bcc820b5ba6441c70159a250157ac



NB-320-O2 and H<sub>2</sub> curve.py

sha256: 06a3d6350deea81473a72eaf51a1df31fd3c06d91e50eb2d9e4d029d332ae65a

Unisense-NB-320Screenshot 2025-11-07 093103.png

sha256: 6d6792790fa62490978f7a55483553df4ba3a286d02c4540fc8b265f647c02f1



NB-320-SrTiO<sub>3</sub>-photocatalytic H<sub>2</sub>O splitting.cdxml

sha256: dc6941c4ba571efe8a190d085ebadee0c5fc432b02f7b226bc107c104191d748

NB-320-SrTiO<sub>3</sub>-photocatalytic H<sub>2</sub>O splitting.png

sha256: f9541119d3b1697940988dad6b96e3acc271559fdf00a3ad36913aa8e3309d61



20251027\_175519-after irradiation.jpg

sha256: 8595db9f876d723b4de4b807102b30ff7c6e2fa5f3b2c6bd29be08a04aa67697



20251027\_172931-after start of irradiation.jpg

sha256: b537ec0799cde01912c25db2c5cef9f57ce4ab2c14b7b17cd39ac2a163c3a9a1



20251027\_160356-degassing of the suspension.jpg

sha256: cc21a56fc2f6ac650faf975c15780bac561fd12f3bf3afe86bcb09c80709332f



20251027\_163001-introducing H<sub>2</sub> sensor under Ar flow.jpg

sha256: 5b9fbfcc391289426327e6f9dd8855fe7abb256303007663f9b6fc07a660bdac



20251027\_171842-after introducing H<sub>2</sub> sensor to the suspension.jpg

sha256: 0de6ea3f880e20c77f05f61fcac095fc8abf43843f9e2fb6b25f51b8fc1bbabf



20251027\_142159-H2 table.jpg

sha256: 71a6a691950adc520713ff0702be0ed5683aad2b0734ac70208c635179981041




20251027\_151031-suspension after vortex.jpg

sha256: 03268bacae5e34d8d440fa479690e0caba0a0736850c3587b817f530f577f194



20251027\_141044-2nd balloon.jpg

sha256: 4bc06aa0af35ae49a8ad1f715a2c63ea0c594deaea3c851c2f19481a1cf8cb5b



20251027\_140423-introducing H2 balloon.jpg

sha256: 553d6b21dbc3ff780c45215931aeefb2b417c034572c88ec4f34937262db389b



20251027\_132610-degassing of water.jpg

sha256: 0f4ecec37e600b8cb28638430ecbd2a5f29ed2f0f29279398a551419015bfd43



20251027\_122208-first attempt of calibration.jpg

sha256: 39996abfc64f186d0d333771c9c0588b2e51f053ce5cbf84c2cf943942858578




NB-320-Logger6-during irradiation.csv

sha256: dfff876c2397631123dd25914a0b1a65e50f03767c36e3e1e106167823bf0eab


NB-320-Logger6-during irradiation.bmp

sha256: 0160b7dd28e6d9dd26938369f7cdd5152ac8e1ea54633345ec7dfddbd324aeda



NB-320-Logger4-calibration curve step.bmp

sha256: 02289e895b3aef33da6029edb29eb7a44275284079da529448467cbdaa048faf



NB-320-Logger4-used calibration.csv

sha256: d23643e73464ba5f5bf39a71128305e9dc6eff310985cfab500ae4329c9b28d4

NB-320-Logger4-calibration info.png

sha256: 483c0e9861a00700ab47e59bbcd1c190147f428b0536857deb46b058fbe6ba99



NB-320-Logger4-calibration 2 points.bmp

sha256: 4325e6859d28ee8e7410bafeb804493481b2c63e25452b80809b9f0964ec4d20




NB-320-Logger3-not saved calibration.csv

sha256: 7888e7d555eb50e87c83bdee06cdb8a8de57303a3ca0fb9bf739dcf7ebdf080a

NB-320-Logger3-calibration curve step.bmp

sha256: 78237887f94bc9cb64655fe06ab27012b56839bcf7747e7edc11841aa2c6658b



NB-320-Logger2.bmp

sha256: 0307b5398c79198ca5119fe370de242ccb1ade4dbe77baee8d2daaaee2264145




NB-320-Logger1.csv

sha256: 214dbe359106663ba23fdf50ed46aa5e5a5dfb6ac9e9d35c74aeeed806ba9aae

NB-320-Logger1.bmp

sha256: 2e7ada12ee7a0a341dd28286803e125d5eed2625475952bc95dcc5a819defd26




2025-10-27\_171920\_NB-320-Ch2-7.txt

sha256: de96e1bdd9f36bcaf5eca4c49ba40aa6062977676d7c214a09215762ab51512c


2025-10-27\_171920\_NB-320-Ch2-7.png

sha256: 825bc9eaaca2094687363472f429a5cb5f8d24e3745890c91243a3637d6a1cdd



2025-10-27\_170003\_NB-320-Ch2-6.png

sha256: e33a1537d39ca1777655632a5e4584ff2c0dc3aaef4d61003fbd289765b81243



2025-10-27\_170003\_NB-320-Ch2-6.txt


sha256: 6fd7b316bee332224c9235388c540c3cf3191ce135ea3fdd0319840f10ae0363

2025-10-27\_163556\_NB-320-Ch2-5.txt

sha256: eb020b90e9b8c2f57c4ca901a62942acf0634e7ae654b80f726dd8966496697b

2025-10-27\_163556\_NB-320-Ch2-5.png

sha256: 36c078fa47a5d32950c138326dd5b5b9a85a5e620885b5cdadd1aa897bfbfa19



2025-10-27\_160525\_NB-320-Ch2-4.txt

sha256: f9c2719b2c3df350f311dbe63dd814f06473208ee5d3c0bd099c7e8223f12179

2025-10-27\_160525\_NB-320-Ch2-4.png

sha256: 362144ada02759497310739d59ab4090946b58e6c368ef826ac125686a9a79fb



2025-10-27\_131923\_NB-320-Ch2-3.txt

sha256: d5ea65ed88a765c2c7f0f759d9323ab35f4da5c246083ba9dcb1b751bca9c9b4

2025-10-27\_131923\_NB-320-Ch2-3.png

sha256: 7262e3ea29bdc3e4bb5904811ef4372b982f6f83f75b9548cf3e88fbaad01d6d



2025-10-27\_112151\_NB-320-Ch2-2.txt

sha256: 0e7bfddd8677db7f44c664ccc16f6a636377d9e51b4a3ad1eeec55ac5bcb16da

2025-10-27\_112151\_NB-320-Ch2-2.png

sha256: c673a5f4db1bb93e859a4ae84d0baa12f24c4a5bf1ad8a0c4dd3b25afee4e5a5




2025-10-27\_111158\_NB-320-Ch2-1.txt

sha256: 55f79c9ff8a90b5bcbd884e21aeae16fa9a3d227b9862aa39ba73a2240d15c0

2025-10-27\_111158\_NB-320-Ch2-1.png

sha256: cdbf3e331e1c8301ab3c050122877558d01b67348055dcbebd57b16c81c41de7



NB-320.ulog

sha256: 66c7c58d62eb69cf90ebda6f243332da31f4636a26ca8e7d4a04ede9c2e13b05



Unique eLabID: 20251027-15f09fafef1c2e5c3939717d4c0df694186b1897

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3262>