

NB-364: Gas phase H₂ and O₂ measurements of Al:SrTiO₃ RhCrO_x (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C, 70 min, degassing (reproduction NB-363)

Date: 2025-11-26

Tags: [O2](#) [Test](#) [Calibration](#) [NB](#) [Firing](#) [Irradiation](#) [O2 sensor](#) [H2](#) [advanced irradiation](#) [setup](#) [Unisense](#) [H2 Sensor](#) [temperature](#) [In situ](#) [Trace range](#) [robust oxygen sensor](#) [photocatalysis](#)

Category: SrTiO₃


Status: Done

Created by: Nadzeya Brezhneva

Objectives

Reproduction of NB-363: simultaneous detection of O₂ and H₂ evolution in gas phase for irradiated suspension of Rh_xCrO_y:Al:SrTiO₃ suspension (EA-358, 0.5 mg/mL), 365 nm LED, 50 mW/cm², 20 °C during 70 min.

Reaction scheme



ChemDraw file linked: [NB-362-SrTiO₃-photocatalytic H₂O splitting.cdxml](#)

Literature/reference experiments

Literature	/
Reproduction	SrTiO₃ - NB-363: Gas phase H₂ and O₂ measurements of Al:SrTiO₃ RhCrO_x (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C, 70 min, degassing (reproduction NB-362)
Similar experiments	SrTiO₃ - NB-362: Gas phase H₂ and O₂ measurements of Al:SrTiO₃ RhCrO_x (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C, 70 min, degassing SrTiO₃ - NB-361: Gas phase H₂ and O₂ measurements of Al:SrTiO₃ RhCrO_x (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C, 15 min, degassing

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	pressure [bar]
milli-Q H ₂ O	/	/	/	/	/	/	/	/	25	/

Al:SrTiO ₃ RhCrO _x (EA-358)	SrTiO ₃ - EA-358: Modification of Al:SrTiO ₃ (EA-354) via deposition of Rh, Cr oxide co-catalyst, 350°C, 1h, Upscaling (3.33x)	/	/	/	12.50	12.47	/	/	/	/
Hydrogen	1333-74-0	/	/	/	/	/	/	/	1 balloon (ca. 1 L)	ca. 1

Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-314: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup](#)

	Name
Used Set-up	Equipment - Advanced irradiation setup V1.0 I
Irradiation setup number	Equipment - Irradiation setup 4 (CEEC II, E002)

	Light Source Name	Power Source Name	Wavelength [nm]	Power Setting [mW]	Analog Setting [0.00 - 10.00]
First light source	Light Source - UHP LED 365 nm-4	Power Sources - BLS-18000-14	365	56	0.19

Used beam combiner [Name or None]	/
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	20
Stirring speed [rpm]	500

Irradiation start: 1. Firing [relative to start log] 2. Unisense	1. 615 s 2. 16:06:16
Irradiation stop: 1. Firing [relative to start log] 2. Unisense	1. 4830 s 2. 17:16:31

O₂/H₂ sensor equipment

	Equipment	Used protocol
Used Firing	Equipment - Firing Fiber-Optic Oxygen Meter 2 Channel (Firing 2)	Protocol - Operation of Firing Fiber-Optic Oxygen Meter 2 Channel Software
Used O ₂ sensor	Equipment - Robust probe for liquid O ₂ measurement	Protocol - In-situ hydrogen and oxygen measurement in H ₂ /O ₂ reactor
Used H ₂ sensor	Equipment - H ₂ UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurement in H ₂ /O ₂ reactor

Procedure/observations

Date	Time	Step	Observations	Pictures/Files
25.11.2025		The experiment was done according to Protocol - In-situ hydrogen and oxygen measurement in H₂/O₂ reactor Important steps and deviations are listed below	/	/
		Important note: replacing 10 mm BOLA fitting to 2 mm BOLA fitting for H ₂ sensor.		
	11:39-13:06	Conditioning of H ₂ sensor	NB-364-Logger1 offset - 2 mV	NB-364.ulong NB-364-Logger1-pre-polarization.csv NB-364-Logger1-pre-polarization.bmp
	13:15-20	Assembling the setup for calibration (25 mL of water was added using graduated cylinder), (LAUDA set to 20 °C) done according to Protocol - Liquid phase calibration of H₂ UniAmp sensor with H ₂ bubbling.	/	

	13:24	Start of O2 logging.	NB-364-Ch2-1	2025-11-26_132450_NB-364-Ch2-1.txt 2025-11-26_132450_NB-364-Ch2-1.png
	13:25	Start of H2 logging.	NB-364-Logger2 offset - 2 mV	NB-364.ulog NB-364-Logger2-calibration.csv NB-364-Logger2-calibration step.bmp
	13:27	Degassing was started.	/	20251126_132753-degassing of water.jpg
	14:11	0 ppm was taken.	/	/
	14:13	H2 bubbling of the reactor was started	/	20251126_141600-H2 bubbling.jpg
	14:18	1.000.000 ppm point was taken and calibration was used.	898 mV, slope 0.009, 97334 Pa	20251126_141816-H2 table.jpg
	14:22	Stop of H2 logging.	/	/
	14:22	Stop of O2 logging.	/	/
	14:30	Deassembling the setup, drying the reactor with acetone and compressed air .	/	/
		Sample preparation		
	14:50	Weighing EA-358 photocatalyst in a 50 mL vial.	Creamy solid.	/
	14:55	Addition of 25 mL H2O to the vial via graduated cylinder.	/	/
	14:59-15:02	The suspension was vortexed for 3 min (Equipment - VWR® VV3, Vortex Mixer , stage 4/6), covered with Al foil before further use.	/	20251126_150223-suspension after vortex.jpg
		Continue in Protocol - In-situ hydrogen and oxygen measurment in H2/O2 reactor from step 6		
	15:05	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	/
	15:10	Assembling the setup.	/	/

	15:18	Start of O2 logging.	NB-364-Ch2-2	2025-11-26_151821_NB-364-Ch2-2.txt 2025-11-26_151821_NB-364-Ch2-2.png
	15:21	The degassing was started	/	20251126_152140-degassing of the suspension.jpg
	15:54	The degassing was stopped by removing the cannula and closing the valve.	/	/
	15:55	Stop of O2 logging.	/	
	15:56	Start of O2 logging.	NB-364-Ch2-3	2025-11-26_155601_NB-364-Ch2-3.txt 2025-11-26_155601_NB-364-Ch2-3.png
	15:56	Start of H2 logging.	NB-364-Logger3	NB-364.ulog NB-364-Logger3-during irradiation.csv NB-364-Logger3-during irradiation.bmp
	15:56-16:06	Equilibration time.	/	/
	16:06	The irradiation was started	/	20251126_160700-after start of irradiation.jpg
	17:16	The irradiation was stopped.	/	/
	17:16-26	Equilibration time.	/	/
	17:26	Stop of O2 and H2 logging.	/	/
	17:30	Deassembling the setup, cleaning the reactor.	Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode).	20251126_175358-after irradiation.jpg

Analysis

Used calibration for Firing: 20250910-BOLA fitting-gas phase-4-neck photoreactor-trace oxygen robust probe-Ch2.ini

Used calibration for UniSense: NB-364-Logger2

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
25.11.2025	11:39	NB-364-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-364.ulong NB-364-Logger1-pre-polarization.csv	/	NB-364-Logger1-pre-polarization.bmp	/	Pre-polarization of the sensor.
	13:25	NB-364-Logger2	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-364.ulong NB-364-Logger2-calibration.csv	/	NB-364-Logger2-2point calibration.bmp NB-364-Logger2-calibration step.bmp	/	Calibration of H2 sensor, 10 ⁶ ppm corresponds to 898 mV
	15:56	NB-364-Logger3	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-364.ulong NB-364-Logger3-during irradiation.csv	NB-364-O2 and H2 curve.py	NB-364-Logger3-during irradiation.bmp NB-364-O2 and H2 curves.png	SrTiO3 - NB-363: Gas phase H2 and O2 measurements of Al:SrTiO3 RhCrOx (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C, 70 min, degassing (reproduction NB-362)	H2 evolution during irradiation.
	13:24	NB-364-Ch2-1	Optical O2 detection	Equipment - Firing Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-26_132450_NB-364-Ch2-1.txt	/	2025-11-26_132450_NB-364-Ch2-1.png	/	Degassing of water followed by calibration of H2 sensor.
	15:18	NB-364-Ch2-2	Optical O2 detection	Equipment - Firing Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-26_151821_NB-364-Ch2-2.txt	/	2025-11-26_151821_NB-364-Ch2-2.png	/	Degassing of the suspension.
	15:56	NB-364-Ch2-3	Optical O2 detection	Equipment - Firing Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-26_155601_NB-364-Ch2-3.txt	NB-364-O2 and H2 curve.py	2025-11-26_155601_NB-364-Ch2-3.png NB-364-O2 and H2 curves.png	SrTiO3 - NB-363: Gas phase H2 and O2 measurements of Al:SrTiO3 RhCrOx (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C, 70 min, degassing (reproduction NB-362)	O2 evolution during irradiation.

Results

Reproduction of NB-363: simultaneous H₂ and O₂ measurements (gas phase) of irradiated suspension of EA-358 (0.5 mg/mL) in O₂/H₂ photoreactor under 365 nm irradiation (50 mW/cm², 20 °C, 70 min) were performed. Problems with stirring were eliminated.

Linked experiments

SrTiO₃ - NB-312: Gas phase H₂ and O₂ measurements with Unisense H₂ sensor, Firesting O₂ robust probe in irradiated Al:SrTiO₃ RhCrOx (NB-289, 0.5 mg/mL), 365 nm, 50 mW, 1 h, degassing (reproduction NB-304)

SrTiO₃ - NB-361: Gas phase H₂ and O₂ measurements of Al:SrTiO₃ RhCrOx (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C, 15 min, degassing

SrTiO₃ - NB-362: Gas phase H₂ and O₂ measurements of Al:SrTiO₃ RhCrOx (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C, 70 min, degassing

SrTiO₃ - NB-363: Gas phase H₂ and O₂ measurements of Al:SrTiO₃ RhCrOx (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C, 70 min, degassing (reproduction NB-362)

Linked resources

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Robust probe for liquid O₂ measurment](#)

Equipment - [Advanced irradiation chamber V1.0 I](#)

Equipment - [H₂ UniAmp Sensor - Normal range - 2.1 x 80 mm needle](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)

Protocol - [Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software](#)

Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)


Protocol - [Gas phase calibration of H₂ UniAmp sensor](#)

Protocol - [In-situ hydrogen and oxygen measurment in H₂/O₂ reactor](#)

Attached files

NB-364-O₂ and H₂ curves.png

sha256: d1ffd0910a69df2e36e179909e31db09141503e5199e4fced4e2c32ea2b775a1



NB-364-O₂ and H₂ curve.py

sha256: b6c90f85944e2f6a1da38793e84801291c5405d9be7ce68de87e2e4d44443fc7

20251126_132753-degassing of water.jpg

sha256: 2926867415848323a468ab5637b1c46272be680333a837b2ddb10041938ca14d



20251126_141816-H2 table.jpg

sha256: 4f8a1df8f9ca432db06382709d8cbbd105ff1d0a2561ca9b2db0ce57ed59da3e



20251126_141600-H2 bubbling.jpg

sha256: 0ee3c850f4fb751fc2c4181e2cde2f3e0190e50837c4c294a0ffb1f2b1136285



20251126_152140-degassing of the suspension.jpg

sha256: f048a04f7b73345991faafd00b6fdde0c01d77ed868f4075ee36ba99c5f20fca



20251126_150223-suspension after vortex.jpg

sha256: a3cfa9c677ca3c3644f71ee2bc58b0586a1b73ec3ae368d452511df1faa7e903



20251126_175358-after irradiation.jpg

sha256: b0d8b23bc0c41d53fd918c30e3734cb0f75612f61b8b5b2a97c7d7c34d1ce59e



20251126_160700-after start of irradiation.jpg

sha256: 8c9678cc254ea99cfbcf39029c5fca1f60c58862d3f93ad591fa1f7d04039644




NB-364.ulog

sha256: 7c2d58dafd04d39060d445bf08bbb742c3dc703a5f11f7634f145104b3f55b27

2025-11-26_155601_NB-364-Ch2-3.png

sha256: a8726b69d2a33eb8d5242f5f00f074d4b9b20996f43f77d4b2d420b82facd830



2025-11-26_155601_NB-364-Ch2-3.txt


sha256: 8acf0a0f902d16b8307516d73c93e6a5528a674373e01ea25d02c4b02c68da7c

2025-11-26_151821_NB-364-Ch2-2.txt

sha256: 7759e93c459ded19449b086090c17477d68d765b500071319855b7d2a087d2ad

2025-11-26_151821_NB-364-Ch2-2.png

sha256: 6a72448bca7da404dd4a71fd9499c1c3ed12bfcc091500484a7e93f82cf992e8




2025-11-26_132450_NB-364-Ch2-1.txt

sha256: 71623dc9082c0e70e2c7b15f3ccfd6ec7743edac65f97203e97686fe6e26b58c

2025-11-26_132450_NB-364-Ch2-1.png

sha256: 72944a4d34f828737732c9d233d02c9c34f0ce109fd848d92f0ae1e4d7d3a703



NB-364-Logger2-calibration.csv

sha256: 11f0ceb3551697aa9ddd74e59e683e6285ab78927df1fadd1a1112a5528eaf3


NB-364-Logger2-calibration step.bmp

sha256: ad766c73a61dd5956574a4bea027fb1a74f40cdd39e089f548f20c2f729d57cc




NB-364-Logger2-2point calibration.bmp

sha256: 5fb5b8bf3b065450574f2d93d3782fed9e6bbef27840eeb047cb8b68862cd2af



NB-364-Logger1-pre-polarization.bmp

sha256: 1bf9e7e6a291f466dbff51ae1e9c8a872a8e699951a01187f5fad8b85db843d2



NB-364-Logger1-pre-polarization.csv


sha256: cd4b2b233600e733030006ada78fe4b957fc2a8251798d2e2998f8512cad1c51

NB-364-Logger3-during irradiation.csv

sha256: f4a3404ae328ea51e739f9009a2d88e66006b9ae36a358fa800cb6fb82d50845

NB-364-Logger3-during irradiation.bmp

sha256: bc5762ed3cad18c3c49c674f5ee6422fe1b02a3a02c7aa3563dccb3b3372aafc



Unique eLabID: 20251126-37d95891eb960ef7be82a4e9ff8be4dea9563e7b

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3563>