

NB-AE-222: Gas phase H2 measurements with Unisense H2 sensor, GC calibration

Date: 2025-05-16
Tags: Test AE Calibration Future NB GC H2 Gas chromatography Unisense H2 Sensor In situ
Category: SrTiO3
Status: Done
Created by: Nadzeya Brezhneva

Literature/reference experiments

Literature	/
Reproduction	/
Similar experiments	SrTiO3 - AE-483: Testing of in-situ electrochemical hydrogen sensor and test H2 evolution reaction with RhCrO:Al:SrTiO3 SrTiO3 - AE-484: Testing of in-situ electrochemical hydrogen sensor and test H2 evolution reaction with RhCrO:Al:SrTiO3 SrTiO3 - NB-AE-215 / AE-485: Gas phase H2 measurements with H2 sensor, SrTiO3:Al loaded with RhxCrYO3 (NB-162-4) suspension (1 mg/ml), 365 nm UHP LED, check with GC measurements

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	pressure [bar]
milli-Q H ₂ O	/	/	/	/	/	/	/	/	ca. 5 (for H2 saturated atmosphere in the RBF)	/
Hydrogen	1333-74-0	/	/	/	/	/	/	/	1 balloon (approx. 2 L)	approx. 1

Procedure/observations

Date	Time	Step	Observations	Pictures
21.05.2025	ca. 10:10	The Device was connected to the computer and the device was connected to the sensor and the PT100. The experiment was started on the device - the PreAmp was set at 1pA=1mV and the polarisation was 100 mV	/	/
	ca. 10:25	The logging was started (sample interval for pre-measurements was set 10 s).	NB-AE-222-Logger1	/
	10:25	Attempt to connect PT 100 sensor.	Problems with recognizing PT 100, by now no use of the sensor.	/
	ca. 10:30	A 10 mL RBF flask was filled with Milli-q water (approx. 5 mL) (flask 1)	/	20250521_104013-RBF-with-H2O.jpg
	ca. 10:45-50	Filling the balloon with H ₂ gas in E005 lab according to the protocol Protocol - Getting hydrogen from hydrogen bottle in CEEC II E014	/	/
	ca. 10:50	Start purging flask 1 with H ₂ from the balloon (two cannulas in the septum) according to the protocol Protocol - Gas phase calibration of H2 UniAmp sensor .	/	20250521_105118-H2-purging-with-balloon.jpg
	ca. 11:05	Important note: removing the second needle to reduce too fast release of H ₂ from the balloon.	/	/
	11:10-12:40	Continuous recording during pre-measurements were continued.	/	Unisense-pre-measurements-Screenshot-2025-05-21-160016.png

	12:40	PT 100 sensor was connected.	Now everything went well, the sensor is recognized.	/
	12:46	The logging was stopped	/	/
	ca. 12:50	The H2 sensor was placed on a GL14/NS14 5 mL Schlenk flask (flask 2) equipped with a PTFE stirring bar inside (Prep work - AE-482: Volume determination of GL14/NS14 5 mL Schlenk flask) using a 2.0 mm BOLA fitting	/	20250521_130448-GL14-NS14-flask-setup.jpg 20250521_130207-with-H2-sensor.jpg
	13:02	Start stirring (400 rpm), LAUDA on (25 °C), PT 100 sensor was placed inside the water bath.		
	13:06	A new logging was started.	NB-AE-222-Logger2	/
	13:12	The first calibration point was taken at 0 ppm H2	Calibration --> Offset --> Set value close to 0 (in current experiment -4 mV).	Unisense-calibration-Screenshot-2025-05-21-160445.png
		Calculations of the amount of saturated H₂ gas that is needed to be injected into a GL14/NS14 flask $V_{total}(GL14/NS14) = 11.66 \text{ mL}$ according to Prep work - AE-482: Volume determination of GL14/NS14 5 mL Schlenk flask $0.05 \% (V_{total}) = 0.0058 \text{ mL} = 5.8 \mu\text{L}$ (6 μL) $0.3 \% (V_{total}) = 0.0348 \text{ mL} = 34.8 \mu\text{L}$ (35 μL) $1 \% (V_{total}) = 0.116 \text{ mL} = \mathbf{116 \mu\text{L}}$		
	13:25	6 uL of H ₂ saturated gas was taken from the head space of flask 1 using 100 μL Pressure lok syringe and injected into flask 2 (via septum on the NS14 neck).	/	/
	13:25-27	Waiting for ca. 1-2 min for stabilizing the signal, afterwards adding point to the calibration curve (corresponded to the signal 2.8 mV).		
	13:29	Taking 100 uL of the sample from the head space of GL14/NS14 flask (preliminary flushing it with the gas phase from the head space) and kept in Pressure lok syringe for further investigations with GC.	NB-AE-222-gas 0.05 %theor H2	Unisense-signal-vs-time-Screenshot-2025-05-21-155932.png
	13:30-35	The sample was transferred to GC lab (314) in TO Building.		
	ca. 13:50	GC measurements of the sample NB-AE-222-gas 0.05 %theor H2 with Equipment - GC chromatograph Shimadzu Nexis GC-2030, Peneva's Lab 304, Physical Chemistry building in accordance with the protocol Protocol - Hydrogen Measurement Using GC .	Change of septum in the GC --> delay in start of the measurements in ca. 15 min	/
	14:07	For 0.3 vol.% H ₂ in total volume, 29 uL H ₂ saturated gas was taken from the head space of flask 1 using Pressure lok syringe and injected into flask 2 (via septum on the NS14 neck).	Since the signal after first injection of H2 sat gas was stable, it was necessary to inject just 29 uL (35 uL-6 uL) instead of the required 35 uL	/
		Important note: adjust the parameters in H ₂ table each time before injecting the sample (check current values of temperature 24.5°C and pressure 991 mbar). The input of the value in ppm will give the result in Pa at the bottom of the panel.	Screenshot-2025-05-27-184655-Unisense calculator.png	/
	14:08	Waiting for 1 min for stabilizing signal, addition of the third point to the calibration curve.	/	/
	ca. 14:10	Taking 100 uL of the sample from the head space of GL14/NS14 flask (preliminary flushing it with the gas phase from the head space) and kept in Pressure lok syringe for further investigations with GC.	NB-AE-222-gas 0.3 %theor H2	Unisense-signal-vs-time-Screenshot-2025-05-21-155932.png
	ca. 14:10-15	The sample was transferred to GC lab (314) in TO Building.	/	/
	ca. 14:20	GC measurements of the sample NB-AE-222-gas 0.3 %theor H2 with Equipment - GC chromatograph Shimadzu Nexis GC-2030, Peneva's Lab 304, Physical Chemistry building in accordance with the protocol Protocol - Hydrogen Measurement Using GC .	/	/
	ca. 14:30	For 1 vol.% H ₂ in total volume, 82 uL H ₂ saturated gas was taken from the head space of flask 1 using Pressure lok syringe and injected into flask 2 (via septum on the NS14 neck).	Since the signal after the second injection of H2 sat gas was stable, it was necessary to inject just 82 uL (117 uL-35 uL) instead of the required 117 uL	/
	ca. 14:32	Waiting for ca. 1-2 min for stabilizing signal, addition of the fourth point to the calibration curve.	/	/

	ca. 14:35	Taking 100 uL of the sample from the head space of GL14/NS14 flask (preliminary flushing it with the gas phase from the head space) and kept in Pressure lok syringe for further investigations with GC.	NB-AE-222-gas 1 %theor H2	Unisense-signal-vs-time-Screenshot-2025-05-21-155932.png
	ca. 14:40	The sample was transferred to GC lab (314) in TO Building.	/	/
	ca. 14:50	GC measurements of the sample NB-AE-222-gas 1 %theor H2 with Equipment - GC chromatograph Shimadzu Nexis GC-2030, Peneva's Lab 304, Physical Chemistry building in accordance with the protocol Protocol - Hydrogen Measurement Using GC.	/	/
	ca. 15:20	The calibration was applied	/	/
	15:22	The logging was stopped and subsequently the experiment was stopped	/	
	ca. 15:40	Deassembling the setup.	/	/
22.05.2025	ca. 15:40	GC calibration using 25 uL, 50 uL, 75 uL, 100 uL H ₂ with 1 % All-in-gas reference gas mixture (1 vol.%) and Pressure lok syringe using Equipment - GC chromatograph Shimadzu Nexis GC-2030, Peneva's Lab 304, Physical Chemistry building in accordance with the protocol Protocol - Hydrogen Measurement Using GC.	/	/

Analysis

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Processed Data	Comparative Data	Interpretation
21.05.2025	10:25-12:46	NB-AE-222-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor	/	NB-222-1.ulong NB-AE-222-log1-pre-measurements.xlsx Calibration-data-for-pre-measurements.csv Comments-pre-measurements.csv Data-1-H2-UniAmp-503157-pre-measurements.csv Devices-pre-measurements.csv Miscellaneous-pre-measurements.csv	Unisense-pre-measurements-Screenshot-2025-05-21-160016.png	/	Pre-measurements are done, the signal is stable during measurements.
	13:06-15:35	NB-AE-222-Logger2	electrochemical H2 detection	Equipment - H2 UniAmp Sensor	/	NB-222-1.ulong NB-AE-222-log2-calibration.xlsx Calibration-data-for-calibration.csv Comments-for-calibration.csv Data-1-H2-UniAmp-503157-for-calibration.csv Devices-for-calibration.csv Miscellaneous-for-calibration.csv	Unisense-calibration-Screenshot-2025-05-21-160445.png Unisense-signal-vs-time-Screenshot-2025-05-21-155932.png Unisense-calculations.xlsx	/	<i>Signal vs. time</i> Step-by-step increase in signal value after gradual injection of H ₂ -containing mixture. <i>4-point calibration</i> R ² = 1, intercept = -1.031, slope = 0.076. After processing Unisense data: theor 0.05 vol.% ----> 0.0517 vol.% theor 0.3 vol.% ----> 0.3001 vol.% theor 1 vol.% ----> 1.0004 vol.%
	13:29	NB-AE-222-gas 0.05 % theor H2	GC	Equipment - GC chromatograph Shimadzu Nexis GC-2030, Peneva's Lab 304, Physical Chemistry building	/	Calibration-Unisense.gcb Blank1.gcd theor-0.05-vol.gcd theor-0.05-vol-H2-GC.txt	theor-0.05-vol.-H2.PNG	/	Rough estimation by LabSolution software (conc vol.%) 0.239, much higher than expected
	14:10	NB-AE-222-gas 0.3 % theor H2	GC	Equipment - GC chromatograph Shimadzu Nexis GC-2030, Peneva's Lab 304, Physical Chemistry building	/	Calibration-Unisense-1.gcb Blank1.1.gcd theor-0.3-vol.gcd theor-0.3-vol-H2-GC.txt	theor-0.3-vol.-H2.PNG	/	Rough estimation by LabSolution software (conc vol.%) 0.569, higher than expected
	14:35	NB-AE-222-gas 1 % theor H2	GC	Equipment - GC chromatograph Shimadzu Nexis GC-2030, Peneva's Lab 304, Physical Chemistry building	/	Calibration-Unisense-2.gcb Blank1.2.gcd theor-1-vol.gcd theor-1-vol-H2-GC.txt	theor-1-vol.-H2.PNG	/	Rough estimation by LabSolution software (conc vol.%) 1.122, a bit higher, but comparable with expected value
22.05.2025	15:40	/	GC	Equipment - GC chromatograph Shimadzu Nexis GC-2030, Peneva's Lab 304, Physical Chemistry building	/	NB-222-calibration.gcb Blank002.gcd 100H2-calibration-GC.gcd 100H2-calibration-GC.txt 75H2-calibration-GC.gcd 75H2-calibration-GC.txt 50H2-calibration-GC.gcd 50H2-calibration-GC.txt 25H2-calibration-GC.gcd 25H2-calibration-GC.txt	100-H2.PNG 75-H2.PNG 50-H2.PNG 25-H2.PNG GC_Data_Analysis_H2-1.py NB_222_Data.csv NB_222_Data-Results.txt NB_222_Data-Analysis.tif	/	After calibration and processing data in Python: according to GC: theor 0.05 vol.% ----> 0.252 vol.% theor 0.3 vol.% ----> 0.588 vol.% theor 1 vol.% ----> 1.151 vol.%

Results

H₂ test measurements with Unisense hydrogen sensor was performed. GC measurements (as well as 4-point calibration with 1 vol.% H₂ gas mixture) was performed as a control.

Four points were tested - 0, 0.05 vol.%, 0.3 vol.%, 1 vol.% (theoretical values)

According to the Unisense data:

theor 0.05 vol.% ---> 0.0517 vol.%

theor 0.3 vol.% ---> 0.3001 vol.%

theor 1 vol.% ---> 1.0004 vol.%

Used as calibration data for Unisense.

Check with GC (after 4 point calibration with 1 vol.% H₂ gas mixture and processing data in Python):

theor 0.05 vol.% ---> 0.252 vol.%

theor 0.3 vol.% ---> 0.588 vol.%

theor 1 vol.% ---> 1.151 vol.%

Unisense data are in a very good agreement with theoretical values, whereas GC data (especially theor 0.05 vol.% and 0.3 vol.%) are not in good agreement with Unisense and theoretical values, high error.

Future recommendations

Old procedure	Problem	Future recommendations
/	/	placing water into Schlenk where H ₂ measurement is performed.
/	/	In future: exporting only data-1 csv file from Unisense Logger software.

Linked experiments

Prep work - [AE-477: Testing of in-situ electrochemical hydrogen sensor and test H₂ evolution reaction with RhCrO:Al:SrTiO₃](#)

Prep work - [AE-482: Volume determination of GL14/NS14 flask](#)

SrTiO₃ - [NB-162: Preparation of SrTiO₃:Al with RhCr₂-yO₃ \(Osterloh route\), I attempt, 1000 C 10h](#)

SrTiO₃ - [AE-483: Testing of in-situ electrochemical hydrogen sensor and test H₂ evolution reaction with RhCrO:Al:SrTiO₃](#)

SrTiO₃ - [AE-484: Testing of in-situ electrochemical hydrogen sensor and test H₂ evolution reaction with RhCrO:Al:SrTiO₃](#)

SrTiO₃ - [NB-AE-215 / AE-485: Gas phase H₂ measurements with H₂ sensor, SrTiO₃:Al loaded with RhxCr_yO₃ \(NB-162-4\) suspension \(1 mg/ml\), 365 nm UHP LED, check with GC measurements](#)

Linked resources

Equipment - [GC chromatograph Shimadzu Nexis GC-2030, Peneva`s Lab 304, Physical Chemistry building](#)

Equipment - [H₂ UniAmp Sensor - Low range - 2.1 x 80 mm needle](#)

General Protocol - [Hydrogen Measurement Using GC \(AG Paneva\)](#)

General Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)

General Protocol - [Gas phase calibration of H₂ UniAmp sensor](#)

Light Source - [UHP LED 365 nm-1](#)

Attached files


Screenshot-2025-05-27-184655-Unisense-calculator.png

sha256: 11af0eaf15a97d7622a53dbe605b30a1eec4fee21b6ab8889b82eb9b634b9ca9



NB_222_Data-Analysis.tiff

sha256: 08def4db49886d8f6e31bd6614187afc20d366414c9d680acee8884e2c1f4db9



NB_222_Data-Results.txt

sha256: a1905d4573b759bdfed9bdaba9376dc76d3cbad8752b16077144d12d0929029b

NB_222_Data.csv

sha256: 228edee5a09c46fb915fe72ed1dc8397c47036e0de52d345198cdde71115dc1d

GC_Data_Analysis_H2-1.py

sha256: a99a287753dc68423fb003ce77237e3aaa2605af26edbae977965923b6d58979

NB_222_Calibration.csv

sha256: 6adf3c54b0ee93ed2865b7bfeaaa6b08801427ef1c3233d93d0ab66964af0396

Unisense-calculations.xlsx


sha256: 13671b8e3812c6472d55584fddf8e2dca1bd5eb99b65906618b3e0821cf58a86

75H2-calibration-GC.gcd

sha256: 82eacb3bc107ae7a82cf64beefc12abc159f79aa9f6602621cf866823ee87f7a

75-H2.PNG

sha256: 93f877d57226d40e1903cb068513592f14ee45663b8cc7a6dfe0f959ede5014d




75H2-calibration-GC.txt

sha256: cf5d2df7500253a78f43b64d56e35c49c40f0929f34e497d4e2f53bc1e95d612

100-H2.PNG

sha256: 9f9ac972a6de30b171c3a4cf7b6fd6c6028b5dce6da36ebc2d6dcf6278ca3745



100H2-calibration-GC.txt

sha256: cfab6c7655907b6f69bb4e6002b223788b86bbb29919c07d6eac8f0914cf075b

Blank002.gcd

sha256: af960c68ca7977641fbf9ce86ddae1bc3b23fc3e31f115188cb5617e671a0346

NB-222-calibration.gcb

sha256: 89acaafd18d403ec47ecbc78efeee8cd19bb5f7b9f11a007aefd9f499bb50f1c

100H2-calibration-GC.gcd


sha256: 9bd5f8a2e2b74626b5cd0bbc13bd9ba6cd9c521f4c69004cdbc8df6d4342be87

25H2-calibration-GC.txt

sha256: cf7679bea0686cad75ed48aa7485e0d6be1f4c91b66d7e3b8540471c8aeb2de

25-H2.PNG

sha256: fae0863460215d3c9d4491ad6af7d8d3553c2322f245f1f09de940d9b10c7390




25H2-calibration-GC.gcd

sha256: 606cfc83bdb614e12fec2921f720323140f0c8aaded0d8a7ab72bb3573384d53

50-H2.PNG

sha256: e5fd5eaa6ca4e9459f62c2d250c356dc191ec121dcf71d6aacf6bc106777c98d



50H2-calibration-GC.gcd

sha256: d7a41c53eaac11de255619ca0c6a25c86509c82b797e2c36ff0660aa7a89c7d7

50H2-calibration-GC.txt

sha256: fc17c43a74e571002207f90038862643d8cf93fe8b96aff404d1af5b33d27083

theor-1-vol-H2-GC.txt

sha256: 5bf7a92b6d98dff15f134e3e0c14d7cc2494368a91fcee06a447b55254318c0

theor-0.05-vol-H2-GC.txt

sha256: d4d75c681474284609c38285d705eac8d5b97a6fe7ee811d7da665bb805191d1

theor-0.3-vol-H2-GC.txt

sha256: bc489dfe2bb3404bc5e71c7b50699c6e4eb3e66b4c0d327d0b6e2cec944ad9f8


20250521_104013-RBF-with-H2O.jpg

sha256: df3505b47836756e6f03ebbc5a88983321944f8791c2d6864e9615b6c2c6133f



20250521_105118-H2-purging-with-balloon.jpg

sha256: 0060a4ba6a514b31192d92ce6c682b7c9c17c65bf97ff5e515be3ee2d28c76d5



20250521_130207-with-H2-sensor.jpg

sha256: 35a20815ead61130ec6ec0c09e34918977481413fc38fbc3c0c111350556ccda



20250521_130448-GL14-NS14-flask-setup.jpg

sha256: 73f969fbd0364e6957737d541ee17807f6632d01324912aeb78c52d95af4fe



Calibration-Unisense-2.gcb

sha256: c9d05f84b1801595ae0c5b8237188319f9f163c5b3d510c14559935f776db8e9

Calibration-Unisense-1.gcb

sha256: e95480a8e6d097b65133930ca3bbdcb77712695d2289e1047c65fe27ce66e903

Calibration-Unisense.gcb

sha256: 170533ffc58615fb627c73980861e5c218887020cea12a5a35c254720974b7dc

Blank1.gcd

sha256: af18dcc9ef8843c5965ff3e40ae20bc266d381cdc37b344bb30a3b18c7ed4a

Blank1.2.gcd


sha256: e68d1c5ade755b77ecef68b0a376fc3dc9c2901fca85387dbcdf13d75d6e1564

Blank1.1.gcd

sha256: ce80e2224ee03b44cb91f86498831342afd864c177ec17890731b11a54b631a4

theor-1-vol.-H2.PNG

sha256: 0ec23f842d42374d6802a2cc3c431b355c587888b2f4431c4264261d53cb7877



theor-1-vol.gcd


sha256: df7cec9952334107aa2a611c75e556545eeaae40fc388d647cf334651824df01

theor-0.05-vol.gcd

sha256: 4cffd205821141a2f4758a12fd904f49f2f6d16522302689430e171f52fc5982


theor-0.05-vol.-H2.PNG

sha256: 038f224b7723e792d394c3f7a92ac4dad57cb8d9415921958efd9aca13e07c30




theor-0.3-vol.gcd

sha256: e2cd3ae76a374c2c5c011ce55324d0e9a0c71d0d7a982e42e4a78f9e79b312e4




theor-0.3-vol.-H2.PNG

sha256: d52677e9e349fab15a7c506cd8f48f115dc16ef4cdd0583a6eed8170f8f5825f



Unisense-calibration-Screenshot-2025-05-21-160445.png

sha256: c6f53137d5709764a159b9a5dd2ff73829de837e07ade4407012cb3716b39e73



Unisense-signal-vs-time-Screenshot-2025-05-21-155932.png

sha256: 50774b4b8911a16cf3ef0ed4e9d66d39069de9ccc0df5e8734305a7384d72029



Unisense-pre-measurements-Screenshot-2025-05-21-160016.png

sha256: 6a03329e452d2011f4b23b5e4a13bc42cb581c491322ec7894fc9c299f3173b8



Miscellaneous-pre-measurements.csv

sha256: 20e130a088a8323888633615ec6438955b3361c0eb6558e65b60ea18f3281297

Miscellaneous-for-calibration.csv

sha256: 20e130a088a8323888633615ec6438955b3361c0eb6558e65b60ea18f3281297

Devices-pre-measurements.csv

sha256: f34bd2c51fc0dd1b10e86d60221b1a50f2580f2d3d47e4398a59cdfa0cfea7e5

Devices-for-calibration.csv

sha256: f34bd2c51fc0dd1b10e86d60221b1a50f2580f2d3d47e4398a59cdfa0cfea7e5

Data-1-H2-UniAmp-503157-pre-measurements.csv

sha256: 931ccffa3e1139d10a35937eb6f07e7d94a94bc6eece2264a8bd3754004e6845

Data-1-H2-UniAmp-503157-for-calibration.csv

sha256: 25b5ef10fde173e62c11ecb54d8b7347d581f1f19053fddde5f9b7f7b01c6428

Comments-pre-measurements.csv

sha256: 4aa81964d195af15a2cb230da0cdf929bf3fa1e8b671e563972a03a18fb98067

Comments-for-calibration.csv

sha256: 4aa81964d195af15a2cb230da0cdf929bf3fa1e8b671e563972a03a18fb98067

Calibration-data-for-calibration.csv

sha256: 568030a0e4c29c04e6cb8120d3f500f2c3ca3e7e492a8d6ae1e3363511a608cb

Calibration-data-for-pre-measurements.csv

sha256: 568030a0e4c29c04e6cb8120d3f500f2c3ca3e7e492a8d6ae1e3363511a608cb

NB-AE-222-log2-calibration.xlsx

sha256: 2181dea31ff813012ec1dbc581294469f17c9ab664bbf76cf4474613bb13e03b

NB-AE-222-log1-pre-measurements.xlsx

sha256: 11649e79b8f56d1e6d12c188eb548fae762688874a7035e9a4d7ef6a2c8ec05e

NB-222-1.ulog

sha256: 139fb27706a12def8b7eb712d787bc97040858e8d1a715257bfe504c3a136cbc

Comments

On 2025-05-27 12:14:38 Jacob Schneidewind wrote:

Notes:

- * Adding "in situ" tag
- * Removing literature
- * Linking H2 sensor calibration protocol (gas phase)
- * Linking protocol for how to fill H2 balloon
- * Specifying where PT100 was placed, ideally also including a picture with Pt100 in the future
- * 13:25: specifying that 100 ul pressure lok syringe is used
- * Linking protocol for measurements of GC samples
- * 14:07: including screenshot of H2 table in software, explaining this in a bit more detail
- * In the future only exporting data-1 csv file
- * Analysis table: splitting data for log 1 and log 2

* Future recommendation: placing water into Schlenk where H2 measurement is performed (also adjusting this in general protocol)

On 2025-05-27 19:05:59 Nadzeya Brezhneva wrote:

Notes:

- * Adding "in situ" tag -> done
- * Removing literature -> done
- * Linking H2 sensor calibration protocol (gas phase) -> done
- * Linking protocol for how to fill H2 balloon -> done
- * Specifying where PT100 was placed, ideally also including a picture with Pt100 in the future -> done, ok, the note will be taken into account
- * 13:25: specifying that 100 ul pressure lok syringe is used -> done
- * Linking protocol for measurements of GC samples -> done
- * 14:07: including screenshot of H2 table in software, explaining this in a bit more detail -> done
- * In the future only exporting data-1 csv file -> done, it will be taken into account
- * Analysis table: splitting data for log 1 and log 2 -> done

* Future recommendation: placing water into Schlenk where H2 measurement is performed (also adjusting this in general protocol) --> done



Unique eLabID: 20250516-d194458a77d779c83a79b1695c2138f354aa148e
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=2173>