

# NB-355: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrO<sub>x</sub>,Al:SrTiO<sub>3</sub> (NB-342, 0.2 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 20 °C

Date: 2025-11-19

Tags: [O2](#) [Test](#) [Calibration](#) [NB](#) [Firing](#)  
[O2 sensor](#) [H2](#) [SrTiO3](#) [Unisense](#)  
[RhCrO3:Al:SrTiO3](#) [H2 Sensor](#)  
[temperature](#) [In situ](#) [Trace range robust](#)  
[oxygen sensor](#) [photocatalysis](#) [Unisense](#)  
[normal range](#)

Category: SrTiO<sub>3</sub>


Status: Done

Created by: Nadzeya Brezhneva

## Objectives

Simultaneous detection of H<sub>2</sub> and O<sub>2</sub> evolution in liquid phase for irradiated suspension of Rh<sub>x</sub>Cr<sub>y</sub>O<sub>3</sub>:Al:SrTiO<sub>3</sub> suspension (NB-342, 0.2 wt% Rh,Cr, 0.5 mg/mL), 365 nm LED, 50 mW/cm<sup>2</sup>, 20 °C.

## Reaction scheme



ChemDraw file linked: [NB-328-SrTiO3-photocatalytic H2O splitting.cdxml](#)

## Literature/reference experiments

Literature	/
Reproduction	/
Similar experiments	<a href="#">SrTiO3 - NB-316: Liquid phase H2 and O2 of RhCrO<sub>x</sub>,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 20 °C</a> <a href="#">SrTiO3 - NB-328: Liquid phase H2 and O2 of RhCrO<sub>x</sub>,Al:SrTiO3 (NB-323, 0.2 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 20 °C</a>

## Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass <sub>theo</sub> [mg]	Mass <sub>exp</sub> [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	Concentration [mM]
milli-Q H <sub>2</sub> O	/	/	/	/	/	/	/	0.998	25	/
Al:SrTiO <sub>3</sub> RhCrO <sub>x</sub> (NB-342)	<a href="#">SrTiO3 - NB-342: Modification of EA-354 (SrTiO3:Al, upscaled batch) with Rh, Cr oxide cocatalyst (0.2 wt%), fresh solutions of RhCl3 and Cr(NO3)3</a>	/	/	/	12.50	12.54	/	/	/	/

# Excel sheet for reagent calculation

/

## Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-314: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup](#)

	Name
Used Set-up	<a href="#">Equipment - Advanced irradiation setup V1.0 I</a>
Irradiation setup number	<a href="#">Equipment - Irradiation setup 4 (CEEC II, E002)</a>

	Light Source Name	Power Source Name	Wavelength [nm]	Power Setting [mW]	Analog Setting [0.00 - 10.00]
First light source	<a href="#">Light Source - UHP LED 365 nm-4</a>	<a href="#">Power Sources - BLS-18000-1 4</a>	365	56	0.19

Used beam combiner [Name or None]	/
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	20
Stirring speed [rpm]	500
Irradiation start: 1. Firesting [relative to start log] 2. Unisense	1. 603 s 2. 22:18:33
Irradiation stop: 1. Firesting [relative to start log] 2. Unisense	1. 1510 s 2. 22:43:38

## O<sub>2</sub>/H<sub>2</sub> sensor equipment

	Equipment	Used protocol
Used Firesting	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)	Protocol - Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software
Used O <sub>2</sub> sensor	Equipment - Robust probe for liquid O <sub>2</sub> measurment	Protocol - In-situ hydrogen and oxygen measurment in H <sub>2</sub> /O <sub>2</sub> reactor
Used H <sub>2</sub> sensor	Equipment - H <sub>2</sub> UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurment in H <sub>2</sub> /O <sub>2</sub> reactor

## Procedure/observations

Date	Time	Step	Observations	Pictures/Files
		Calibration from SrTiO <sub>3</sub> - NB-353: Liquid phase H <sub>2</sub> and O <sub>2</sub> of RhCrOx,Al:SrTiO <sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm <sup>2</sup> , 20 °C VIII (reproduction NB-316) was used.	/	/
19.11.2025		<b>Sample preparation</b>		
	21:00	Weighing photocatalyst in a 50 mL vial (the vial was covered with Al foil before further use) (preliminary, the lumps of the solid were broken with Smartspatula inside the vial).	Slightly grey solid	/
	21:05	Addition of 25 mL H <sub>2</sub> O to the vial via graduated cylinder.	/	/
	21:14-17	The suspension was vortexed for 3 min ( Equipment - VWR® VV3, Vortex Mixer, stage 4/6), covered with Al foil before further use.	/	20251119_212245-suspension after vortex.jpg
		Continue in Protocol - In-situ hydrogen and oxygen measurment in H <sub>2</sub> /O <sub>2</sub> reactor from step 6		
	21:25	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	/
	21:30	Assembling the setup.	Currently, stopper instead of H <sub>2</sub> sensor, PT100, PT1000 and O <sub>2</sub> robust probe are inside the reactor immersed in the liquid phase	/
	21:33	Start of O <sub>2</sub> logging.	<b>NB-355-Ch2-1</b>	2025-11-19_213326_NB-355-Ch2-1.txt 2025-11-19_213326_NB-355-Ch2-1.png

	21:36	The degassing was started	/	<a href="#">20251119_213729-degassing of the suspension.jpg</a>
	22:09	Cannula was transferred to gas phase, above the suspension.	/	/
	22:12	H <sub>2</sub> sensor was added in Ar counterflow.	/	/
	22:15	The degassing was stopped by removing the cannula and closing the valve.	/	/
	22:17	<a href="#">Stop of O2 logging.</a>	/	/
	22:18	<a href="#">Start of O2 logging.</a>	<b>NB-333-Ch2-2</b>	<a href="#">2025-11-19_221828_NB-355-Ch2-2.txt</a> <a href="#">2025-11-19_221828_NB-355-Ch2-2.png</a>
	22:18	<a href="#">Start of H2 logging.</a>	<b>NB-333-Logger1</b> ca. 22:26 - very slight movement of the holder to fit in the setup box completely	<a href="#">NB-355.ulong</a> <a href="#">NB-355-Logger1.csv</a> <a href="#">NB-355-Logger1.bmp</a>
	22:18-28	Equilibration time.	/	/
	22:28	The irradiation was started	/	<a href="#">20251119_222851-after start of irradiation.jpg</a>
	22:43	The irradiation was stopped.	/	/
	22:43-53	Equilibration time.	/	/
	22:53	<a href="#">Stop of O2 and H2 logging.</a>	/	/
	ca. 23:00	Deassembling the setup, cleaning the reactor.	Tips of the sensors and reactor were covered with attached photocatalyst particles. Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode).	<a href="#">20251119_231806-after irradiation.jpg</a>

## Analysis

Used calibration for Firesting: [20250611-BOLA-fitting-liquid-phase-trace-oxygen-sensor-H2-O2 reactor.ini](#)

Used calibration for UniSense: NB-353-Logger2

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
19.11.2025	22:18	NB-355-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-355.ulong NB-355-Logger1.csv	NB-355-O2 and H2 curve.py	NB-355-Logger1.bmp NB-355-O2 and H2 curves.png	/	H2 evolution during irradiation.
	21:33	NB-355-Ch2-1	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-19_213326_NB-355-Ch2-1.txt	/	2025-11-19_213326_NB-355-Ch2-1.png	/	Degassing of the suspension.
	22:18	NB-355-Ch2-2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-19_221828_NB-355-Ch2-2.txt	NB-355-O2 and H2 curve.py	2025-11-19_221828_NB-355-Ch2-2.png NB-355-O2 and H2 curves.png	/	O2 evolution during irradiation.

## Results

Simultaneous H<sub>2</sub> and O<sub>2</sub> measurements of irradiated suspension of NB-342 (0.5 mg/mL, 0.2 wt% RhCr) in O<sub>2</sub>/H<sub>2</sub> photoreactor under 365 nm irradiation (50 mW/cm<sup>2</sup>, 20 °C, 15 min) were performed.

## Linked experiments

SrTiO<sub>3</sub> - NB-316: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrOx,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 20 °C

SrTiO<sub>3</sub> - NB-318: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrOx,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm<sup>2</sup>, 20 °C

SrTiO<sub>3</sub> - NB-322: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrOx,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm<sup>2</sup>, 20 °C (reproduction NB-318)

SrTiO<sub>3</sub> - NB-323: Modification of EA-354 (SrTiO<sub>3</sub>:Al, upscaled batch) with Rh, Cr oxide cocatalyst (0.2 wt%)

SrTiO<sub>3</sub> - NB-325: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrOx,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm<sup>2</sup>, 20 °C

SrTiO<sub>3</sub> - NB-326: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrOx,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm<sup>2</sup>, 20 °C (reproduction NB-325)

SrTiO<sub>3</sub> - NB-327: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrOx,Al:SrTiO<sub>3</sub> (NB-321, 0.05 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 20 °C

SrTiO<sub>3</sub> - NB-328: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrOx,Al:SrTiO<sub>3</sub> (NB-323, 0.2 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 20 °C

SrTiO<sub>3</sub> - NB-331: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrOx,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 20 °C (reproduction NB-316) I I I

SrTiO<sub>3</sub> - NB-333: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrOx,Al:SrTiO<sub>3</sub> (NB-323, 0.2 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 20 °C (reproduction NB-328)

SrTiO<sub>3</sub> - NB-342: Modification of EA-354 (SrTiO<sub>3</sub>:Al, upscaled batch) with Rh, Cr oxide cocatalyst (0.2 wt%), fresh solutions of RhCl<sub>3</sub> and Cr(NO<sub>3</sub>)<sub>3</sub>

SrTiO<sub>3</sub> - NB-353: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrOx,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm<sup>2</sup>, 20 °C VIII (reproduction NB-316)

## Linked resources

Equipment - [VWR® VV3, Vortex Mixer](#)

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)

Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)

Protocol - [In-situ hydrogen and oxygen measurment in H2/O2 reactor](#)


## Attached files

NB-355-O2 and H2 curve.py

sha256: 989c9819cc27cb310a4a32e6832dc616a9430c9a8e59f37e6bc1ea51768b1e1c

NB-355-O2 and H2 curves.png

sha256: 682cafea68a4d453a8d547f54d4c83936c1dc6e8b513c8edeef09607afc3d679



20251119\_222851-after start of irradiation.jpg

sha256: 86a7c323b95b5cf33d9d23f9c7f36947f16144a155a085295974c22c3b270405



20251119\_231806-after irradiation.jpg

sha256: 3e4fbfc58bd27f984ef97677af0023bf5116cb8a26f72c8602da9e828d2f2a02



20251119\_213729-degassing of the suspension.jpg

sha256: c024d8a3ec5e773357b0b70e35084410f00dfa9a6ea4980aee2fc26be0a5c1ea




20251119\_212245-suspension after vortex.jpg

sha256: 6bdf08b99dab4b3b8d4a7835b13b9e49ed5f89dea7fcbfd28ac346e2192d6b39



2025-11-19\_221828\_NB-355-Ch2-2.png

sha256: a2853bd9d4cf96bb4780ac7a1f2cbb3757d4803169d1c2bb2114ad4fffc95fde



2025-11-19\_221828\_NB-355-Ch2-2.txt

sha256: 57e29f29f558594eb7275dbea404274842d24e1ea831e5df5bca72eb1ced629a

2025-11-19\_213326\_NB-355-Ch2-1.png

sha256: 91d6a8506e8c2bc64517ca63f3a8febb8dbcccbf0bfbcbbe811ea90ca6f9122e




2025-11-19\_213326\_NB-355-Ch2-1.txt

sha256: 64a2dd79929ed7eef2be9518d87bfd71da54980c4ef4ff8277da1b6df88d6964

NB-355-Logger1.bmp

sha256: 56bdb37443f8c4a6bdef45e856692b3f70ab5d2e2a8a40b5a30565efc106095d



NB-355-Logger1.csv

sha256: 6dfe450cc884040c2bd6ee7e70d5eff2b9de1c4bb225304e79f250086ae8a200

NB-355.ulog

sha256: eca66793b5b7b0b42edbc14fe744dc27892dd9b9891d764b5aaa64136edb5bfa



Unique eLabID: 20251119-f7dd5a411637b0160f0c812827ae6390ad1ba721  
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3502>