

NB-344: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 150 mW/cm², 20 °C

Date: 2025-11-12

Tags: O₂ Test Calibration NB Firing
O₂ sensor H₂ SrTiO₃ Unisense
RhCrO₃:Al:SrTiO₃ H₂ Sensor
temperature In situ Trace range robust
oxygen sensor photocatalysis Unisense
normal range

Category: SrTiO₃


Status: Done

Created by: Nadzeya Brezhneva

Objectives

Simultaneous detection of H₂ and O₂ evolution in liquid phase for irradiated suspension of Rh_xCr_yO₃:Al:SrTiO₃ suspension (EA-358, 0.5 mg/mL), 365 nm LED, **150 mW/cm²**, 20 °C (changing intensity relatively to reference experiment).

Reaction scheme



ChemDraw file linked: [NB-344-SrTiO3-photocatalytic H2O splitting.cdxml](#)

Literature/reference experiments

Literature	/
Reproduction	/
Similar experiments	SrTiO₃ - NB-318: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm², 20 °C SrTiO₃ - NB-316: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	Concentration [mM]
milli-Q H ₂ O	/	/	/	/	/	/	/	0.998	25.00	/
Al:SrTiO ₃ RhCrO _x (EA-358)	SrTiO₃ - EA-358: Modification of Al:SrTiO₃ (EA-354) via deposition of Rh, Cr oxide co-catalyst, 350°C, 1h, Upscaling (3.33x)	/	/	/	12.50	12.44	/	/	/	/

Excel sheet for reagent calculation

/

Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-343: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup I I I](#)

	Name
Used Set-up	Equipment - Advanced irradiation setup V1.0 I
Irradiation setup number	Equipment - Irradiation setup 4 (CEEC II, E002)

	Light Source Name	Power Source Name	Wavelength [nm]	Power Setting [mW]	Analog Setting [0.00 - 10.00]
First light source	Light Source - UHP LED 365 nm-4	Power Sources - BLS-18000-1 4	365	170	0.61

Used beam combiner [Name or None]	/
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	20
Stirring speed [rpm]	500
Irradiation start: 1. Firesting [relative to start log] 2. Unisense	1. 601 s 2. 22:08:43
Irradiation stop: 1. Firesting [relative to start log] 2. Unisense	1. 1530 s 2. 22:24:13

O₂/H₂ sensor equipment

	Equipment	Used protocol
Used Firesting	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)	Protocol - Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software
Used O ₂ sensor	Equipment - Robust probe for liquid O ₂ measurment	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor
Used H ₂ sensor	Equipment - H ₂ UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor

Procedure/observations

Date	Time	Step	Observations	Pictures/Files
11.11.2025		Calibration from experiment SrTiO ₃ - NB-339: Liquid phase H ₂ and O ₂ of RhCrO _x ,Al:SrTiO ₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm ² , 20 °C VI (reproduction NB-316) was used.		
		Sample preparation		
	20:35	Weighing EA-358 photocatalyst in a 50 mL vial.	Creamy solid	/
	20:40	Addition of 25 mL H ₂ O to the vial via graduated cylinder, covering with Al foil before further use.	/	/
	20:42-45	The suspension was vortexed for 3 min (Equipment - VWR® VV3, Vortex Mixer, stage 4/6), covered with Al foil before further use.	/	20251111_204556-suspension after vortex.jpg
		Continue in Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor from step 6		
	20:50	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	/
	20:55	Assembling the setup.	Currently, stopper instead of H ₂ sensor, PT100, PT1000 and O ₂ robust probe are inside the reactor immersed in the liquid phase	
	21:00	Start of O ₂ logging.	NB-344-Ch2-1	2025-11-11_210025_NB-344-Ch2-1.txt 2025-11-11_210025_NB-344-Ch2-1.png

	21:03	The degassing was started	/	20251111_210354-degassing of the suspension.jpg
	21:45	Cannula was transferred to gas phase, above the suspension.	/	/
	21:47	H ₂ sensor was added in Ar counterflow.	/	/
	21:51	The degassing was stopped by removing the cannula and closing the valve.	/	/
	21:53	Stop of O2 logging.	/	/
	21:53	Start of O2 logging.	NB-344-Ch2-2	2025-11-11_215318_NB-344-Ch2-2.txt 2025-11-11_215318_NB-344-Ch2-2.png
	21:53	Start of H2 logging.	NB-344-Logger1	NB-344.ulong
	21:58	Stop of H2 and O2 loggings.	Spike in H2 sensor signal --> stop logging, adjusting the cabel position	
	21:58	Start of O2 logging.	NB-344-Ch2-3	2025-11-11_215842_NB-344-Ch2-3.txt 2025-11-11_215842_NB-344-Ch2-3.png
	21:58	Start of H2 logging.	NB-344-Logger2	NB-344.ulong NB-344-Logger2-during irradiation.csv NB-344-Logger2-during irradiation.bmp
	21:58-22:08	Equilibration time.	/	/
	22:08	The irradiation was started	/	20251111_220941-after start of irradiation.jpg
	22:24	The irradiation was stopped.	/	/
	22:24-34	Equilibration time.	/	/
	22:34	Stop of O2 and H2 logging.	/	/

	22:40	Deassembling the setup, cleaning the reactor.	Tips of the sensors and reactor were covered with attached photocatalyst particles. Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode).	20251111_223525-after irradiation.jpg
--	-------	---	--	---

Analysis

Used calibration for Firesting: [20250611-BOLA-fitting-liquid-phase-trace-oxygen-sensor-H2-O2 reactor.ini](#)

Used calibration for UniSense: NB-339-Logger2

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
11.11.2025	21:53	NB-344-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-344.ulog	/	/	/	Spike in H2 signal in pre-reaction baseline, stop logging, adjusting cable position.
	21:58	NB-344-Logger2	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-344.ulog NB-344-Logger2-during irradiation.csv	NB-344-O2 and H2 curve.py	NB-344-Logger2-during irradiation.bmp NB-344-O2 and H2 curve.png	/	H2 evolution during irradiation, up to ca.124 uM at the end of irradiation
	21:00	NB-344-Ch2-1	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-11_210025_NB-344-Ch2-1.txt	/	2025-11-11_210025_NB-344-Ch2-1.png	/	Degassing of the suspension.
	21:53	NB-344-Ch2-2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-11_215318_NB-344-Ch2-2.txt	/	2025-11-11_215318_NB-344-Ch2-2.png	/	Interrupted logging to eliminate spike in H2 signal.
	21:58	NB-344-Ch2-3	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-11_215842_NB-344-Ch2-3.txt	NB-344-O2 and H2 curve.py	2025-11-11_215842_NB-344-Ch2-3.png NB-344-O2 and H2 curve.png	/	O2 evolution during irradiation, up to ca. 100 uM at the end of irradiation.

Results

Simultaneous H₂ and O₂ measurements of irradiated suspension of EA-358 (0.5 mg/mL) in O₂/H₂ photoreactor under 365 nm irradiation (150 mW/cm², 20 °C, 15 min) were performed .

H₂ level at the end of irradiation - ca. 124 uM umol/L, O₂ level - ca. 100 umol/L.

Linked experiments

SrTiO₃ - NB-316: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

SrTiO₃ - NB-318: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm², 20 °C

SrTiO₃ - NB-339: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C VI (reproduction NB-316)

Linked resources

Equipment - [VWR® VV3, Vortex Mixer](#)

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)

Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)

Protocol - [Liquid phase calibration of H₂ UniAmp sensor](#)

Protocol - [In-situ hydrogen and oxygen measurment in H₂/O₂ reactor](#)

Attached files

NB-344-SrTiO₃-photocatalytic H₂O splitting.cdxml

sha256: 066d8885f0f62fd2d240c615c34b8fc60edc7cef53430afbae30ef02502c4d70


NB-344-SrTiO₃-photocatalytic H₂O splitting.png

sha256: f7cb3993b0a0873eb4a59f6e519dda7242a17c0c757b03f75590d3ecf0dd015d



NB-344-O₂ and H₂ curve.png

sha256: e76a7c66e1a8acf4fd3847c5d577e9fff1b50ed13cfcac0c679b50ed051a4429



NB-344-O₂ and H₂ curve.py

sha256: 5da9d286ea5dfccb109f5940d283534f1c88946820b2d8d193162fb70a7563ff

20251111_220941-after start of irradiation.jpg

sha256: 9ba701cc0f7354144207875eb72e70e0e6f80cd0bfe8871fe5652ac43343af9c



20251111_223525-after irradiation.jpg

sha256: 139cb44506b7cf7d071a3c7e99e60017abb7ffd353892d5ccefb50f9cf39ce74



20251111_210354-degassing of the suspension.jpg

sha256: e2188e517aabb82b999887583f9412aeea4f6f13bffad23d51ee11a9ba3a0780



20251111_204556-suspension after vortex.jpg

sha256: d1dc8311e12eda430fa1757e87cf893901cab7380c45f89b358ed8927a07b887




NB-344-Logger2-during irradiation.csv

sha256: 11b266c01a615e321f4630950b5a48f2bed52ccfd434b852c6eb8b66656006ba

NB-344-Logger2-during irradiation.bmp

sha256: 3dc8a66a3a2552b803e4bb722db8048f646eacb3f43818b6a09f63111092075e



NB-344.ulog


sha256: 7ab4447dc6c6373bfe807f02d155bb3672417bec31fcd7ed0d120dbeb5376e5f

2025-11-11_215842_NB-344-Ch2-3.txt

sha256: 0f6955dc1f2805b361c541a7e94f92171bc365831b7bb1de7b520ecc2d3a3d81

2025-11-11_215842_NB-344-Ch2-3.png

sha256: e9e2ee71aaada35649f89aeb08abba443c015a636ee660253e029ba9b8898ee7




2025-11-11_215318_NB-344-Ch2-2.txt

sha256: 839a4659f21386e5ff2ef08cbbac88cb3001c4b53b8a097f0e08a62e690dce37

2025-11-11_215318_NB-344-Ch2-2.png

sha256: 214eeba7c3af2e746420391db069f67c11550ce5437edccca82a659798cc648e



2025-11-11_210025_NB-344-Ch2-1.png

sha256: e565ae5a610aa94e42dedd010a227b9474cb5bb8beade5601b27d01dcc5fb9e5



2025-11-11_210025_NB-344-Ch2-1.txt

sha256: 9c51132a6f26ae15364dddfef97bc561cde0b3bff8007c7279c26e59a9023e3d



Unique eLabID: 20251112-f4219ad6aa0289c0b5e5cbafb2f3845a0969a47b

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3425>