

NB-362: Gas phase H2 and O2 measurements of Al:SrTiO3 RhCrOx (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C, 70 min, degassing

Date: 2025-11-26

Tags: O2 Test Calibration Future NB
Firesting Irradiation O2 sensor H2
advanced irrad setup Unisense H2
Sensor temperature In situ Trace range
robust oxygen sensor photocatalysis

Category: SrTiO3


Status: Done

Created by: Nadzeya Brezhneva

Objectives

Simultaneous detection of O₂ and H₂ evolution in gas phase for irradiated suspension of Rh_xCrO_x:Al:SrTiO₃ suspension (EA-358, 0.5 mg/mL), 365 nm LED, 50 mW/cm², 20 °C during 70 min.

Reaction scheme



ChemDraw file linked: [NB-362-SrTiO3-photocatalytic H2O splitting.cdxml](#)

Literature/reference experiments

Literature	/
Reproduction	/
Similar experiments	SrTiO3 - NB-361: Gas phase H2 and O2 measurements of Al:SrTiO3 RhCrOx (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C, 15 min, degassing

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	pressure [bar]
milli-Q H ₂ O	/	/	/	/	/	/	/	/	25	/
Al:SrTiO3 RhCrOx (EA-358)	SrTiO3 - EA-358: Modification of Al:SrTiO3 (EA-354) via deposition of Rh, Cr oxide co-catalyst, 350°C, 1h, Upscaling (3.33x)	/	/	/	12.50	12.61	/	/	/	/

Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-314: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup](#)

	Name
Used Set-up	Equipment - Advanced irradiation setup V1.0 I
Irradiation setup number	Equipment - Irradiation setup 4 (CEEC II, E002)

	Light Source Name	Power Source Name	Wavelength [nm]	Power Setting [mW]	Analog Setting [0.00 - 10.00]
First light source	Light Source - UHP LED 365 nm-4	Power Sources - BLS-18000-1 4	365	56	0.19

Used beam combiner [Name or None]	/
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	20
Stirring speed [rpm]	500
Irradiation start: 1. Firesting [relative to start log] 2. Unisense	1. 610 s 2. 19:09:08
Irradiation stop: 1. Firesting [relative to start log] 2. Unisense	1. 4830 s 2. 20:19:29

O₂/H₂ sensor equipment

	Equipment	Used protocol
Used Firesting	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)	Protocol - Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software
Used O ₂ sensor	Equipment - Robust probe for liquid O ₂ measurment	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor
Used H ₂ sensor	Equipment - H ₂ UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor

Procedure/observations

Date	Time	Step	Observations	Pictures/Files
25.11.2025		Calibration from experiment SrTiO ₃ - NB-361: Gas phase H ₂ and O ₂ measurements of Al:SrTiO ₃ RhCrOx (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm ² , 20 °C, 15 min, degassing was used.		
		Sample preparation		
	17:50	Weighing EA-358 photocatalyst in a 50 mL vial.	Creamy solid.	/
	17:52	Addition of 25 mL H ₂ O to the vial via graduated cylinder.	/	/
	17:55-58	The suspension was vortexed for 3 min (Equipment - VWR® VV3, Vortex Mixer, stage 4/6), covered with Al foil before further use.	/	20251125_180013-suspension after vortex.jpg
		Continue in Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor from step 6		
	18:00	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	/
	18:05	Assembling the setup.	/	/
	18:16	Start of O ₂ logging.	NB-362-Ch2-1	2025-11-25_181629_NB-362-Ch2-1.txt 2025-11-25_181629_NB-362-Ch2-1.png
	18:17	The degassing was started	/	20251125_181835-degassing of the suspension.jpg

	18:54	The degassing was stopped by removing the cannula and closing the valve.	/	/
	18:55	Stop of O2 logging.	/	/
	18:58	Start of O2 logging.	NB-362-Ch2-2	2025-11-25_185858_NB-362-Ch2-2.txt 2025-11-25_185858_NB-362-Ch2-2.png
	18:58	Start of H2 logging.	NB-362-Logger1	NB-362.ulog NB-362-Logger1.csv NB-362-Logger1.bmp
	18:58-19:09	Equilibration time.	/	/
	19:09	The irradiation was started	After 20:10 the O ₂ and H ₂ curves reached plateau	20251125_191001-after start of irradiation.jpg
	20:19	The irradiation was stopped.	/	/
	20:19-30	Equilibration time.	/	/
	20:30	Stop of O2 and H2 logging.	/	/
	ca. 20:40	Deassembling the setup, cleaning the reactor.	When the front lid of the setup was opened, it was noticed that the stirring bar was not stirred properly. Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode).	/

Analysis

Used calibration for Firesting: [20250910-BOLA fitting-gas phase-4-neck photoreactor-trace oxygen robust probe-Ch2.ini](#)

Used calibration for UniSense: NB-361-Logger2

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
25.11.2025	18:58	NB-362-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-362.ulong NB-362-Logger1.csv	NB-362-O2 and H2 curve.py	NB-362-Logger1.bmp NB-362-O2 and H2 curves.png	/	H2 evolution during irradiation
	18:16	NB-362-Ch2-1	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-25_181629_NB-362-Ch2-1.txt	/	2025-11-25_181629_NB-362-Ch2-1.png	/	Degassing of the suspension.
	18:58	NB-362-Ch2-2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-25_185858_NB-362-Ch2-2.txt	NB-362-O2 and H2 curve.py	2025-11-25_185858_NB-362-Ch2-2.png NB-362-O2 and H2 curves.png	/	O2 evolution during irradiation

Results

Simultaneous H₂ and O₂ measurements (gas phase) of irradiated suspension of EA-358 (0.5 mg/mL) in O₂/H₂ photoreactor under 365 nm irradiation (50 mW/cm², 20 °C, 70 min) were performed. After 20:10 change in slope in both O2 and H2 values was observed.

After opening the lid, the problems with stirring were observed (most probably, it changed the slope in O2 and H2 curves during irradiation).

Futture recommendations

Old procedure	Problem	Suggested new procedure
/	Problems with stirring were noticed (the problem was found out after opening the lid of the setup).	Check the stirring bar or replace it with a new one if necessary.

Linked experiments

SrTiO₃ - NB-312: Gas phase H₂ and O₂ measurements with Unisense H₂ sensor, Firesting O₂ robust probe in irradiated Al:SrTiO₃ RhCrOx (NB-289, 0.5 mg/mL), 365 nm, 50 mW, 1 h, degassing (reproduction NB-304)

SrTiO₃ - EA-358: Modification of Al:SrTiO₃ (EA-354) via deposition of Rh, Cr oxide co-catalyst, 350°C, 1h, Upscaling (3.33x)

SrTiO₃ - NB-361: Gas phase H₂ and O₂ measurements of Al:SrTiO₃ RhCrOx (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C, 15 min, degassing

Linked resources

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Robust probe for liquid O2 measurment](#)

Equipment - [Advanced irradiation chamber V1.0 I](#)

Equipment - [H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)

Protocol - [Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software](#)

Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)

Protocol - [Gas phase calibration of H2 UniAmp sensor](#)

Protocol - [In-situ hydrogen and oxygen measurment in H2/O2 reactor](#)


Attached files

NB-362-O2 and H2 curve.py

sha256: 39781bd2f52857fb4c0ee55166eaa582caead52749a2b1ab734a304c6aa76c56

NB-362-O2 and H2 curves.png

sha256: 2ae68adb34943fd21ff94c6571851e8ba211157233620ddcad923ec14dec4361



NB-362-SrTiO3-photocatalytic H2O splitting.cdxml

sha256: 0ebd7d0ec0bdaf4edf2efc3efc5e472a05e8d1f4fffa835ee58268b9a6014b5b

NB-362-SrTiO3-photocatalytic H2O splitting.png

sha256: 45230fe7f6dc671214ae5b87efb2149b695488dc88aac47e2264cfd7711fd5a8



20251125_181835-degassing of the suspension.jpg

sha256: a1642b4cdb4e847f866a20ac9a3845086a6a82530c4bd0e3eb1d7ff139e5e74f



20251125_191001-after start of irradiation.jpg

sha256: 6ca3b61c7f44b2eddd58c3dea5f5ad068823440fa04114658ea942b3ba8ca86b



20251125_180013-suspension after vortex.jpg

sha256: 4e66ad9e6d4281e538ad79aacaec232dc8ee5026c67695aa4d750c3dc72824b4



2025-11-25_185858_NB-362-Ch2-2.png

sha256: df1aee26d69aa454f6245e420f0dfdb7d3283017662993d695bbb05ff25df5




2025-11-25_185858_NB-362-Ch2-2.txt

sha256: 737093603907cdf95ffc1e808122c696ab09c60e841efa77146805aa6e12519

NB-362-Logger1.bmp

sha256: 0e6be08d111cf24c7ad3f0b4ac57f59b7259c47bba8bfb86b8665d7d212722da



NB-362.ulog

sha256: 4732e6c397537e05407c518b478adbed5dbf6f981ef04a6950e6714fb0a686d6

NB-362-Logger1.csv


sha256: df699cf295ccb35b67067fc7d6e82af08fde2c3ecd520e9d691811583b47db71

2025-11-25_181629_NB-362-Ch2-1.txt

sha256: fc7f132c47e2ca1e0fb055290c868ebd387ef6eeba1813ca040e482d6e4c61df

2025-11-25_181629_NB-362-Ch2-1.png

sha256: 47e3d6dab0ca88977830b59ad08e87e324e1c1b3c718ff91e333e639e5a1abe8



Unique eLabID: 20251126-09f4bed7e97631e3105f4b7924d5a33257a1a118
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3557>