

H2/O2 reactor

Date: 2025-08-26

Tags: O2 AE H2 H2/O2 reactor


Category: Equipment

Created by: Alexander Eith

Titel

H2/O2 reactor

Table

Device name	/	
Model number	V1.0	
Manual	/	
Location	CEEC II E004	
Name in Instruments Calendar	/	
Introduction	<i>not needed in theory, but for all known use cases introduction needed, so yes</i>	
Responsible Person	Alexander Eith	
Protocol on usage of the device	Protocol - In-situ hydrogen and oxygen measurment in H2/O2 reactor	
Additional Information	<i>outer diameter 40 mm, 3 NS14 connctions on same hight, 1 NS14 connection lower than others</i> <i>Volume: Prep work - AE-529: Volume determination of H2/O2 reactor</i>	
Status from	26/08/25	

Linked experiment

Prep work - [AE-529: Volume determination of H2/O2 reactor](#)

Linked resource

Protocol - [In-situ hydrogen and oxygen measurment in H2/O2 reactor](#)

Attached files

20250826_135312.jpg

sha256: 66b791966cd6da6a228489d7ae219b37a6dfd3e3c6ea56e969d7318d27ef9069



20250826_135312.jpg

sha256: 66b791966cd6da6a228489d7ae219b37a6dfd3e3c6ea56e969d7318d27ef9069



Unique eLabID: 20250826-af67ca706fd144f5f18b3a9653a1dce52bc1fc13

Link: <https://elab.water-splitting.org/database.php?mode=view&id=260>