

# NB-337: Liquid phase H<sub>2</sub> and O<sub>2</sub> of RhCrO<sub>x</sub>,Al:SrTiO<sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm<sup>2</sup>, 20 °C (reproduction NB-325) I

Date: 2025-11-10

Tags: O<sub>2</sub> Test Calibration Future NB  
Firesting O<sub>2</sub> sensor H<sub>2</sub> SrTiO<sub>3</sub> Unisense  
RhCrO<sub>3</sub>:Al:SrTiO<sub>3</sub> H<sub>2</sub> Sensor  
temperature In situ Trace range robust  
oxygen sensor photocatalysis Unisense  
normal range

Category: SrTiO<sub>3</sub>


Status: Done

Created by: Nadzeya Brezhneva

## Objectives

Reproduction of NB-325: simultaneous detection of H<sub>2</sub> and O<sub>2</sub> evolution in liquid phase for irradiated suspension of Rh<sub>x</sub>CrO<sub>y</sub>:Al:SrTiO<sub>3</sub> suspension (EA-358, 0.5 mg/mL), 365 nm LED, 20 mW/cm<sup>2</sup>, 20 °C.

## Reaction scheme



ChemDraw file linked: [NB-325-SrTiO3-photocatalytic H2O splitting.cdxml](#)

## Literature/reference experiments

Literature	/
Reproduction	SrTiO <sub>3</sub> - NB-326: Liquid phase H <sub>2</sub> and O <sub>2</sub> of RhCrO <sub>x</sub> ,Al:SrTiO <sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm <sup>2</sup> , 20 °C (reproduction NB-325) SrTiO <sub>3</sub> - NB-325: Liquid phase H <sub>2</sub> and O <sub>2</sub> of RhCrO <sub>x</sub> ,Al:SrTiO <sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm <sup>2</sup> , 20 °C
Similar experiments	SrTiO <sub>3</sub> - NB-316: Liquid phase H <sub>2</sub> and O <sub>2</sub> of RhCrO <sub>x</sub> ,Al:SrTiO <sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm <sup>2</sup> , 20 °C

## Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass <sub>theo</sub> [mg]	Mass <sub>exp</sub> [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	Concentration [mM]
milli-Q H <sub>2</sub> O	/	/	/	/	/	/	/	0.998	25	/
Al:SrTiO <sub>3</sub> RhCrO <sub>x</sub> (EA-358)	SrTiO <sub>3</sub> - EA-358: Modification of Al:SrTiO <sub>3</sub> (EA-354) via deposition of Rh, Cr oxide co-catalyst, 350°C, 1h, Upscaling (3.33x)	/	/	/	12.50	12.68	/	/	/	/

# Excel sheet for reagent calculation

/

## Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-324: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup I I](#)

	Name
Used Set-up	<a href="#">Equipment - Advanced irradiation setup V1.0 I</a>
Irradiation setup number	<a href="#">Equipment - Irradiation setup 4 (CEEC II, E002)</a>

	Light Source Name	Power Source Name	Wavelength [nm]	Power Setting [mW]	Analog Setting [0.00 - 10.00]
First light source	<a href="#">Light Source - UHP LED 365 nm-4</a>	<a href="#">Power Sources - BLS-18000-1 4</a>	365	22	0.07

Used beam combiner [Name or None]	/
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	20
Stirring speed [rpm]	500
Irradiation start: 1. Firesting [relative to start log] 2. Unisense	1. 602 s 2. 21:41:42
Irradiation stop: 1. Firesting [relative to start log] 2. Unisense	1. 1527 s 2. 22:07:07

## O<sub>2</sub>/H<sub>2</sub> sensor equipment

	Equipment	Used protocol
Used Firesting	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)	Protocol - Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software
Used O <sub>2</sub> sensor	Equipment - Robust probe for liquid O <sub>2</sub> measurment	Protocol - In-situ hydrogen and oxygen measurment in H <sub>2</sub> /O <sub>2</sub> reactor
Used H <sub>2</sub> sensor	Equipment - H <sub>2</sub> UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurment in H <sub>2</sub> /O <sub>2</sub> reactor

## Procedure/observations

Date	Time	Step	Observations	Pictures/Files
10.11.2025		Calibration from SrTiO <sub>3</sub> - NB-336: Liquid phase H <sub>2</sub> and O <sub>2</sub> of RhCrOx,Al:SrTiO <sub>3</sub> (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm <sup>2</sup> , 20 °C V (reproduction NB-316) was used		
		<b>Sample preparation</b>		
	20:43	Weighing EA-358 photocatalyst in a 50 mL vial.	Creamy solid.	/
	20:45	Addition of 25 mL H <sub>2</sub> O to the vial via graduated cylinder.	Breaking lumps manually using Smartspatula inside the vial.	/
	20:46-49	The suspension was vortexed for 3 min (Equipment - VWR® VV3, Vortex Mixer, stage 4/6), covered with Al foil before further use.	/	20251110_204943-suspension after vortex.jpg
		Continue in Protocol - In-situ hydrogen and oxygen measurment in H <sub>2</sub> /O <sub>2</sub> reactor from step 6		
	20:50-55	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	/
	20:55-21:00	Assembling the setup.	Currently, stopper instead of H <sub>2</sub> sensor, PT100, PT1000 and O <sub>2</sub> robust probe are inside the reactor immersed in the liquid phase	/

	21:00	Start of O2 logging.	<b>NB-337-Ch2-1</b>	2025-11-10_210048_NB-337-Ch2-1.txt 2025-11-10_210048_NB-337-Ch2-1.png
	21:02	The degassing was started	/	20251110_210323-degassing of the suspension.jpg
	21:32	Cannula was transferred to gas phase, above the suspension.	/	/
	21:34	H <sub>2</sub> sensor was added in Ar counterflow.	/	/
	21:39	The degassing was stopped by removing the cannula and closing the valve.	/	/
	ca. 21:40	Stop of O2 logging.	/	/
	21:41	Start of O2 logging.	<b>NB-337-Ch2-2</b>	2025-11-10_214137_NB-337-Ch2-2.txt 2025-11-10_214137_NB-337-Ch2-2.png
	21:41	Start of H2 logging.	<b>NB-337-Logger1</b>	NB-337.ulog NB-337-Logger1.csv NB-337-Logger1.bmp
	21:41-51	Equilibration time.	/	/
	21:51	The irradiation was started	/	20251110_215343-after start of irradiation.jpg
	22:07	The irradiation was stopped.	/	/
	22:07-16	Equilibration time.	/	/
	22:16	Stop of O2 and H2 logging.	/	/

	ca. 22:20	Deassembling the setup, cleaning the reactor.	Tips of the sensors and reactor were covered with attached photocatalyst particles. Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode).	/
--	-----------	---	--	---

## Analysis

Used calibration for Firesting: [20250611-BOLA-fitting-liquid-phase-trace-oxygen-sensor-H2-O2 reactor.ini](#)

Used calibration for UniSense: NB-336-Logger2

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
10.11.2025	21:41	NB-337-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-337.ulong NB-337-Logger1.csv	NB-337-O2 and H2 curve.py	NB-337-O2 and H2 curve.png	SrTiO3 - NB-325: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm2, 20 °C	Sharp increase/decrease in H2 signal after start and stop of irradiation.
	21:00	NB-337-Ch2-1	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-10_210048_NB-337-Ch2-1.txt	/	2025-11-10_210048_NB-337-Ch2-1.png	/	Degassing of the suspension.
	21:41	NB-337-Ch2-2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-10_214137_NB-337-Ch2-2.txt	NB-337-O2 and H2 curve.py	NB-337-O2 and H2 curve.png 2025-11-10_214137_NB-337-Ch2-2.png	SrTiO3 - NB-325: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm2, 20 °C	Change in O2 value during irradiation.

## Results

Reproduction of NB-325 was performed: simultaneous H<sub>2</sub> and O<sub>2</sub> measurements of irradiated suspension of EA-358 (0.5 mg/mL) in O<sub>2</sub>/H<sub>2</sub> photoreactor under 365 nm irradiation (20 mW/cm<sup>2</sup>, 20 °C, 15 min) were performed.

Sharp increase/decrease in H2 signal after start and end of irradiation.

## Future recommendations

Old procedure	Problem	Suggested new procedure
---------------	---------	-------------------------

/	Sharp increase/decrease in H2 signal.	Adjust the position of H2 sensor.
---	---------------------------------------	-----------------------------------

## Linked experiments

SrTiO3 - NB-316: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C

SrTiO3 - NB-318: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm2, 20 °C

SrTiO3 - NB-322: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm2, 20 °C (reproduction NB-318)

SrTiO3 - NB-325: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm2, 20 °C

SrTiO3 - NB-326: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm2, 20 °C (reproduction NB-325)

SrTiO3 - NB-336: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C V (reproduction NB-316)

## Linked resources

Equipment - [VWR® VV3, Vortex Mixer](#)

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)

Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)

Protocol - [In-situ hydrogen and oxygen measurment in H2/O2 reactor](#)

## Attached files

NB-325-SrTiO3-photocatalytic H2O splitting.cdxml

sha256: 1e25f76a703054a4e3a1c95ca0b8c004e2e67c7a646758e715ac13662fa83a8c

NB-325-SrTiO3-photocatalytic H2O splitting.png

sha256: 84aba7b2eaa7e7c6c75ec4bc4edf260ee29035205566e821b254e3f00db612f8




NB-337-O2 and H2 curve.py

sha256: 66671b92adc4bf8ce737737ece90324c770119725fa7e5022c4d8bab65dbb8df


NB-337-O2 and H2 curve.png

sha256: c1abe77942b636a343fd0f4bf6eb59f05b5a2c931fc8ac0a9f99240db3ef58b5



20251110\_204943-suspension after vortex.jpg

sha256: d0c660ca98bb33962e966bd1ea55203679edc46a371f3410543b6e0e53652ae5



20251110\_210323-degassing of the suspension.jpg

sha256: 8b08fa19557ea5de2730a500e39276d3a05d630f5aa6130c99aaf5706b1c1c0a



20251110\_215343-after start of irradiation.jpg

sha256: cfb9df95414d4d3d0a6e2165cd5de642543621a86d7c783a886ac4af8a3c816f



NB-337.ulog

sha256: 3e7ed960e8f5f58fb834d5ea3eb919364c5938d855e6646b8658ab9ab636ece4

NB-337-Logger1.bmp

sha256: 49d79917497ceb3daf74d297f3015e1748fd412eee5d8daf36691cf479383053




NB-337-Logger1.csv

sha256: b72b27db88cbbde0de3c29b470e8667650d4246d71811e7e4f9cdb0dc30d9bd1

2025-11-10\_214137\_NB-337-Ch2-2.png

sha256: 11890cae9b1eb906720b5aadd5b339fad86bdc2c8a0dc377be8799007da8fe3



2025-11-10\_214137\_NB-337-Ch2-2.txt

sha256: 2b8de6fe78c9ae3907966657cdd1a239afcd710d3a97c534ab89fc502893f38c

2025-11-10\_210048\_NB-337-Ch2-1.png

sha256: 2bdfef4eef00854992de09dd4b63fd431c5cefcf45d4e2f446e94aa6273860d8



2025-11-10\_210048\_NB-337-Ch2-1.txt

sha256: 56dc9c9b92e612bd45994101a5cd6130074441db89f4bace1d00d4671382f13e



Unique eLabID: 20251110-c8a60acd2179c4733063111fc2404b72ab1e5f8b  
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3405>