

NB-322: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm², 20 °C (reproduction NB-318)

Date: 2025-10-27

Tags: O₂ Test Calibration NB Firing
O₂ sensor H₂ SrTiO₃ Unisense
RhCrO₃:Al:SrTiO₃ H₂ Sensor
temperature In situ Trace range robust
oxygen sensor photocatalysis Unisense
normal range

Category: SrTiO₃


Status: Done

Created by: Nadzeya Brezhneva

Objectives

Reproduction NB-318: Simultaneous detection of H₂ and O₂ evolution in liquid phase for irradiated suspension of Rh_xCrO_y:Al:SrTiO₃ suspension (EA-358, 0.5 mg/mL), 365 nm LED, **100 mW/cm²**, 20 °C (changing intensity relatively to reference experiment).

Reaction scheme



ChemDraw file linked: [NB-322-SrTiO3-photocatalytic H2O splitting.cdxml](#)

Literature/reference experiments

Literature	/
Reproduction	SrTiO₃ - NB-318: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm², 20 °C
Similar experiments	SrTiO₃ - NB-316: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	Concentration [mM]
milli-Q H ₂ O	/	/	/	/	/	/	/	0.998	25	/
Al:SrTiO ₃ RhCrO _x (EA-358)	SrTiO₃ - EA-358: Modification of Al:SrTiO₃ (EA-354) via deposition of Rh, Cr oxide co-catalyst, 350°C, 1h, Upscaling (3.33x)	/	/	/	12.50	12.37	/	/	/	/

Excel sheet for reagent calculation

/

Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-317: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup I](#)

	Name
Used Set-up	Equipment - Advanced irradiation setup V1.0 I
Irradiation setup number	Equipment - Irradiation setup 4 (CEEC II, E002)

	Light Source Name	Power Source Name	Wavelength [nm]	Power Setting [mW]	Analog Setting [0.00 - 10.00]
First light source	Light Source - UHP LED 365 nm-4	Power Sources - BLS-18000-1 4	365	113	0.39

Used beam combiner [Name or None]	/
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	20
Stirring speed [rpm]	500
Irradiation start: 1. Firesting [relative to start log] 2. Unisense	1. 613 s 2. 22:07:52
Irradiation stop: 1. Firesting [relative to start log] 2. Unisense	1. 1536 s 2. 22:23:16

O₂/H₂ sensor equipment

	Equipment	Used protocol
Used Firesting	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)	Protocol - Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software
Used O ₂ sensor	Equipment - Robust probe for liquid O ₂ measurment	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor
Used H ₂ sensor	Equipment - H ₂ UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor

Procedure/observations

Date	Time	Step	Observations	Pictures/Files
27.10.2025		Calibration from experiment SrTiO ₃ - NB-320: Liquid phase H ₂ and O ₂ of RhCrO _x ,Al:SrTiO ₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm ² , 20 °C (reproduction NB-316) I was used.		
		Sample preparation		
	20:25	Weighing EA-358 photocatalyst in a 50 mL vial.	Creamy solid	/
	20:27	Addition of 25 mL H ₂ O to the vial via graduated cylinder, covering with Al foil before further use.	/	/
	20:48-51	The suspension was vortexed for 3 min (Equipment - VWR® VV3, Vortex Mixer, stage 4/6), covered with Al foil before further use.	/	20251027_205218-suspension after vortex.jpg
		Continue in Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor from step 6		
	21:00	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	/
	21:05	Assembling the setup.	Currently, stopper instead of H ₂ sensor, PT100, PT1000 and O ₂ robust probe are inside the reactor immersed in the liquid phase	20251027_210553-before degassing.jpg

	21:06	Start of O2 logging.	NB-322-Ch2-1	2025-10-27_210616_NB-322-Ch2-1.txt 2025-10-27_210616_NB-322-Ch2-1.png
	21:10	The degassing was started	/	/
	21:51	Cannula was transferred to gas phase, above the suspension.	/	/
	21:53	H ₂ sensor was added in Ar counterflow.	/	/
	21:55	The degassing was stopped by removing the cannula and closing the valve.	/	/
	21:57	Stop of O2 logging.	/	/
	21:57	Start of O2 logging.	NB-322-Ch2-2	2025-10-27_215740_NB-322-Ch2-2.txt 2025-10-27_215740_NB-322-Ch2-2.png
	21:57	Start of H2 logging.	NB-322-Logger1	NB-322.ulog NB-322-Logger1.csv NB-322-Logger1.bmp
	21:57-22:08	Equilibration time.	/	20251027_220036-after introducing H2 sensor.jpg
	22:08	The irradiation was started	/	20251027_220824-after start of irradiation.jpg
	22:23	The irradiation was stopped.	/	/
	22:23-33	Equilibration time.	/	/
	22:33	Stop of O2 and H2 logging.	/	/

	ca. 22:40	Deassembling the setup, cleaning the reactor.	Tips of the sensors and reactor were covered with attached photocatalyst particles. Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode).	20251027_223530-after irradiation.jpg
--	-----------	---	--	---

Analysis

Used calibration for Firesting: [20250611-BOLA-fitting-liquid-phase-trace-oxygen-sensor-H2-O2 reactor.ini](#)

Used calibration for UniSense: NB-320-Logger6

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
27.10.2025	21:57	NB-322-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-322.ulong NB-322-Logger1.csv	NB-322-O2 and H2 curve.py	NB-322-Logger1.bmp NB-322-O2 and H2 curve.png	SrTiO3 - NB-318: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm2, 20 °C	Clean response of H2 sensor, ca. 90 μ M H2 at the end of irradiation.
	21:06	NB-322-Ch2-1	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-10-27_210616_NB-322-Ch2-1.txt	/	2025-10-27_210616_NB-322-Ch2-1.png	/	Degassing of the suspension.
	21:57	NB-322-Ch2-2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-10-27_215740_NB-322-Ch2-2.txt	NB-322-O2 and H2 curve.py	2025-10-27_215740_NB-322-Ch2-2.png NB-322-O2 and H2 curve.png	SrTiO3 - NB-318: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm2, 20 °C	O2 curve during irradiation, O2 value ca. 79 μ M at the end of irradiation.

Results

Simultaneous H₂ and O₂ measurements of irradiated suspension of EA-358 (0.5 mg/mL) in O₂/H₂ photoreactor under 365 nm irradiation (100 mW/cm², 20 °C, 15 min) were performed (reproduction NB-318).

H₂ level at the end of irradiation - 90 μ mol/L, O₂ level - 79 μ mol/L.

Linked experiments

SrTiO₃ - NB-313: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (NB-301, 0.5 mg/mL), 365 nm, 50 mW

SrTiO₃ - NB-316: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

SrTiO₃ - NB-318: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm², 20 °C

SrTiO₃ - NB-320: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C (reproduction NB-316) I

Linked resources

Equipment - [VWR® VV3, Vortex Mixer](#)

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)

Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)


Protocol - [Liquid phase calibration of H₂ UniAmp sensor](#)

Protocol - [In-situ hydrogen and oxygen measurment in H₂/O₂ reactor](#)

Attached files

NB-322-O₂ and H₂ curve.png

sha256: 930dc1109ddf62019ee30046b82569a6d6dcb9e96549b7f7eb23aa8c0d1abb95



NB-322-O₂ and H₂ curve.py

sha256: a9cddd56779a0113205ffdd530c7cfdc48c3b6a3d528534ca0dd56da5e0166df

NB-322-SrTiO₃-photocatalytic H₂O splitting.png

sha256: 7408e07503f6e62c3e9b953bcd65664c695e64e9805ca7b6def802182c54f077



NB-322-SrTiO₃-photocatalytic H₂O splitting.cdxml

sha256: 1e25f76a703054a4e3a1c95ca0b8c004e2e67c7a646758e715ac13662fa83a8c

20251027_205218-suspension after vortex.jpg

sha256: ec784f748818ef46d0c101062c7336289469caaf149c49f9eaec2e0e84b09e63



20251027_220824-after start of irradiation.jpg

sha256: d039558180791237c57b9640fb1832ad6a970bf3294851ba2f3a6ba5d836d94c



20251027_223530-after irradiation.jpg

sha256: 0a57b958a2d9743127730f2286a832408a2760ab7e52d223c650ee205796af70



20251027_220036-after introducing H2 sensor.jpg

sha256: cfa56c9d1df9c165ed0c64cb19788b4d7fd8a00446befca498f7028742534057




20251027_210553-before degassing.jpg

sha256: c4fb8a736050f8805e2a5864b92d6a126536332795ab0929f2fe74d417d0d4ac



NB-322-Logger1.bmp

sha256: d3a7a76bbdd9640487a39677d792f14335d00d4a4d01dbb8a409bdacba809f44




NB-322-Logger1.csv

sha256: 696b07f8d9ff31c9f207f9567855c43fb203a92cd486c68aef5d3d6eeaa5427d

2025-10-27_215740_NB-322-Ch2-2.png

sha256: 048cef4d2875c24589fdf1b933b5efd579e79f30277a4f4835cbdc6fce811795



2025-10-27_215740_NB-322-Ch2-2.txt

sha256: 7e4bf5dccba1069cd71adfc5fb471d592c302ec78f94cde409e02c23827fd7c6

2025-10-27_210616_NB-322-Ch2-1.png

sha256: d8f3fc3cffbcc2aacf1b72a385efe7a63199cb6ac25a44007657418b38be3345



2025-10-27_210616_NB-322-Ch2-1.txt

sha256: 8b98b0a41feeea4e8a72e84769d4521e4406c9979686284538469487f134ae9c

NB-322.ulog

sha256: e8577baf3cf70f9844623ce0ae0c90fd8223bcf6d313ee116f659b968deb62ba



Unique eLabID: 20251027-577cf695b74e1c8be49932204898d947a6c16256
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3264>