

NB-316: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

Date: 2025-10-20

Tags: O₂ Test Calibration Future NB
Firing O₂ sensor H₂ SrTiO₃
troubleshooting Unisense
RhCrO₃:Al:SrTiO₃ H₂ Sensor
temperature In situ Trace range robust
oxygen sensor photocatalysis Unisense
normal range

Category: SrTiO₃


Status: Done

Created by: Nadzeya Brezhneva

Objectives

Simultaneous detection of H₂ and O₂ evolution in liquid phase for irradiated suspension of Rh,CrO_x:Al:SrTiO₃ suspension (EA-358 sample, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C, 15 min (reference conditions).

Reaction scheme



ChemDraw file linked: [NB-316-SrTiO₃-photocatalytic H₂O splitting.cdxml](#)

Literature/reference experiments

| | |
|---------------------|--|
| Literature | / |
| Reproduction | / |
| Similar experiments | SrTiO₃ - NB-313: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (NB-301, 0.5 mg/mL), 365 nm, 50 mW SrTiO₃ - NB-315: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (NB-301, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C |

Reagents

| Name | CAS Number / Experiment Number | Inventory number | Amount [mmol] | Equivalents | Mass _{theo} [mg] | Mass _{exp} [mg] | Molar mass [g/mol] | Density (g/ml) | Volume [ml] |
|--------------------------|--------------------------------|------------------|---------------|-------------|---------------------------|--------------------------|--------------------|----------------|-------------|
| milli-Q H ₂ O | / | / | / | / | / | / | / | 0.998 | 25.00 |

| | | | | | | | | | |
|---------------------------------|---|---|---|---|------|-------|---|---|---|
| Al:SrTiO3 RhCrOx (EA-358) | SrTiO3 - EA-358: Modification of Al:SrTiO3 (EA-354) via deposition of Rh, Cr oxide co-catalyst, 350°C, 1h, Upscaling (3.33x) | / | / | / | 12.5 | 12.86 | / | / | / |
|---------------------------------|---|---|---|---|------|-------|---|---|---|

Excel sheet for reagent calculation

/

Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-314: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup](#)

| | |
|--------------------------|---|
| | Name |
| Used Set-up | Equipment - Advanced irradiation setup V1.0 I |
| Irradiation setup number | Equipment - Irradiation setup 4 (CEEC II, E002) |

| | Light Source Name | Power Source Name | Wavelength [nm] | Power Setting [mW] | Analog Setting [0.00 - 10.00] |
|--------------------|---|---|-----------------|--------------------|-------------------------------|
| First light source | Light Source - UHP LED 365 nm-4 | Power Sources - BLS-18000-1 4 | 365 | 56 | 0.19 |

| | |
|-----------------------------------|-----|
| Used beam combiner [Name or None] | / |
| Irradiation distance [cm] | 6.5 |
| Thermostat temperature [°C] | 20 |
| Stirring speed [rpm] | 500 |

| | |
|---|--------------------------|
| Irradiation start: 1. Firing [relative to start log] 2. Unisense | 1. 605 s 2. 21:15:23 |
| Irradiation stop: 1. Firing [relative to start log] 2. Unisense | 1. 1524 s 2. 21:30:42 |

O₂/H₂ sensor equipment

| | Equipment | Used protocol |
|----------------------------|--|--|
| Used Firing | Equipment - Firing Fiber-Optic Oxygen Meter 2 Channel (Firing 2) | Protocol - Operation of Firing Fiber-Optic Oxygen Meter 2 Channel Software |
| Used O ₂ sensor | Equipment - Robust probe for liquid O ₂ measurement | Protocol - In-situ hydrogen and oxygen measurement in H ₂ /O ₂ reactor |
| Used H ₂ sensor | Equipment - H ₂ UniAmp Sensor - Normal range - 2.1 x 80 mm needle | Protocol - In-situ hydrogen and oxygen measurement in H ₂ /O ₂ reactor |

Procedure/observations

| Date | Time | Step | Observations | Pictures/Files |
|------------|-----------|--|--------------|---|
| 20.10.2025 | | Calibration from SrTiO ₃ - NB-315: Liquid phase H ₂ and O ₂ of RhCrOx,Al:SrTiO ₃ (NB-301, 0.5 mg/mL), 365 nm, 50 mW/cm ² , 20 °C was used | | |
| | | Sample preparation | | |
| | 19:05 | Weighing EA-358 photocatalyst in a 50 mL vial. | Creamy solid | 20251020_190334-weighed solid.jpg |
| | ca. 19:08 | Addition of 25 mL H ₂ O to the vial via graduated cylinder. | / | / |
| | 19:10-13 | The suspension was vortexed for 3 min (Equipment - VWR® VV3, Vortex Mixer, stage 4/6), covered with Al foil before further use. | / | 20251020_191321-suspension after vortex.jpg |
| | | Continue in Protocol - In-situ hydrogen and oxygen measurement in H ₂ /O ₂ reactor from step 6 | | |
| | 19:30 | The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) . | / | / |

| | | | | |
|--|-----------|--|--|--|
| | 19:35 | Assembling the setup. | Currently, stopper instead of H ₂ sensor, PT100, PT1000 and O ₂ robust probe are inside the reactor immersed in the liquid phase | 20251020_193516-before degassing and irradiation.jpg |
| | 19:39 | Start of O2 logging. | NB-316-Ch2-1 | 2025-10-20_193930_NB-316-Ch2-1.txt 2025-10-20_193930_NB-316-Ch2-1.png |
| | 19:48 | The degassing was started | / | / |
| | ca. 20:20 | Cannula was transferred to gas phase, above the suspension. | / | / |
| | 20:21 | H ₂ sensor was added in Ar counterflow. | / | 20251020_202318-introducing H2 sensor.jpg |
| | 20:26 | The degassing was stopped by removing the cannula and closing the valve. | / | / |
| | 20:29 | Stop of O2 logging. | / | / |
| | 20:30 | Start of H2 logging. | NB-316-Logger1 | NB-316.ulong NB-316-Logger1.bmp |
| | 20:31 | Start of O2 logging. | NB-316-Ch2-2 | 2025-10-20_203104_NB-316-Ch2-2.txt 2025-10-20_203104_NB-316-Ch2-2.png |
| | 20:35 | Stop O2 logging. | Too high leakage rate - 20:36 - 4.58 uM O ₂ | / |
| | 20:42 | Stop H2 logging. | Problem that caused leakage - BOLA fitting for PT100 was not tight, Deassembling the setup, eliminating problem --> solved | / |
| | 20:46 | Restart of degassing. | For degassing procedure H ₂ sensor was removed from the reactor, it was replaced with NS14 stopper | / |
| | 20:59 | Introducing H2 sensor under Ar flow (cannula was in gas phase above the liquid). | / | / |
| | ca. 21:03 | Removing cannula, closing the valve. | / | / |
| | 21:05 | Start of O2 logging. | NB-316-Ch2-3 | 2025-10-20_210518_NB-316-Ch2-3.txt 2025-10-20_210518_NB-316-Ch2-3.png |
| | 21:05 | Start of H2 logging. | NB-316-Logger2 | NB-316.ulong NB-316-Logger2-during irradiation.csv NB-316-Logger2-during irradiation.bmp |

| | | | | |
|--|-----------|---|--|--|
| | 21:05-15 | Equilibration time. | / | / |
| | 21:15 | The irradiation was started | / | 20251020_211617-after start of irradiation.jpg |
| | 21:30 | The irradiation was stopped. | / | / |
| | 21:30-41 | Equilibration time. | / | / |
| | 21:41 | Stop of O2 and H2 logging. | / | / |
| | ca. 21:45 | Deassembling the setup, cleaning the reactor. | Tips of the sensors and reactor were covered with attached photocatalyst particles. Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode). | 20251020_214330-after irradiation.jpg 20251020_214344-after irradiation-closer view.jpg |

Analysis

Used calibration for Firesting: [20250611-BOLA-fitting-liquid-phase-trace-oxygen-sensor-H2-O2 reactor.ini](#)

Used calibration for UniSense: NB-315-Logger3

| Date | Time | Sample name | Analysis method | Analytical device | Solvent | Raw Data | Python script | Processed Data | Comparative Data | Interpretation |
|------------|-------|----------------|------------------------------|--|---------|--|-------------------------------------|--|------------------|---|
| 20.10.2025 | 20:30 | NB-316-Logger1 | electrochemical H2 detection | Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle | water | NB-316.ulog | / | NB-316-Logger1.bmp | / | High O2 leakage rate was found - logging was stopped, needs to be redone. |
| | 21:05 | NB-316-Logger2 | electrochemical H2 detection | Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle | water | NB-316.ulog NB-316-Logger2-during irradiation.csv | NB-316-O2 and H2.py | NB-316-Logger2-during irradiation.bmp NB-316-H2 and O2 curves.png | / | Clean H2 response, H2 value ca. 52 uM at the end of irradiation. |
| | 19:39 | NB-316-Ch2-1 | Optical O2 detection | Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel | water | 2025-10-20_193930_NB-316-Ch2-1.txt | / | 2025-10-20_193930_NB-316-Ch2-1.png | / | Degassing of the suspension. |

| | | | | | | | | | |
|-------|--------------|----------------------|--|-------|------------------------------------|---------------------|---|---|--|
| 20:31 | NB-316-Ch2-2 | Optical O2 detection | Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel | water | 2025-10-20_203104_NB-316-Ch2-2.txt | / | 2025-10-20_203104_NB-316-Ch2-2.png | / | High O2 leakage rate. |
| 21:05 | NB-316-Ch2-3 | Optical O2 detection | Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel | water | 2025-10-20_210518_NB-316-Ch2-3.txt | NB-316-O2 and H2.py | 2025-10-20_210518_NB-316-Ch2-3.png NB-316-H2 and O2 curves.png | / | Change in O2 value during irradiation, O2 value ca. 40 uM at the end of irradiation. |

Results

Simultaneous H₂ and O₂ measurements of irradiated suspension of EA-358 (0.5 mg/mL) in O₂/H₂ photoreactor under 365 nm irradiation (50 mW/cm², 20 °C, 15 min, reference conditions) were performed.

H₂ level at the end of irradiation - 52 umol/L, O₂ level - 40 umol/L.

Future recommendations

| Old procedure | Problem | Suggested new procedure |
|---------------|---------------------------------------|--|
| / | High O2 leakage rate after degassing. | Check BOLA fitting to each sensor before experiment, should fit tightly. |

Linked experiments

SrTiO₃ - [NB-313: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 \(NB-301, 0.5 mg/mL\), 365 nm, 50 mW](#)

SrTiO₃ - [NB-315: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 \(NB-301, 0.5 mg/mL\), 365 nm, 50 mW/cm², 20 °C](#)

Linked resources

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)

Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)

Protocol - [Liquid phase calibration of H2 UniAmp sensor](#)


Protocol - [In-situ hydrogen and oxygen measurment in H2/O2 reactor](#)

Attached files

NB-316-O2 and H2.py
sha256: ba3195fc32aa9261b7ae02f4ac926dd5bf69825c6dc1c9c137fd3283a6672193

NB-316-H2 and O2 curves.png

sha256: c5925f567c35b5673158cdfe50114a41f6e58254e1472bc684d484cb89a13bfe



Unisense-NB-316-Screenshot 2025-11-07 083140.png

sha256: e3b3adce7a35ccee765a73e3bbc60a2d1e17184be789033010065d04880da780



NB-316-SrTiO3-photocatalytic H2O splitting.png

sha256: f9541119d3b1697940988dad6b96e3acc271559fdf00a3ad36913aa8e3309d61



NB-316-SrTiO3-photocatalytic H2O splitting.cdxml


sha256: 492252c43cd5886d865bfe0134680d27dce9e3bbfbf36d2890b07a732289313c

NB-316-Logger2-during irradiation.csv

sha256: 584b70d220a548571176fcc50ec042c82db9c41b01271091b8b6a3a4e73dc52f

NB-316-Logger2-during irradiation.bmp

sha256: 46e84a1bf03f7f2329b7132e0115f00c8f5b397d768c96d5d7679007f26ce34f



20251020_214344-after irradiation-closer view.jpg

sha256: 6f3f8dbe4173771b997881df1882ffc75994278ede43e3f4876f9c7509605f5



20251020_214330-after irradiation.jpg

sha256: 60b9842c24fd7f18a3cf8ef02f86e90c33d6493a7552829b41a7a1ccc166e7ef



20251020_202318-introducing H2 sensor.jpg

sha256: 79ae5b5ef12b6094a06644d795ce9fb40c40fa44696be1a445cfe34c2bf65b0c



20251020_211617-after start of irradiation.jpg

sha256: 8e09dcb02772d9203e1b3b9cc2e83050c4c9022875ae83596d5d16dd268fefb0



20251020_193516-before degassing and irradiation.jpg

sha256: dba50468ab236eb833f172973da212525cba9bff59b47da010f6d72281a69ac7



20251020_191321-suspension after vortex.jpg

sha256: 4285bc1fe5be4cff80adadf3abed434934f19b3f9a2a31581708eaefe1737238



20251020_190334-weighed solid.jpg

sha256: d1050f4236d68eebf8f64019e2ef2dab713f85bb0b938e3340aa2d7dd8ec0304




2025-10-20_210518_NB-316-Ch2-3.txt

sha256: ca813324d36186e99c665a2dd771cd90777fb7738642e2d2eb5312a60355a5d9

2025-10-20_210518_NB-316-Ch2-3.png

sha256: 2415c3aa2680ca667c3f27fab6b190f82192a9fa7a097cee1f88d5393f9fa361




2025-10-20_203104_NB-316-Ch2-2.txt

sha256: ff763d076d3a6cf949ff737b7504bce9107df4345529ee7fbceaae7f6fbb5e92

2025-10-20_203104_NB-316-Ch2-2.png

sha256: 91a71707699b975617e8c20c403b5364edc5795bd8e5ce7cc7929d43c0be0aac




2025-10-20_193930_NB-316-Ch2-1.txt

sha256: 7a451240f18f8c6b3a9e3865e4a72e079545a847e6ace576ac7eca48c8bbe01c

2025-10-20_193930_NB-316-Ch2-1.png

sha256: 5727296e6a1e9b04cde57dbec77e07436d12dc44fb65573394cc737e4a291b6f



NB-316-Logger1.bmp

sha256: f8e3d83369aca084df2ba11f7c0777d6bc9e5a232e06e284aabbcd11c8397a5a



NB-316.ulog

sha256: 369162fb004bd173d2c95039c0b837b6410ef9bc6b8253ffc0125b033a1fd315



Unique eLabID: 20251020-447a9663b9794ef9311885a8d19530c6afee6244
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3238>