

NB-360: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 10 °C (reproduction NB-351) I

Date: 2025-11-24

Tags: O₂ Test Calibration NB Firing
O₂ sensor H₂ SrTiO₃ Unisense
RhCrO₃:Al:SrTiO₃ H₂ Sensor
temperature In situ H₂/O₂ reactor Trace
range robust oxygen sensor
photocatalysis Unisense normal range

Category: SrTiO₃


Status: Done

Created by: Nadzeya Brezhneva

Objectives

Second reproduction of NB-359: simultaneous detection of H₂ and O₂ evolution in liquid phase for irradiated suspension of Rh_xCrO_y:Al:SrTiO₃ suspension (EA-358, 0.5 mg/mL), 365 nm LED, 50 mW/cm², **10 °C** (changing temperature relatively to reference conditions).

Reaction scheme



ChemDraw file linked: [NB-351-SrTiO3-photocatalytic H2O splitting.cdxml](#)

Literature/reference experiments

Literature	/
Reproduction	SrTiO₃ - NB-351: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 10 °C SrTiO₃ - NB-359: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 10 °C (reproduction NB-351)
Similar experiments	SrTiO₃ - NB-316: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	Concentration [mM]
milli-Q H ₂ O	/	/	/	/	/	/	/	0.998	2	/
Al:SrTiO ₃ RhCrO _x (EA-358)	SrTiO₃ - EA-358: Modification of Al:SrTiO₃ (EA-354) via deposition of Rh, Cr oxide co-catalyst, 350°C, 1h, Upscaling (3.33x)	/	/	/	12.50	12.48	/	/	/	/

Excel sheet for reagent calculation

/

Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-314: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup](#)

	Name
Used Set-up	Equipment - Advanced irradiation setup V1.0 I
Irradiation setup number	Equipment - Irradiation setup 4 (CEEC II, E002)

	Light Source Name	Power Source Name	Wavelength [nm]	Power Setting [mW]	Analog Setting [0.00 - 10.00]
First light source	Light Source - UHP LED 365 nm-4	Power Sources - BLS-18000-1 4	365	56	0.19

Used beam combiner [Name or None]	/
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	10
Stirring speed [rpm]	500
Irradiation start: 1. Firesting [relative to start log] 2. Unisense	1. 604 s 2. 20:08:24
Irradiation stop: 1. Firesting [relative to start log] 2. Unisense	1. 1535 s 2. 20:23:56

O₂/H₂ sensor equipment

	Equipment	Used protocol
Used Firesting	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)	Protocol - Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software
Used O ₂ sensor	Equipment - Robust probe for liquid O ₂ measurment	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor
Used H ₂ sensor	Equipment - H ₂ UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor

Procedure/observations

Date	Time	Step	Observations	Pictures/Files
24.11.2025		Calibration from experiment SrTiO₃ - NB-359 : Liquid phase H ₂ and O ₂ of RhCrO _x ,Al:SrTiO ₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm ² , 10 °C (reproduction NB-351) was used.		
		Sample preparation		
	18:45	Weighing EA-358 photocatalyst in a 50 mL vial.	Creamy solid.	/
	18:49	Addition of 25 mL H ₂ O to the vial via graduated cylinder.	/	/
	18:50-53	The suspension was vortexed for 3 min (Equipment - VWR® VV3, Vortex Mixer , stage 4/6), covered with Al foil before further use.	/	20251124_185332-suspension after vortex.jpg
		Continue in Protocol - In-situ hydrogen and oxygen measurment in H₂/O₂ reactor from step 6		
	18:55	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	/
	19:00	Assembling the setup.	Currently, stopper instead of H ₂ sensor, PT100, PT1000 and O ₂ robust probe are inside the reactor immersed in the liquid phase	/
	19:00-11	Keeping suspension at 10 °C, T control by PT1000.	/	/
	19:12	Start of O₂ logging.	NB-360-Ch2-1	2025-11-24_191210_NB-360-Ch2-1.txt 2025-11-24_191210_NB-360-Ch2-1.png

	19:15	The degassing was started	/	20251124_191601-degassing of the suspension.jpg
	19:47	Cannula was transferred to gas phase, above the suspension.	/	/
	19:49	H ₂ sensor was added in Ar counterflow.	/	/
	19:52	The degassing was stopped by removing the cannula and closing the valve.	offset for H2 sensor - 1mV	/
	19:58	Stop of O2 logging.	/	/
	19:58	Start of O2 logging.	NB-360-Ch2-2	2025-11-24_195821_NB-360-Ch2-2.txt 2025-11-24_195821_NB-360-Ch2-2.png
	19:58	Start of H2 logging.	NB-360-Logger1	NB-360.ulog NB-360-Logger1.csv NB-360-Logger1.bmp
	19:58-20:08	Equilibration time.	/	/
	20:08	The irradiation was started	/	20251124_200900-after start of irradiation.jpg
	20:23	The irradiation was stopped.	/	/
	20:23-33	Equilibration time.	/	/
	20:33	Stop of O2 and H2 logging.	/	/
	ca. 20.40	Deassembling the setup, cleaning the reactor.	Tips of the sensors and reactor were covered with attached photocatalyst particles. Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode).	20251124_204350-after irradiation.jpg

Analysis

Used calibration for Firesting: [20250611-BOLA-fitting-liquid-phase-trace-oxygen-sensor-H2-O2 reactor.ini](#)

Used calibration for UniSense: NB-359-Logger2

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
24.11.2025	19:58	NB-360-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-360.ulong NB-360-Logger1.csv	NB-360-O2 and H2 curve.py	NB-360-Logger1.bmp	SrTiO3 - NB-351: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 10 °C	H2 evolution during irradiation, weird increase in signal during irradiation, change in baseline.
	19:12	NB-360-Ch2-1	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-24_191210_NB-360-Ch2-1.txt	/	2025-11-24_191210_NB-360-Ch2-1.png	/	Degassing of the suspension.
	19:58	NB-360-Ch2-2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-24_195821_NB-360-Ch2-2.txt	NB-360-O2 and H2 curve.py	2025-11-24_195821_NB-360-Ch2-2.png	SrTiO3 - NB-351: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm2, 10 °C	O ₂ evolution during irradiation.

Results

Second reproduction of NB-351, Simultaneous H₂ and O₂ measurements of irradiated suspension of EA-358 (0.5 mg/mL) in O₂/H₂ photoreactor under 365 nm irradiation (50 mW/cm², 10 °C, 15 min) were performed. Strong increase in H2 signal during irradiation --> data not applicable for analysis.

Linked experiments

SrTiO₃ - NB-316: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

SrTiO₃ - NB-318: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm², 20 °C

SrTiO₃ - NB-322: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm², 20 °C (reproduction NB-318)

SrTiO₃ - NB-325: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm², 20 °C

SrTiO₃ - NB-330: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 30 °C

SrTiO₃ - NB-334: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 30 °C (reproduction NB-330)

SrTiO₃ - NB-351: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 10 °C

SrTiO₃ - NB-359: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 10 °C (reproduction NB-351)

Linked resources

Equipment - [VWR® VV3, Vortex Mixer](#)

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Advanced power measurment chamber V1.0 I](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)

Light Source - [UHP LED 365 nm-4](#)

Power Sources - [BLS-18000-1 4](#)

Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)

Protocol - [Liquid phase calibration of H₂ UniAmp sensor](#)

Protocol - [In-situ hydrogen and oxygen measurment in H₂/O₂ reactor](#)


Attached files

NB-360-O2 and H2 curve.py

sha256: c8f782a92a3c63c2ebe3d7974aed7d3e83e2c4ef36fa1284790fce888e96bf3f

NB-360-O2 and H2 curves.png

sha256: 4dcbc875b97d1f6531414166460aa31f8e381d91b52426261b1583666c5227ca



20251124_200900-after start of irradiation.jpg

sha256: 03a6ecad5ea7084097a1e8a27555eb1b915b7e703f4c26752b9ec25b4d8b3a02




20251124_204350-after irradiation.jpg

sha256: b1a827bd9db8f5135d75c6244e0e4645ef56b299f28cc99d52e2206ec7d99509



20251124_185332-suspension after vortex.jpg

sha256: 22b00d5aeb623b989e53e6184631f2deff1e73dcb66bbbf77e22d5518a38a90



20251124_191601-degassing of the suspension.jpg

sha256: 96a98a56a6231d9a73d09a31e8cb9c0e93c1c8d3fef1a2055f583580319e9a02



NB-360-Logger1.csv


sha256: e52cab7c87910caefc720300069bee65c846920d6feee073ed1cfa3f508e9814

NB-360.ulog

sha256: 7403996e3d592b28a1d414eac39bb9bbbd2a5eb90ef4a249e7745e733c42051c

NB-360-Logger1.bmp

sha256: 347b2c9803e88d01500973ead0f98e6bf22b6effaf945e6e4674e991c91699f7




2025-11-24_195821_NB-360-Ch2-2.txt

sha256: 1a7580759adf5354ce949101c4d0775bcb77d32e145d0bac94c08b875ea21fcc

2025-11-24_195821_NB-360-Ch2-2.png

sha256: a2d881571a2be1dd984f0a2cc199050f064ed06f71eb4451dbcc5bd292ec283b



2025-11-24_191210_NB-360-Ch2-1.png

sha256: 1bedfbcebb647ec0af86dd74fb330c1d30cd95b51b17e6f381b17835f94c576f



2025-11-24_191210_NB-360-Ch2-1.txt

sha256: 84726997e332984dc33165f5a4e5408027bcbe2f094b2566e4eb128626152162



Unique eLabID: 20251124-06d37f1eb406b35265537de84353999dc98f845f
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3546>