

NB-357: Liquid phase H₂ and O₂ of RhCrO_x,Al:SrTiO₃ (NB-341, 0.05 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C (reproduction NB-354)

Date: 2025-11-21

Tags: [O2](#) [Test](#) [Calibration](#) [NB](#) [Firing](#) [O2 sensor](#) [H2](#) [SrTiO3](#) [Unisense](#) [RhCrO3:Al:SrTiO3](#) [H2 Sensor](#) [temperature](#) [In situ](#) [Trace range robust](#) [oxygen sensor](#) [photocatalysis](#) [Unisense](#) [normal range](#)

Category: SrTiO₃


Status: Done

Created by: Nadzeya Brezhneva

Objectives

Reproduction of NB-354: simultaneous detection of H₂ and O₂ evolution in liquid phase for irradiated suspension of Rh_xCr_yO₃:Al:SrTiO₃ suspension (NB-341, 0.05 wt% Rh,Cr, 0.5 mg/mL), 365 nm LED, 50 mW/cm², 20 °C.

Reaction scheme



ChemDraw file linked: [NB-327-SrTiO3-photocatalytic H2O splitting.cdxml](#)

Literature/reference experiments

Literature	/
Reproduction	/
Similar experiments	SrTiO3 - NB-327: Liquid phase H2 and O2 of RhCrO_x,Al:SrTiO₃ (NB-321, 0.05 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C SrTiO3 - NB-316: Liquid phase H2 and O2 of RhCrO_x,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	Concentration [mM]
milli-Q H ₂ O	/	/	/	/	/	/	/	0.998	25.00	/

Al:SrTiO3 RhCrOx (NB-341)	SrTiO3 - NB-341: Modification of EA-354 (SrTiO3:Al, upscaled batch) with Rh, Cr oxide cocatalyst (0.05 wt%), fresh solutions of RhCl3 and Cr(NO3)3	/	/	/	12.50	12.50	/	/	/	/
---------------------------------	--	---	---	---	-------	-------	---	---	---	---

Excel sheet for reagent calculation

/

Irradiation Parameters

Power measurement was performed using [Power Meter - 843-R-USB + 919P-020-12](#) in [Equipment - Advanced power measurment setup V1.0 I](#)

Power measurement was performed in experiment [Prep work - NB-314: Measuring power output of UHP-365 nm #4 with 18A-4 in advanced irradiation setup](#)

	Name
Used Set-up	Equipment - Advanced irradiation setup V1.0 I
Irradiation setup number	Equipment - Irradiation setup 4 (CEEC II, E002)

	Light Source Name	Power Source Name	Wavelength [nm]	Power Setting [mW]	Analog Setting [0.00 - 10.00]
First light source	Light Source - UHP LED 365 nm-4	Power Sources - BLS-18000-1 4	365	56	0.19

Used beam combiner [Name or None]	/
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	20
Stirring speed [rpm]	500

Irradiation start: 1. Firesting [relative to start log] 2. Unisense	1. 605 s 2. 2:27:16
Irradiation stop: 1. Firesting [relative to start log] 2. Unisense	1. 1515 s 2. 2:42:26

O₂/H₂ sensor equipment

	Equipment	Used protocol
Used Firesting	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel (Firesting 2)	Protocol - Operation of Firesting Fiber-Optic Oxygen Meter 2 Channel Software
Used O ₂ sensor	Equipment - Robust probe for liquid O ₂ measurment	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor
Used H ₂ sensor	Equipment - H ₂ UniAmp Sensor - Normal range - 2.1 x 80 mm needle	Protocol - In-situ hydrogen and oxygen measurment in H ₂ /O ₂ reactor

Procedure/observations

Date	Time	Step	Observations	Pictures/Files
		Calibration from SrTiO ₃ - NB-356: Liquid phase H ₂ and O ₂ of RhCrO _x ,Al:SrTiO ₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm ² , 20 °C IX (reproduction NB-316) was used.	/	/
21.11.2025		Sample preparation		
	1:05	Weighing photocatalyst in a 50 mL vial.	Slightly creamy solid	/
		Important note: before weighing, the solid was dispersed inside the vial manually with Smartspatula - to break lumps manually, but to avoid mortaring.		
	1:06	Addition of 25 mL H ₂ O to the vial via graduated cylinder.	/	/

1:09-12	The suspension was vortexed for 3 min (Equipment - VWR® VV3, Vortex Mixer , stage 4/6), covered with Al foil before further use.	/	20251121_011342-suspension after vortex.jpg
	Continue in Protocol - In-situ hydrogen and oxygen measurment in H2/O2 reactor from step 6		
1:15	The suspension was transferred to the reactor using glass pipette (preliminary the vial was manually shaken ca. 15 s) .	/	/
1:20	Assembling the setup.	Currently, stopper instead of H ₂ sensor, PT100, PT1000 and O ₂ robust probe are inside the reactor immersed in the liquid phase	/
1:29	Start of O2 logging.	NB-357-Ch2-1	2025-11-21_012927_NB-357-Ch2-1.txt 2025-11-21_012927_NB-357-Ch2-1.png
1:32	The degassing was started	/	20251121_013310-degassing of the suspension.jpg
2:07	Cannula was transferred to gas phase, above the suspension.	/	/
2:10	H ₂ sensor was added in Ar counterflow.	/	/
2:13	The degassing was stopped by removing the cannula and closing the valve.	/	/
2:17	Stop of O2 logging.	/	/
2:17	Start of O2 logging.	NB-357-Ch2-2	2025-11-21_021711_NB-357-Ch2-2.txt 2025-11-21_021711_NB-357-Ch2-2.png
2:17	Start of H2 logging.	NB-357-Logger1	NB-357.ulong NB-357-Logger1.csv NB-357-Logger1.bmp
2:17-27	Equilibration time.	/	/
2:27	The irradiation was started	/	20251121_022807-after start of irradiation.jpg
2:42	The irradiation was stopped.	/	/

	2:42-52	Equilibration time.	/	/
	2:52	Stop of O2 and H2 logging.	/	/
	ca. 3:00	Deassembling the setup, cleaning the reactor.	Tips of the sensors and reactor were covered with attached photocatalyst particles, bubbles attached to the surface of the sensors tips. Tip: After preliminary cleaning with sticks, wipes, the residual particles attached to the walls of the reactor could be removed by sonication - fill the reactor with water and place it in ultrasonic bath for ca. 20 s (Eco mode).	20251121_025513-after irradiation.jpg

Analysis

Used calibration for Firesting: [20250611-BOLA-fitting-liquid-phase-trace-oxygen-sensor-H2-O2 reactor.ini](#)

Used calibration for UniSense: NB-356-Logger2

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
21.11.2025	2:17	NB-357-Logger1	electrochemical H2 detection	Equipment - H2 UniAmp Sensor - Normal range - 2.1 x 80 mm needle	water	NB-357.ulog NB-357-Logger1.csv	NB-357-O2 and H2 curve.py	NB-357-Logger1.bmp NB-357-O2 and H2 curves.png	SrTiO3 - NB-354: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (NB-341, 0.05 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C	H2 evolution during irradiation.
	1:29	NB-357-Ch2-1	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-21_012927_NB-357-Ch2-1.txt	/	2025-11-21_012927_NB-357-Ch2-1.png	/	Degassing of the suspension.
	2:17	NB-357-Ch2-2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	water	2025-11-21_021711_NB-357-Ch2-2.txt	NB-357-O2 and H2 curve.py	2025-11-21_021711_NB-357-Ch2-2.png NB-357-O2 and H2 curves.png	SrTiO3 - NB-354: Liquid phase H2 and O2 of RhCrOx,Al:SrTiO3 (NB-341, 0.05 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm2, 20 °C	O2 evolution during irradiation

Results

Reproduction of NB-354: simultaneous H₂ and O₂ measurements of irradiated suspension of NB-341 (0.05 wt% RhCr, 0.5 mg/mL) in O₂/H₂ photoreactor under 365 nm irradiation (50 mW/cm², 20 °C, 15 min) were performed.

Linked experiments

SrTiO₃ - NB-316: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

SrTiO₃ - NB-318: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm², 20 °C

SrTiO₃ - NB-321: Modification of EA-354 (SrTiO₃:Al, upscaled batch) with Rh, Cr oxide cocatalyst (0.05 wt%)

SrTiO₃ - NB-322: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 100 mW/cm², 20 °C (reproduction NB-318)

SrTiO₃ - NB-325: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm², 20 °C

SrTiO₃ - NB-326: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 20 mW/cm², 20 °C (reproduction NB-325)

SrTiO₃ - NB-327: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (NB-321, 0.05 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

SrTiO₃ - NB-331: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C (reproduction NB-316) I I I

SrTiO₃ - NB-332: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (NB-321, 0.05 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C (reproduction NB-327)

SrTiO₃ - NB-341: Modification of EA-354 (SrTiO₃:Al, upscaled batch) with Rh, Cr oxide cocatalyst (0.05 wt%), fresh solutions of RhCl₃ and Cr(NO₃)₃

SrTiO₃ - NB-354: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (NB-341, 0.05 wt% Rh, Cr, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C

SrTiO₃ - NB-356: Liquid phase H₂ and O₂ of RhCrOx,Al:SrTiO₃ (EA-358, 0.5 mg/mL), 365 nm, 50 mW/cm², 20 °C IX (reproduction NB-316)

Linked resources

Equipment - [VWR® VV3, Vortex Mixer](#)

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 2\)](#)

Equipment - [Irradiation setup 4 \(CEEC II, E002\)](#)

Protocol - [Getting hydrogen from hydrogen bottle in CEEC II E014](#)

Protocol - [Liquid phase calibration of H2 UniAmp sensor](#)

Protocol - [In-situ hydrogen and oxygen measurment in H2/O2 reactor](#)


Attached files

NB-357-O2 and H2 curve.py

sha256: ab7d6282ebef9e8e357affcbe44bf93e8a681aeb8198d9a8378988d34b86e12b

NB-357-O2 and H2 curves.png

sha256: c923b1633996e7e9371284955b2e533b8e62d03b2521ab649f639504c72a6f4f



20251121_013310-degassing of the suspension.jpg

sha256: 4c1a4c9cee193d3f89a569986594e7916683ecaa6f72447b9118618760519443



20251121_022807-after start of irradiation.jpg

sha256: ed3910537e8d3c27f6c3cab6ab21580079a46b70235fa17aaba8191a6e3c43ec



20251121_025513-after irradiation.jpg

sha256: 83e5c580f9469c5846260c16f6c32007060a04a8f074b0e244cb2ef09100cf3b




20251121_011342-suspension after vortex.jpg

sha256: c003d3d2b7f7b2996f9c708902a0d808bc3e35bda6c786a34547017eba6b897c



2025-11-21_021711_NB-357-Ch2-2.png

sha256: f77d3ef1187e59f5032e53cedd961465c07b3d0fd1cfc8700de1b8e8c55f403a



2025-11-21_021711_NB-357-Ch2-2.txt

sha256: 93f4a60a4c1233b4c85814d40be9be4af237782aa5f956789b67532e66dc201b

2025-11-21_012927_NB-357-Ch2-1.txt

sha256: 7f93d357eb083c2f34ca54745cf4fd9bf434b853b24516d130d0620471eec874


2025-11-21_012927_NB-357-Ch2-1.png

sha256: 4841957a4151d08424d465d553abb2695d0d3679058486a578de01424bd7291c



NB-357-Logger1.bmp

sha256: 530abc1664d2ed3b5dbeac25bf160575bc1cec732527a944c8969ed439c7d375



NB-357.ulog

sha256: 22f983992f4776bd6d5c928a18f94b4290c74a5055ef1b69c2bda813aed6a68f

NB-357-Logger1.csv

sha256: 8e97c71f4ecac0b53ef8e69de843af5db12baaf508cc9e4b97660b3f94516fb6



Unique eLabID: 20251121-e353e205b4bf5ac34b1eacb79a36f6e9f20c17d9

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3522>