

# NB-342: Modification of EA-354 (SrTiO<sub>3</sub>:Al, upscaled batch) with Rh, Cr oxide cocatalyst (0.2 wt%), fresh solutions of RhCl<sub>3</sub>\*3H<sub>2</sub>O and Cr(NO<sub>3</sub>)<sub>3</sub>\*9H<sub>2</sub>O

Date: 2025-11-11


Tags: NB synthesis Furnace Muffle  
Furnace SrTiO<sub>3</sub> RhCl<sub>3</sub> Osterloh  
RhxCr<sub>y</sub>O<sub>3</sub> Al:SrTiO<sub>3</sub> RhCrO<sub>3</sub>:Al:SrTiO<sub>3</sub>  
Cr(NO<sub>3</sub>)<sub>3</sub>

Category: SrTiO<sub>3</sub>

Status: Done

Created by: Nadzeya Brezhneva

## Reaction scheme/sample structure



ChemDraw File (linked): [NB-189-RhCr-oxide-loaded-Al-SrTiO3.cdxml](#)

## Literature/reference experiments

|                     |  |
|---------------------|--|
| Literature          | <a href="https://doi.org/10.1039/C9EE00310J">https://doi.org/10.1039/C9EE00310J</a>  |
| Reproduction        | /  |
| Similar experiments | <a href="#">SrTiO<sub>3</sub> - NB-323: Modification of EA-354 (SrTiO<sub>3</sub>:Al, upscaled batch) with Rh, Cr oxide cocatalyst (0.2 wt%)</a> |

## Reagents

| Name  | CAS Number / Experiment Number   | Inventory number | Amount    | Equivalents               | Mass <sub>theo</sub> [mg] | Mass <sub>exp</sub> [mg] | Molar mass [g/mol] | Density [g/ml] | Volume [ul]         |
|---|--|------------------|-----------|---------------------------|---------------------------|--------------------------|--------------------|----------------|---------------------|
| Al:SrTiO <sub>3</sub> (EA-354)  | <a href="#">SrTiO<sub>3</sub> - EA-354: Preparation of Al:SrTiO<sub>3</sub> (using EA-352-SrTiO<sub>3</sub>-molten-salt and SrCl<sub>2</sub> as flux) at 1000°C, 10 h, Osterloh route, upscaling x15, filtration with PVDF</a> | /                | 0.83 mmol | 1                         | 150                       | 152.05                   | 183.49             | /              | /                   |
| RhCl <sub>3</sub> solution (9.87 mM in H <sub>2</sub> O)                  | <a href="#">Prep work - NB-KRA-340: Preparation of RhCl<sub>3</sub> and Cr(NO<sub>3</sub>)<sub>3</sub> stock solutions</a>   | /                | 2.92 umol | 0.2 wt% Rh <sup>[a]</sup> | /                         | 295                      | 209.26             | /              | 295                 |
| Cr(NO <sub>3</sub> ) <sub>3</sub> solution (28.81 mM in H <sub>2</sub> O) | <a href="#">Prep work - NB-KRA-340: Preparation of RhCl<sub>3</sub> and Cr(NO<sub>3</sub>)<sub>3</sub> stock solutions</a>   | /                | 5.76 umol | 0.2 wt% Cr <sup>[b]</sup> | /                         | 200                      | 238.011            | /              | 200                 |
| milli-Q water   | /  | /                | 0.19 mol  | /                         | /                         | /                        | 18                 | 0.998          | 3505 <sup>[c]</sup> |

[a] - Rh content in the final product



- [b] - Cr content in the final product
- [c] - for suspension preparation

## Excel sheet for reagent calculation

NB-342-calculations.xlsx

## Furnace Parameters

Equipment - Muffelofen L3/11/P320, Nabertherm GmbH, Lab 106, CEEC I, (Matilda)

### Temperature/time parameters

|  |           |
|--|-----------|
| Used zone or charge sensor                                     | Zone      |
| Used delayed start   | /         |
| Used automatic/manual/extended holdback                        | Automatic |
| The temperature band entered for manual/extended holdback (°C) | /         |
| End time [min], relative to start of program                   | /         |

### Segments

|                | Target Temperature (°C) | Duration (h) | Rate (°C/h) | Temperature band (°C) | Description of the segment | Observations |
|----------------|-------------------------|--------------|-------------|-----------------------|----------------------------|--------------|
| First segment  | 350                     | 0:33         | 600         | /                     | Increase                   | /            |
| Second segment | 350                     | 1:00         | /           | /                     | Hold                       | /            |
| End segment    | 0                       | /            | /           | /                     | Natural cooling down       | /            |



## Procedure/observations

For transfer of precise liquid amount, Eppendorf pipettes were used (for 10 - 100  $\mu$ L: 100  $\mu$ L Eppendorf, above 100  $\mu$ L: 1000  $\mu$ L Eppendorf).

| Date       | Time        | Step  | Observations                         | Pictures  |
|------------|-------------|---|--------------------------------------|---|
| 11.11.2025 | 17:35       | Weighing <b>EA-354</b> , transfer to a 25 mL Schott glass beaker.   | Slightly creamy solid.               | /   |
|            | 17:40       | Placing the glass stirring bar (8 mm) inside the beaker.  | /                                    | /   |
|            | 17:50       | Addition of 3505 $\mu$ L of $\text{H}_2\text{O}$ ( $3 \times 1000 + 505$ ) to the beaker.                       | Creamy, slightly pinkish suspension. | /   |
|            | 17:51       | Placing the beaker with the suspension into the water bath (currently, at room T).                              | /                                    | /   |
|            | 17:52       | Switching the stirring (400 rpm).   | /                                    | /   |
|            | 17:53       | Drop-by-drop addition of 295 $\mu$ L of 9.87 mM $\text{RhCl}_3$ solution to a beaker under stirring.            | /                                    | /   |
|            | 17:54       | Drop-by-drop addition of 200 $\mu$ L of 28.9 mM $\text{Cr}(\text{NO}_3)_3$ solution to a beaker under stirring. | /                                    | /   |
|            | 17:55       | Switching the heating mode on (70 $^\circ\text{C}$ , precise mode).   | /                                    | /   |
|            | 17:55-19:30 | Evaporating of water under constant stirring at 70 $^\circ\text{C}$ .   | /                                    | <a href="#">20251111_175528-start of evaporation of water.jpg</a> |



|            |       |  |  |  |
|------------|-------|--|--|--|
|            | 19:30 | Stop stirring, since all water has been evaporated, removing the stirring bar from the beaker.   | In general solid seems quite homogeneous   | <a href="#">20251111_193204-after evaporation of water.jpg</a>   |
|            | 19:40 | Transfer of the solid to the quartz crucible using Smartspatula, covering the crucible with lid and Al foil before further use.  | /  | <a href="#">20251111_195700-crucibles with dried samples-top 0.2 wt% loading bottom 0.05 wt% loading.jpg</a> |
|            | 20:00 | Loading the crucible with the dried sample (together with the crucible with the sample <a href="#">SrTiO<sub>3</sub> - NB-341: Modification of EA-354 (SrTiO<sub>3</sub>:Al, upscaled batch)</a> with Rh, Cr oxide cocatalyst (0.05 wt%), fresh solutions of <a href="#">RhCl<sub>3</sub></a> and <a href="#">Cr(NO<sub>3</sub>)<sub>3</sub></a> in a muffle furnace Equipment - Muffelofen L3/11/P320, Nabertherm GmbH, Lab 106, CEEC I, (Matilda). | /  | <a href="#">20251111_200317-crucibles inside the furnace.jpg</a>   |
|            | 20:03 | Start of the program (350 °C, 1 h, 10 °/min).  | /  | /  |
| 12.11.2025 | 0:30  | Removal of the samples from the furnace.   | T = 163 °C when opening the lid<br>Grey solid (more intense grey colour in comparison with <a href="#">SrTiO<sub>3</sub> - NB-321: Modification of EA-354 (SrTiO<sub>3</sub>:Al, upscaled batch)</a> with Rh, Cr oxide cocatalyst (0.05 wt%)). | /  |
|            | 15:50 | Collecting the sample from the beaker, transfer to a 5 mL vial with a snap cap, weighing.<br>m = 142.57 mg   | <b>NB-342-0.2 wt% Rh,Cr-oxide loaded Al:SrTiO<sub>3</sub></b><br>Creamy greyish solid.   | <a href="#">20251112_160416-final product.jpg</a>  |

## Product characterization

| Sample  | Mass [mg] | Purity | Mass <sub>pure</sub> [mg] | Amount [mol] | Yield [%] | Description           | Image   | Storage location                |
|---|-----------|--------|---------------------------|--------------|-----------|-----------------------|---|---------------------------------|
| <b>NB-342-0.2 wt% Rh,Cr-oxide loaded Al:SrTiO<sub>3</sub></b> | 142.57    | /      | /                         | /            | 93.76     | Creamy greyish solid. | <a href="#">20251112_160416-final product.jpg</a> | Lab E004 (CEEC II), SSC cabinet |

## Results

Modification of Al:SrTiO<sub>3</sub> sample (**EA-354**, upscaled batch) with Rh, Cr oxide cocatalyst (0.2 wt% Rh and Cr in final product) was performed. The final product (**NB-342-0.2 wt% Rh,Cr-oxide loaded Al:SrTiO<sub>3</sub>**) is



represented with grey solid, m = 142.57 mg, yield 93.76 %.

## Linked experiments

Prep work - [NB-340: Preparation of  \$\text{RhCl}\_3 \cdot 3\text{H}\_2\text{O}\$  and  \$\text{Cr}\(\text{NO}\_3\)\_3 \cdot 9\text{H}\_2\text{O}\$  stock solutions](#)

SrTiO<sub>3</sub> - [NB-162: Preparation of SrTiO<sub>3</sub>:Al with RhyCr<sub>2</sub>-yO<sub>3</sub> \(Osterloh route\), I attempt, 1000 C 10h](#)

SrTiO<sub>3</sub> - [NB-175: Preparation of SrTiO<sub>3</sub>:Al \(1100 C, 3 h, IAAC furnace\) \(Osterloh route\)](#)

SrTiO<sub>3</sub> - [NB-285: Preparation of SrTiO<sub>3</sub>:Al \(from self-made SrTiO<sub>3</sub>, NB-283, batch V, Osterloh \(no Al<sub>2</sub>O<sub>3</sub>\), upscaling x3\), 1000 C, 10 h, filtration with PVDF](#)

SrTiO<sub>3</sub> - [NB-289: Modification of NB-285 \(SrTiO<sub>3</sub>:Al \(from self-made SrTiO<sub>3</sub>, Osterloh, no Al<sub>2</sub>O<sub>3</sub>, PVDF filter\) 1000 C, 10 h\) with Rh, Cr oxide cocatalyst](#)

SrTiO<sub>3</sub> - [NB-301: Modification of NB-285 \(SrTiO<sub>3</sub>:Al \(from self-made SrTiO<sub>3</sub>, Osterloh, no Al<sub>2</sub>O<sub>3</sub>, PVDF filter\) 1000 C, 10 h\) with Rh, Cr oxide cocatalyst I](#)

SrTiO<sub>3</sub> - [EA-354: Preparation of Al:SrTiO<sub>3</sub> \(using EA-352-SrTiO<sub>3</sub>-molten-salt and SrCl<sub>2</sub> as flux\) at 1000°C, 10 h, Osterloh route , upscaling x15, filtration with PVDF](#)

SrTiO<sub>3</sub> - [EA-359: Modification of Al:SrTiO<sub>3</sub> \(EA-354\) via deposition of Rh, Cr oxide co-catalyst, 350°C, 1h](#)

SrTiO<sub>3</sub> - [NB-321: Modification of EA-354 \(SrTiO<sub>3</sub>:Al, upscaled batch\) with Rh, Cr oxide cocatalyst \(0.05 wt%\)](#)

SrTiO<sub>3</sub> - [NB-323: Modification of EA-354 \(SrTiO<sub>3</sub>:Al, upscaled batch\) with Rh, Cr oxide cocatalyst \(0.2 wt%\)](#)

SrTiO<sub>3</sub> - [NB-341: Modification of EA-354 \(SrTiO<sub>3</sub>:Al, upscaled batch\) with Rh, Cr oxide cocatalyst \(0.05 wt%\), fresh solutions of  \$\text{RhCl}\_3 \cdot 3\text{H}\_2\text{O}\$  and  \$\text{Cr}\(\text{NO}\_3\)\_3 \cdot 9\text{H}\_2\text{O}\$](#)

## Linked resources

Equipment - [Muffelofen L3/11/P320, Nabertherm GmbH, Lab 106, CEEC I, \(Matilda\)](#)

Equipment - [Muffle furnace Nabertherm LT 15/11/P330 \(AWZ 304 lab\)](#)

Protocol - [Heat treatment using Muffelofen L3/11/P320, Nabertherm GmbH, Lab 106, CEEC I](#)

Protocol - [Muffle furnace Nabertherm GmbH LT 15/11/P330 \(Lab AWZ 304\)](#)



## Attached files

NB-342-calculations.xlsx

sha256: 43615633e81e5fc207f524517a456ba718f3b9c7e5a303b75ce8dd3c8abf2de3

20251112\_160416-final product.jpg

sha256: 5f20259581609a2d9ac5907be5640940415dd5137e4e692776fde454a751b3c2



20251111\_175528-start of evaporation of water.jpg

sha256: 9f27b47ac30d281dfd813cf1372af0af681504898eef97a35e48c06c98befbf3



20251111\_193204-after evaporation of water.jpg

sha256: 7646c895d1ce788dfa07845b2aa855c73d1946bf8e7852fbe4991910a25676c9



20251111\_200317-crucibles inside the furnace.jpg

sha256: c98bb38e9e64b141bec82adca182e2da1f5ecbc491268eacd21840a206c5dacd



20251111\_195700-crucibles with dried samples before calcination-top 0.2 wt% loading bottom 0.05 wt% loading.jpg

sha256: 6040d4408319e3b764c1d2c72439667a93b17fc239b3f64af137ef3772fc68ab



Unique eLabID: 20251111-34c23151f470164172068e0e106c0106ed79c314  
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3422>